

## Installation Guide

### Philips MP2 Post Rail Clamp (PRC) Mounting Kit


The purpose of this guide is to:

1. Describe mounting of MP2 on post (pages 1 - 2)
2. Describe mounting of MP2 on rail (page 3).

#### Parts Reference

The following parts and hardware are included with this installation kit (hardware not shown):


Item #	Description	Qty
1	Post Rail Clamp (PRC)	2
2	MP2 Interface Plate (Post Mounting Only)	1
3	M4 x 16 mm Flat Head Machine Screw (FHMS)	2
4	M4 x 14 mm Pan Head Machine Screw (PHMS)	2
5	M4 x 8 mm PHMS	4



**Tools Required:** Phillips screwdriver (not provided).

#### Mounting Docking Station and Power Supply on Post

1. Fasten Interface Plate to docking station with two (2) M4 x 14 mm FHMS as shown below.




## Mounting Docking Station and Power Supply on Post

1. Fasten PRC to Interface Plate with two (2) M4 x 8 mm PHMS as shown below.



1. Fasten PRC to power supply with two (2) M4 x 8 mm PHMS as shown below.





## Mounting Devices on Post

1. Clamp devices onto post as shown and tighten Clamping Knobs. Mount monitor according to Philips instructions.

## Mounting Docking Station and Power Supply on Rail

1. Fasten PRC to docking station with two (2) M4 x 14 mm PHMS as shown below left.
2. Fasten PRC to power supply with two (2) M4 x 8 mm PHMS as shown below right.



## Mounting Devices on Rail

1. Clamp devices onto rail as shown and tighten Clamping Knobs. Mount monitor according to Philips instructions.



## **Routine Maintenance**

Periodically check all mounting hardware. Tighten as necessary for optimal operation and safety.

## **Cleaning the Mounting Assembly**

**CAUTION:** GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted instruments or accessory equipment, refer to the specific instructions delivered with those products.

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.
3. Do not use steel wool or other abrasive material to clean the mounting assembly.
4. Damage caused by the use of unapproved substances or processes will not be covered by warranty. We recommended that you test any cleaning solution on a small area of the mounting assembly that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the mounting assembly. Wipe any cleaning agents off the mounting assembly immediately using a water-dampened cloth. Dry all mounting assemblies thoroughly after cleaning.