About the "Children's Respiration Monitor"

The "Children's Respiration Monitor" system comprises of an accelerometer based device, with a belt to attach the device to the child. There are two variants, Community health worker model called CHW and a Professional Health worker model called PHW. This Philips product has to be used and operated only in accordance with the safety procedures and operating instructions provided in this "Instructions for Use," and for the purposes for which it is designed. This product is intended for the below mentioned purposes. However, nothing stated in this "Instructions For Use" reduces the responsibilities of the operator for sound clinical judgment and use of the best clinical procedures.

Intended Use

The "Children's Respiration Monitor" measures the respiration rate in children under five years old and automatically classify fast breathing rate according to the IMCI guidelines set by the "World Health Organization". This device is to be used by health workers at the community level in low-resource settings, clinical officers, nurses, midwives, clinicians (Professional Health Workers) at primary or secondary care facilities. This device is neither intended to provide automated treatment decisions, nor to be used as a substitute for professional healthcare judgment. All patient medical diagnosis and treatment has to be performed under the supervision and oversight of an appropriate healthcare professional.

Conventions Used in this Guide

This guide uses the following conventions

- "CAUTION"** to be mandatorily exercised.
- (P) Applicable only to Professional Health Worker model (PHW).
- (C) Applicable only to Community Health Worker model (CHW).

Safety Summary


Symbols Marked on the Device and Packaging

Symbol	Description
	The Legal manufacturing address is next to this symbol.
	Input DC voltage rating.
	Indicates the need for separate collection for electrical and electronic equipment in compliance with the Waste Electrical and Electronic Equipment (WEEE) directive.
	Date of manufacture.
	For Medical Device directive. See the Children's Respiration Monitor "Instructions for Use".
	Type B Applied Part.
	Dust-protected, Protected against splashing water.
	The number next to this symbol is the serial number.
Pxxxxxx	Serial number for PHW.
Cxxxxxx	Serial number for CHW.

Important Safety Information

- This Device has to be used by trained and authorized personnel only.
- This system should only be used on FULL TERM babies (i.e. babies completed 37 weeks of gestation).
- The system shall not be used on Children wearing any supportive device at the area of the chest/belly.
- The system does not need direct skin contact for measurement, the user is recommended to attach directly on cloth. If used without a cloth use on intact skin only.
- Use only the belt supplied with the system, use of any other belt is not allowed.
- Be cautious while attaching the Belt to ensure that no injuries are caused to the child.
- Ensure that the belt is not too tightly attached to the child.
- Do not connect the Charging Cable to the device when it is attached to the child.
- Charge the PHW device to full before first use. Always charge the device to full battery state for correct battery indication and functioning.
- Do not leave the system unattended while not in use.
- Ensure the system is not exposed to unhygienic conditions.
- The device or the belt should not be sterilized in any sterilizer equipment.
- The concerned device user is responsible for proper disposal of the device past their expected life span or when no longer be used for its intended purpose, dispose in accordance with applicable waste disposal laws and local regulations.

Getting Started

Frequently used functions

- Attach the device to the child using the belt (Refer to the "Quick Reference Guide" (QRG), Fig- 1).
- Power "ON" the Device (Refer to the QRG Fig- 2).
- Select an Age Group (Refer to the QRG Fig-3).
- Start a Measurement (Refer to the QRG Fig-4)
- Power "OFF" the Device (Refer to the QRG Fig-2).
- Charge the Device (P) (Refer to the QRG for a pictorial representation).
- Download the Data (P) Refer to this Link (<http://www.philips.com/a-w/about/sustainability/healthy-people/supporting-communities/fabric-of-africa/programs/childrens-automated-respiratory-monitor.html>) under Tools and Resources to find the detailed instructions.

Instructions for Disinfection

1. Turn off the device before disinfection.
2. Wipe the patient contact surface of the device and the belt (**as shown below**) before every use with a clean cloth moistened with 70% Isopropyl alcohol.

Routine Cleaning Instructions

1. Turn off the device before cleaning.
2. Wipe the external surfaces of the device and the belt with a cloth dampened in mild soap solution or water.

Troubleshooting Guide

- If the unit does not power up after repeated attempts, and you just see a blank screen, contact your local Philips representative.
- Displayed at end of measurement indicates that either there was excess motion or loosely attached belt. Ensure that the measurement conditions are followed for accurate results.
- Appearing on the screen indicates that the accelerometer is faulty, contact your local Philips representative.
- This icon indicates that the battery is low please charge the device in case of PHW, order a new device in case of CHW.
- If the "RED LED" on the device appeared and the unit powered down while turning on, this is due to excess temperature. Allow the unit to cool and try powering ON after 30 minutes. If the issue persists, contact your local Philips representative. If the RED LED appeared after power ON or while initiating a measurement, it is due to critically low battery. Charge the device to full in case of PHW. Order a new device in case of CHW.

Technical Specification

Meter Type: Breathing Rate Measurement.

Principle: Uses a 3D accelerometer with advanced signal processing to measure the specific abdominal motions related to breathing.

Display: Liquid Crystal Display.

Accuracy: ± 2 Breaths per minute when measured under the recommended conditions.

Measurement Range: 10 to 150 Breaths per minute.

Completed Measurement: Beep sound and LED indication.

Mode of Operation: Non-Continuous.

Type of Protection: Class II (while charging Mode) (P).

Degree of Protection: Type B Applied Part.

Power Specifications for Charging Adapter (P): Input Voltage: 100-240V@50/60 Hz, Output Voltage: 5VDC Output Current: > 0.5A.

Battery Type: Alkaline for CHW (C) and Rechargeable NiMH for PHW (P).

Battery Life: 2000 measurements over two year usage (C).

4380 measurements over three years of usage with 200 measurements per full charge cycle (P).

Automatic Shut-off: Approximately after Two minutes of inactivity.

Memory: 300 measurement records (P).

Degree of Ingress Protection: IP54 rated.

Operating Conditions: 0 °C to 40 °C, 10% to 90% RH.

Atm Pressure (Operating / Storage): 700 hPa – 1060 hPa.

Storage Conditions (Packaged): -20 °C to 50 °C, 10% to 95% RH.

Dimensions (approx.): 65 x 65 x 28 mm (without belt).

Weight (approx.): 95 gms including belt.

Safety Compliance: IEC 60601-1-Ed3.1, IEC 60601-1-2 Ed 3, IEC 60601-1-11 Ed1.

Guidance & Manufacturer's Declaration, according to IEC 60601-1-2:2007.

Do not use mobile RF communication equipment, which generate strong electrical or electromagnetic fields, near the medical device. The user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF equipment (transmitters) and the Children's Respiration Monitor as recommended below.

Rated maximum output power of Transmitter	Minimum Safe Distance
0.01 Watt	0.23 meter
0.1 Watt	0.73 meter
1 Watt	2.3 meter
10 Watt	7.3 meter
100 Watt	23 meter

The CE mark on the product denotes compliance with all applicable EU Directives. Note that the Notified Body number does not apply to RoHS Directive (restriction of the use of certain hazardous substances in electrical and electronic equipment).

REACH Requirements

REACH requires Philips Healthcare (PH) to provide chemical content information for Substances of Very High Concern (SVHC) if they are present above 0.1% of the product weight. Components within electric and electronic equipment may contain phthalates above the threshold (e.g.bis (2-ethylhexyl)phthalate), CAS nr.: 117-81-7). The SVHC list is updated on a regular basis. Please refer to the following Philips REACH website for the most up-to-date list of products containing SVHC above the threshold: <http://www.philips.com/about/sustainability/reach.page>

www.philips.com/healthcare

healthcare@philips.com

Philips Medical Systems Nederland B.V.
Veenpluis 4-6 5684 PC Best, The Netherlands

4598 010 77026 / August 2017

Manufactured by

Philips India Limited, Plot No.B-79, MIDC, Phase - II, Chakan, Taluka-Khed Village-Savardari, District: Pune - 410501, Maharashtra, India

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