

PHILIPS

www.philips.com

Transitioning Mechanical Ventilation across the care continuum: The Hospital to Home story

Robert S Campbell, RRT FAARC

HRC

5-22-2019

innovation  you



Objectives

- Describe the Why, Who, and How of Home Ventilatory Support
- Understand the clinical applications of AVAPS
- Understand the clinical and hospital value of AVAPS-AE
- Understand our technologies and the COPD patient journey



Why Provide Home Ventilatory Support

- Economics
 - Up to 70% cost savings compared to MV in ICU
 - Cost of home MV ranges from \$3,500-8,500 per month
- Improved Quality of Life indicators
- Enhanced social interactions
- Reduced infection rate
- Alternate sites include: LTACH, Skilled nursing facilities, Home sites



Who is a candidate for home MV?

- COPD patients (typically Gold stage III and IV)
- Neuromuscular patients (ALS, DMD, SMA, etc.)
- Chronic Hypoventilation syndrome patients



Transition key components


- Device Selection
- Interface Selection
- Mode Selection
- Settings optimization/titration
- Alarm management
- Oxygen needs and delivery options



Transition key components (cont.)

- Humidification
- Secretion management
- Patient and caregiver education/training/troubleshooting
- Case management to incorporate and organize multidisciplinary approach
- Home provider
- Reimbursement issues

Trilogy100 Circuit Options – Passive



- Utilizes an Exhalation Port
 - Integrated into a mask
 - Whisper Swivel II
- Passive Porting Block
- Choose Passive in Set up Menu



Failed approach



ED - ICU

General ward

Home

Repeat



Example: Hospital impacted by readmissions

Balance Sheet | Income Statement | Financial Ratios | Charges/Costs by Dept | Medicare Incentive Programs **New** | Payor Mix | Beds/Components | Open Payments

Export Medicare Incentive Programs to Excel

Medicare Incentive Programs ⓘ

REVENUES	FY2013	FY2014	FY2015	FY2016	FY2017	FY2018
Est. Revenue Loss Due to Readmission Penalty	\$0	\$0	(\$141,881)	(\$706,018)	(\$714,843)	(\$405,960)
Est. Revenue Adjustment Due to Value Based Purchasing	\$204,103	\$180,598	(\$144,317)	(\$105,068)	(\$436,821)	(\$358,201)
Est. Revenue Loss Due to Hospital Acquired Condition Penalty			\$0	\$0	(\$882,523)	(\$882,523)



Transitional care approach

Readmission prevention begins in the hospital




AVAPS-AE



AVAPS-AE

Shown to reduce readmission.* **Start early in the hospital.**

COPD Transitional Care Story



Reduce or eliminate unnecessary readmission



Patient journey

AECOPD presents to the ED



Presentation: High work of breathing, shortness of breath, discomfort, tachycardia, hypotension or hypertension, anxious, sense of 'suffocation' or 'drowning,' erratic/chaotic breathing

Place on V60, S/T or AVAPs

Assessment is vital. Rapid changes to NIV settings are based on:

- ABG results
- Blood pressure
- Patient comfort and WOB
- Sudden changes in compliance (V_T)
- May need to make changes to IPAP, EPAP, FiO_2 , Rise time



Benefits of Philips V60 in this scenario:

- Clinicians can make changes to S/T settings quickly
- Large screen for monitoring and alarm management

If hemodynamically stable, transfer to ICU

Example: Transferred to ICU with V60, S/T



Upon patient and NIV stabilization consider AVAPS



- Vt
- EPAP
- Rate
- FiO₂
- Min P and Max P
- Rise time

Mode:	S/T
IPAP:	15 cmH ₂ O
EPAP:	7 cmH ₂ O
V _T measured:	~450mL
Rate:	10 bpm

Mode:	AVAPS
Min P:	12 cmH ₂ O
Max P:	22-25 cmH ₂ O
EPAP:	7 cmH ₂ O
V _T (set):	6-8 mL/Kg
Rate:	10 bpm



Example



Start-up algorithm

AVAPS will choose from the highest of these three:

- 1 IPAP Minimum: In this example, IPAP minimum = $15\text{cmH}_2\text{O}$
- 2 $(V_T/60\text{ml/cmH}_2\text{O}) + \text{EPAP}$: $(450\text{mL}/60\text{ml/cmH}_2\text{O}) + 7 = 14.5\text{cmH}_2\text{O}$
- 3 $8\text{cmH}_2\text{O} + \text{EPAP}$: $8\text{cmH}_2\text{O} + 7\text{cmH}_2\text{O} = 15\text{cmH}_2\text{O}$



Benefits of AVAPS in ICU



Clinical

- Lowest pressure support to achieve target volume
- Improves patient comfort and compliance¹
- Less sedation¹
- Improved Health-Related Quality of Life (HRQL)³
- Improves gas exchange^{1,2,3}
 - **Small** changes to PaCO₂ lead to **significant** outcomes
- Improved sleep quality³


Workflow

- Smart algorithm with ongoing monitoring and adjustments
- Incremental IPAP changes (up to 2.5 cmH₂O/min) to minimize sudden mask leaks and sleep disturbances
- Clinician is in control; AVAPS does the work

1. Canpolat G, Ozgultekin A, Turan G, Iskender A, Adiyek E, Ekinci O. Does average volume-assured pressure support make any difference compared with BIPAP? *Critical Care*. 2014;18(Suppl 1):P265. doi:10.1186/cc13455.
2. Windisch, Wolfram. "Average Volume-Assured Pressure Support in Obesity Hypoventilation." *CHEST* 130 (2006): 815. With respiratory insufficiency patients diagnosed with Obesity Hypoventilation Syndrome
3. Murphy, PBTorax thoraxjnl-2011-201081: Published Online First: 1 March 2012 doi:10.1136/thoraxjnl-2011-201081 Produces results comparable to sleep lab titration of PS

V30 is the cornerstone in this story

For a successful transition home Philips solutions enable consistent modes, features and patient experience for an easy return home.



Reduce or eliminate unnecessary readmission



Clinical application of AVAPS-AE

Treatment of COPD



COPD and readmission: AVAPS-AE

- Maintain stable and secure ventilation
- Lowest pressure support to achieve target volume
- Adjusting EPAP to maintain open airway
- Identifies and retains PS requirements

Obstructive sleep apnea: Auto-CPAP and Auto-bilevel

- Reverse upper airway obstruction

Prevent respiratory failure: AVAPS-AE

Transitional care approach

Readmission prevention begins in the hospital



AVAPS-AE



AVAPS-AE

Shown to reduce readmission.* **Start early in the hospital.**

Change the cycle of readmissions

Introduce AVAPS-AE in the hospital

Patient education

Pulmonary rehabilitation

Multidisciplinary care team

Medication reconciliation

Follow-up plan prior to discharge

Device / therapy reconciliation

Remote ventilator monitoring with

Care Orchestrator



Chronic Obstructive Pulmonary Diseases: Journal of the COPD Foundation



National COPD Readmissions Summit

The 2nd National COPD Readmissions Summit and Beyond: From Theory to Implementation

Kristen S. Willard, MS¹ Jamie B. Sullivan, MPH¹ Byron M. Thomas, MD²
Catherine S. Jones, PhD, APRN, ANP-C³ Leonard Fromer, MD⁴ Barbara P. Fawcett, MD, MSc⁵
Alpesh Amin, MD, MBA, MACP⁶ Jean M. Rommes, PhD⁷ Rhonda L. ...


Experts in attendance identified the best available evidence-based approaches to improving care delivery across the continuum, including ongoing device reconciliation

Abstract

Chronic obstructive pulmonary disease (COPD) hospitalizations and readmissions adversely impact the health and quality of life of COPD patients. Under the Hospital Readmissions Reduction Program, the Centers for Medicare & Medicaid Services reduce payments to those hospitals exceeding expected rates of COPD readmissions within 30 days of hospital discharge. It was within this climate that the COPD Foundation held its 2nd COPD Readmissions Summit in March 2015. Experts in attendance: (1) categorized challenges to optimal COPD care, (2) analyzed the state of care delivery and readmissions reduction strategies and (3) identified the best available evidence-based approaches to improving care delivery across the continuum including early diagnosis via spirometry, ongoing device, oxygen and medication reconciliation, treatment that addresses comorbidities and preventive care, robust patient education, prompt post-acute follow up, home health services and pulmonary rehabilitation. Results of this collaborative event formed the basis for PRAXIS, the COPD Foundation's initiative to improve COPD care across the health continuum and to reduce readmissions.

A culture of readmission reduction

Medication reconciliation



Family



Physician



Respiratory therapist



RN



Discharge planner

Therapy reconciliation



Family



Physician



Respiratory therapist




RN



Discharge planner

Example



AVAPS

AVAPS-AE

Mode	AVAPS-AE
Trigger Type	Auto-Trak
AVAPS Rate	5.0 cmH ₂ O/min
Tidal Volume	450 ml
Maximum Pressure	30.0 cmH ₂ O
Pressure Support Max	11.0 cmH ₂ O
Pressure Support Min	5.0 cmH ₂ O
EPAP Max Pressure	12.0 cmH ₂ O
EPAP Min Pressure	7.0 cmH ₂ O
Breath Rate	10 BPM
Inspiratory Time	1.0 seconds

Buttons: Finish, Navigate, Modify

Transitioning home



Increasing **Pressure Support Max** allows algorithm autonomy

After the initial session the device will remember the last effective PS and start at it

▲ Mode	AVAPS-AE
Trigger Type	Auto-Trak
AVAPS Rate	5.0 cmH2O/min
Tidal Volume	450 ml
Maximum Pressure	30.0 cmH2O
Pressure Support Max	11.0 cmH2O
Pressure Support Min	5.0 cmH2O
EPAP Max Pressure	12.0 cmH2O
EPAP Min Pressure	7.0 cmH2O
Breath Rate	10 BPM
▼ Inspiratory Time	1.0 seconds

Finish Navigate ▲▼ Modify

Transitioning home



▲ Mode	AVAPS-AE	
Trigger Type	Auto-Trak	
AVAPS Rate	5.0 cmH2O/min	
Tidal Volume	450 ml	
Maximum Pressure	30.0 cmH2O	
Pressure Support Max	11.0 cmH2O	
Pressure Support Min	5.0 cmH2O	
EPAP Max Pressure	12.0 cmH2O	
EPAP Min Pressure	7.0 cmH2O	
Breath Rate	10 BPM	
▼ Inspiratory Time	1.0 seconds	
Finish	Navigate ▲▼	Modify



▲ Mode	AVAPS-AE	
Trigger Type	Auto-Trak	
AVAPS Rate	5.0 cmH2O/min	
Tidal Volume	450 ml	
Maximum Pressure	30.0 cmH2O	
Pressure Support Max	11.0 cmH2O	
Pressure Support Min	5.0 cmH2O	
EPAP Max Pressure	12.0 cmH2O	
EPAP Min Pressure	7.0 cmH2O	
Breath Rate	10 BPM	
▼ Inspiratory Time	1.0 seconds	
Finish	Navigate ▲▼	Modify

Transition home

- AVAPS-AE with Trilogy
- Care Orchestrator for wireless remote ventilator monitoring
- RT driven disease management programs



Versatile, flexible, proven





Discharge Planning

- Discharge planning starts on admission!
- Case manager is central to organizing care and education (COPD Navigator)
- Multidisciplinary approach
- Transition within hospital under controlled environment, but with home equipment and support
- Patient self-care education and disease management is also key to success
- Remote monitoring and follow up care

