



**PHILIPS**


Complex PCI

Your partner in complex PCI:  
In-stent restenosis (ISR)



# Your partner in complex PCI: In-stent restenosis (ISR)

Philips provides a portfolio of specialty coronary diagnostic and therapy devices that enable safe and effective treatment of a wide variety of the most complex coronary lesions types and morphologies, included in-stent restenosis.





**ISR is >50% diameter stenosis at stent's inside or edges, with different ISR morphologies predicting different TLR rates at 1-year<sup>1,2</sup>**


Several factors with both biological and/or procedural causes contribute to the nature of ISR<sup>3</sup>

**Biologic causes**

- Reaction to metal or polymer
- Drug resistance
- Thrombosis

**Procedural causes**

- Stent under-expansion/mal-apposition
- Stent fracture
- Edge trauma
- Geographical miss




Malapposition (\*) occurs at the junction of calcified and noncalcified plaque and in localized areas of reference segment ectasia at the stent edges. However, the stent is fully (albeit not symmetrically) expanded since the MSA of 6.0 mm<sup>2</sup> matches the proximal and distal reference segment lumen areas.

**Incidence of ISR (meant PCI for DES restenosis) quite consistent over either National Registries, All-comers RCTs and Registries, in the range of ≈5%**


Source	PCI	Design	Incidence	Notes
Cassese et al. <sup>4</sup>	10.004	Routine angiographic surveillance after unrestricted use of newer-generation devices	30,1% BMS 14,6% 1st gen DES 12,2% 2nd gen DES	Angio-binary restenosis
UK BCIS Audit <sup>5</sup>	96.143	Registry of the British Cardiovascular Intervention Society	4,8%	PCI for restenosis, overall incidence (85% DES)
Spanish Registry <sup>6</sup>	67.671	National Registry of Coronary and Structural Interventions (2010-2015)	4,7%	PCI for restenosis, could include DES, BMS, multi-layers
Norstent Study <sup>7</sup>	4.504	Multicentric RCT 1:1 DES Vs. BMS	4,6%	PCI for DES restenosis
RESOLUTE All-Comers <sup>8</sup>	2.292	Patients randomly assigned to R-ZES (n = 1,140) or EES (n = 1,152)	7,0% ZES, 6,5% EES	Clinically driven TLR
SCAAR Registry <sup>9</sup>	>25.000	SCAAR/SWEDHEART records consecutive patients from all centres (n=29) performing PCI in Sweden	5,6%	PCI treated in-stent restenoses

# Physiology outperforms angiography in predicting functional significance of ISR, IVUS helps identify location and etiology of ISR to quickly determine and confirm best treatment options




IVUS helps you verify and optimize your results


Eagle Eye Platinum digital IVUS  
Refinity ST rotational IVUS catheters



SyncVision with iFR Co-reg




AngioSculptX drug-coated scoring balloon



AngioSculpt PTCA scoring balloon



SyncVision with IVUS Co-reg



### **Verrata Plus pressure guide wire**

Plan your procedure using iFR Co-Registration with SyncVision providing physiologic guidance, discriminating focal Vs. diffuse disease.

### **Eagle Eye Platinum digital IVUS – Refinity ST rotational IVUS catheters**

With SyncVision with IVUS Co-Registration easily assess for geographic miss and edge complications, confirm stent apposition and optimal expansion for luminal gain.

### **AngioSculpt PTCA Scoring Balloon**

Lesion preparation with AngioSculpt outperforms POBA for the treatment of DES-ISR with drug-coated balloons.<sup>11</sup>

### **ELCA coronary laser atherectomy**

Effective and safe plaque de-bulking for greater vessel and stent expansion.<sup>12</sup>

By modifying the plaque even behind the struts laser makes it more amendable to further stent expansion.<sup>12</sup>

### **AngioSculptX drug-coated scoring balloon**

AngioSculptX is the first and only Treatment Solution combining Plaque Scoring and Drug Delivery in a Single Device<sup>13,14</sup> for a safe and effective treatment of ISR.



1. Mehran R, Dangas G, Abizaid AS, et al. Angiographic patterns of in-stent restenosis: classification and implications for long-term outcome. *Circulation* 1999;100:1872–8.
2. Solinas E, Dangas G, Kirtane AJ, et al. Angiographic patterns of drug eluting stent restenosis and one-year outcomes after treatment with repeated percutaneous coronary intervention. *Am J Cardiol* 2008;102:311–5.
3. Mintz GS. “Clinical Utility of Intravascular Imaging and Physiology in Coronary Artery Disease”. *J Am Coll Cardiol* 2014;64:207–22.
4. Cassese S, Byrne RA, Tada T, et al. Incidence and predictors of restenosis after coronary stenting in 10,004 patients with surveillance angiography. *Heart* 2014;100:153–9.
5. BCIS Audit 2014 for web. <https://www.bcis.org.uk/education/bcis-auditpresentation-audit-returns-2014/>
6. Garcia del Blanco B, et al. Spain: coronary and structural heart interventions from 2010 to 2015. *EuroIntervention* 2017;13:Z64–Z69.
7. Bona KH, et al. Drug eluting or bare metal stents for coronary artery disease. *NEJM* 2016; 375:1242–1252.
8. Taniwaki M, et al. 4-year clinical outcomes and predictors of repeat revascularization in patients treated with new-generation drug-eluting stents: a report from the RESOLUTE All-Comers trial. *J. Am. Coll. Cardiol* 2014;63(16):1617–1625
9. <http://www.ucr.uu.se/swedeheart/arsrapport-2017/swedeheart-annualreport-2017>.
10. Dangas GD, et al. In-Stent Restenosis in the Drug-Eluting Stent Era. *J Am Coll. Cardiol* 2010;56:1897–907.
11. Kufner S, et al. “Neointimal Modification With Scoring Balloon and Efficacy of Drug-Coated Balloon Therapy in Patients With Restenosis in Drug-Eluting Coronary Stents: A Randomized Controlled Trial”. *JACC Cardiovasc Interv.* 2017 Jul 10;10(13):1332–1340. doi: 10.1016/j.jcin.2017.04.024.
12. Mehran R, et al. “Treatment of In-Stent Restenosis With Excimer Laser Coronary Angioplasty”. *Circulation* 1997 Oct 7;96(7):2183–9.
13. Scheller B, Fontaine T, Mangner N, et al. A Novel Drug-Coated Scoring Balloon for the Treatment of Coronary In-Stent Restenosis: Results From the Multi-Center Randomized Controlled PATENT-C First in Human Trial. *Cath and Cardiovasc Interv.* 2016; 88:51–59.
14. Scheller B, et al. “A Novel Drug-Coated Scoring Balloon for the Treatment of Coronary In-Stent Restenosis Two-Year Results from the PATENT-C First in- Human Trial”. Presented at TCT 2015.

©2018 Koninklijke Philips N.V. All rights reserved.  
Some or all products manufactured by Spectranetics, a Philips company.  
Approved for external distribution. D046657-00 102018