


**PHILIPS**

IGT Devices



# Coronary portfolio

# Decide, Guide, Treat and Confirm

## Decide

### Physiology




iFR and FFR  
modalities



OmniWire  
pressure guide wire




Verrata Plus  
pressure guide wire



iFR Co-registration

### Imaging



Intravascular  
ultrasound (IVUS)




Eagle Eye Platinum  
digital IVUS catheters




Refinity ST  
rotational IVUS catheter

IntraSight Interventional application  
platform (Integrated and Mobile)

## Guide




Intravascular  
ultrasound (IVUS)



Eagle Eye Platinum  
digital IVUS catheters



Refinity ST  
rotational IVUS catheter



iFR and IVUS  
Co-registration



OmniWire  
pressure guide wire

IntraSight Interventional application  
platform (Integrated and Mobile)


## Treat



AngioSculpt RX PTCA  
scoring balloon



ELCA laser  
atherectomy catheter



Quick-Cross  
support catheters\*

## Confirm


### Imaging



Intravascular  
ultrasound (IVUS)



Eagle Eye Platinum  
digital IVUS catheters



Refinity ST  
rotational IVUS catheter

### Physiology



iFR Co-registration



OmniWire  
pressure guide wire



Verrata Plus  
pressure guide wire

Philips Laser System

CVX-300 Excimer Laser System


IntraSight Interventional application  
platform (Integrated and Mobile)

\* Coronary indication for the device is limited and subject to change.  
Products subject to country availability. Please contact your local sales representative.





See beyond the angiogram to help provide superior care



## Co-registration

Combine iFR and IVUS data with the angiogram using Philips exclusive iFR and IVUS Co-registration<sup>6</sup> technology.

4



Configurations/Features	Series 3	Series 5	Series 7
<b>IntraSight Integrated - interventional application platform</b> (includes IntraSight CPU with Windows 10 OS, 19" monitor kit, mouse, keyboard and cabling kit)	●	●	●
<b>Imaging (IVUS) license</b> (includes digital, rotational and ChromaFlo IVUS)	●	●	●
<b>Physiology license</b> (includes iFR Spot/Pullback - Hyperemia free lesion assessment modalities and FFR modality)		●	●
<b>Touch Screen Module (TSM)</b>	●	●	●
<b>Philips Remote Services</b>	●	●	●
<b>IVUS and iFR co-registration/tri-registration*</b> (includes SyncVision CPU, monitor, joystick, mouse, keyboard and cabling kit)			●
<b>Device Detection*</b>			●
<b>Quantitative Coronary Analysis*</b>			●
<b>Vessel Enhancement*</b>			●

\* Applicable for coronary only. Note: SpinVision and printer are optional extra parts

		IntraSight series 3, 5, 7	SyncVision (with IntraSight 7)
Power requirements	System input	100, 120 v, 220, 240 VAC, 50/60 Hz, 1000 VA	100 V-120 V, 50/60 Hz, 220-240 V, 50/60 Hz, 600 VA
	Workstation	100-240 V, 50/60 Hz, 825 VA	100-240 V, 50/60 Hz, 250 VA
	Monitor	100V-240 V 50/60 Hz, 39 W	100-240 V, 50/60 Hz, 93 VA
Dimensions	Workstation	H=43.18 cm, W=25.4 cm, D=41.91 cm, -20 kg	H=41.91 cm, W=17.15 cm, D=54 cm, -15.87kg
	Touch Screen Module (TSM) with articulating tableside mount	H=17.78 cm, W=30.23 cm, D=22.86 cm, -3.6 kg (articulating arm extends to a max depth of 16.5" and/or 20" above the top of the bedrail)	n/a
	Monitor	H=38.1-48.26 cm (adjustable stand), W=40.13 cm, D=24.64 cm, -5.9 kg	H=38.1-48.26 cm (adjustable stand), W=40.13 cm, D=24.64 cm, -5.9 kg
	Joystick	n/a	H=3.81 cm, W=10.67 cm, D=7.62 cm, -0.9 kg
	Connection box	H=25.1 cm, W=7.5 cm, D=19.68 cm, -2.7 kg	n/a
Processing and data storage	Processor	1 CPU with 2.3 GHz (maximum turbo frequency of 3.2 GHz). 12 core total. 2400 MHz BUS	1 GPU P5000 1 CPU Intel E5-1600/E5-2600 Series Processor
	Memory	32 GB RAM	16 GB RAM
	Hard drive capacity	128 GB NVMe SSD, 1 TB SATA SSD	120 GB SSD SATA + 480 GB SSD SATA
	Digital archiving capabilities	Local, DVD/Blu-ray, DICOM Network (includes Worklist management, DICOM Store)	n/a
	USB export files	.jpg	.jpg/avi

## Ordering information

IntraSight 3	IntraSight 5	IntraSight 7
797403 (series 3)	797403 (series 5)	797403/797406 (series 7)

# Intrasight Mobile

Providing versatility, IntraSight can now be experienced on an easy to maneuver mobile cart.

## Designed for all lab settings

The Philips IntraSight on mobile is ideally suited for acute and non-acute settings. Customize your platform, select the best-in-class imaging and physiology tools that are right for your coronary or peripheral vascular patients.

Seamless mobile integration with any interventional suite enabling the use of Philips interventional precision tools.

Mobile 5 Series





## Configurations/Features

### Series 3

### Series 5

**IntraSight Mobile - interventional application platform** (includes IntraSight Panel PC with 19" touch screen monitor, mobile cart base with storage bin, DICOM Network Connection and Windows 10 OS)

**Imaging (IVUS) license** (includes digital, rotational and ChromaFlo IVUS)

**Physiology licence** (includes iFR Spot/Pullback - Hyperemia free lesion assesment modalities and FFR modality)

**Touch Screen Module (TSM)**

Optional

**Philips Remote Services**

Note: SpinVision and printer are optional extra parts

### Series 3

### Series 5

Power requirements	System input	100 V-240 VAC, 50/60 Hz, 250 W	100 V-240 VAC, 50/60 Hz, 250 W
Dimensions	Overall system	H=160.17 cm, W=55.07 cm, D=66.32 cm, 56.5 kg (includes cart, panel PC, IVUS PIM and all necessary cabling)	H=160.17 cm, W=55.07 cm, D=66.32 cm, 66.32 kg (includes cart, panel PC, IVUS PIM, FM-PIM, TSM and all necessary cabling)
	Display	19" diagonal, 1280 x 1024 resolution	19" diagonal, 1280 x 1024 resolution
Processing and data storage	Processor	1 CPU Intel Core i7-7820EQ 3.0 GHz Quad Core (maximum turbo frequency of 3.7 GHz)	1 CPU Intel Core i7-7820EQ 3.0 GHz Quad Core (maximum turbo frequency of 3.7 GHz)
	Memory	16 GB RAM	16 GB RAM
	Hard drive capacity	256 GB NVMe SSD, 1 TB SATA SSD	256 GB NVMe SSD, 1 TB SATA SSD
	Digital archiving capabilities	Local, DVD/Blu-ray, DICOM Network (includes Worklist management, DICOM Store)	Local, DVD/Blu-ray, DICOM Network (includes Worklist management, DICOM Store)
	USB export files	.jpg	.jpg

## Ordering information

### IntraSight Mobile 3

797415 (series 3)

### IntraSight Mobile 5

797415 (series 5)



# Disposables compatibility

	Applications	IVUS compatible catheters	Compatible pressure guide wires
<b>IntraSight 3</b> Integrated and Mobile	Digital IVUS	Eagle Eye Platinum	NA
	Rotational IVUS	Eagle Eye Platinum ST	
	ChromaFlo IVUS	Refinity ST	
	Philips Remote Services		
<b>IntraSight 5</b> Integrated and Mobile	Digital IVUS	Eagle Eye Platinum	OmniWire
	Rotational IVUS	Eagle Eye Platinum ST	Verrata Plus
	ChromaFlo IVUS	Refinity ST	
	iFR FFR		
	Philips Remote Services		
<b>IntraSight 7*</b> Integrated only	Digital IVUS	Eagle Eye Platinum	OmniWire
	Rotational IVUS	Eagle Eye Platinum ST	Verrata Plus
	ChromaFlo IVUS	Refinity ST	
	iFR FFR		
	IVUS and iFR Co-registration IVUS, iFR and angio Tri-registration		
	Vessel Enhancement Quantitative Coronary Analysis Device Detection		
	Philips Remote Services		

\* Configuration via SyncVision

# Core

## Integrated and Mobile systems



		Core Integrated	Core Mobile
Power requirements	System input	100, 120 V, 220, 240 VAC, 50/60 Hz, 1000 VA	100, 120 V, 220, 240 VAC, 50/60 Hz, 1000 VA
	Workstation	100–240 V, 50/60 Hz, 825 VA	n/a
	Monitor	100V–240 V 50/60 Hz, 39 W	n/a
Dimensions	Core mobile cart	-	H=157.48 cm, W=55.88 cm, D=83.82 cm, ~95.25 kg
	Workstation	H=43.2 cm, W=25.4 cm, D=41.6 cm, ~22.6 kg	-
	Controller	H=12 cm, W=36 cm, D=24 cm, ~3 kg	-
	Monitor	H=42 cm, W=43.2 cm, D=24 cm, ~15 kg	-
Processing and data storage	Processor	1 CPU processor 2.53 GHz, 8 core total, 1366 MHz BUS	1 CPU processor 2.53 GHz, 8 core total, 1366 MHz BUS
	Memory	8 GB RAM	8 GB RAM
	Hard drive capacity	1 TB 7200 RPM SATA	1TB 7200 RPM SATA
	Digital archiving capabilities	Local, DVD, DICOM Network	Local, DVD, DICOM Network
	Dicom services supported	DICOM worklist management, DICOM store	DICOM worklist management, DICOM store

### Ordering information

#### Core Integrated


797402

#### Core Mobile

797414


# SyncVision

## Making a difference with Live Image Guidance



### iFR Co-registration

Provides localization of physiology measurements on the angiogram during the decide and confirm stages of a PCI procedure.



### IVUS Co-registration

Provides localization of the IVUS image on the angiogram during the decide, guide and confirm stages of a PCI procedure. It can additionally be synchronized with physiology measurements, providing a Tri-Registration that together shows the localization of IVUS and physiology on the angiogram image.



### Angio+ Enhancement

Enhancement and quantification of angiographic images during the decide, guide and confirm stages of a PCI procedure.



# An innovative suite of PCI applications

At Philips, we are committed to delivering technology that makes a difference – to clinical workflows, medical procedures, and to patients' and clinicians' lives. Our portfolio of Live Image Guidance solutions integrates imaging technologies and real-time patient information to help you decide, guide and confirm the right therapy for your patient.

Live Image Guidance includes dedicated solutions that help you increase the efficiency of percutaneous coronary interventions (PCI) – in both the diagnostic and therapeutic phase. Philips solutions offer reliable guidance that helps you make informed decisions and enhance clinical workflow.

## A new perspective

One of our innovative PCI suite application is SyncVision, a comprehensive solution designed to address imaging and physiological challenges in the cath lab. SyncVision allows you to plan your PCI strategy providing guidance and mapping of intravascular ultrasound and physiological information into the angiogram. Moreover, it helps to streamline lesion assessment, simplify vessel sizing, and support precise therapy delivery.

## Key features and benefits

- Intuitive user interface
- Physiology Co-registration
- IVUS Co-registration
- Physiology-IVUS Tri-registration
- Angio+ enhanced angiography featuring Vessel Enhancement, Quantitative Coronary Analysis and Device Detection
- Seamless integration into your PCI workflow

### SyncVision

Power requirements	System input	100 V-120 V, 50/60 Hz. 220-240 V, 50/60 Hz, 600 VA
	Workstation	100-240 V, 50/60 Hz, 250 VA
	Monitor	100-240 V, 50/60 Hz, 93 VA
Dimensions	Workstation	H=41.91 cm, W=17.15 cm, D=54 cm, -15.87kg
	Monitor	H=38.1-48.26 cm (adjustable stand), W=40.13 cm, D=24.64 cm, -5.9 kg
	Joystick	H=3.81 cm, W=10.67 cm, D=7.62 cm, -0.9 kg
Processing and data storage	Processor	1 GPU P5000 1 CPU Intel E5-1600/E5-2600 Series Processor
	Memory	16 GB SD RAM
	Hard drive capacity	120 GB SSD SATA + 480GB SSD SATA
	Digital archiving capacity	n/a
Compatible disposables	Catheters	Eagle Eye Platinum Eagle Eye Platinum ST
	Pressure wires	OmniWire Verrata Plus

Ordering information

SyncVision

797406



# What is **IVUS**

## Intravascular ultrasound


- When electricity stimulates a crystal, sound waves are released
- The sound hits tissue and is reflected to a varying degree
- Crystals then receive those reflected signals and code them back into electrical impulses
- Those impulses can construct a sophisticated image using computing power

**Intima**

**Dead zone**  
Catheter

**Adventitia**  
White ring

**Media**  
Black ring




# Two ways to acquire IVUS signals

## Digital IVUS

### Eagle Eye Plainum and Eagle Eye Plainum ST

A series of fixed 64 ultrasound devices that scan sections of the vessel in turn, whilst computers to build the image, 20MHz greyscale.




## Rotational IVUS

### Refinity ST and Revolution

A single IVUS catheter that spins within a sheath and uses computers to build the image in real-time from the sections that come in bit by bit, 45MHz greyscale.

## MHz – megahertz

MHz refers to the frequency of the soundwave. A higher frequency will increase the resolution but reduce the depth penetration.








## Digital IVUS

Eagle Eye Platinum and  
Eagle Eye Platinum ST





## Rotational IVUS

Refinity, Revolution and SpinVision



## ChromaFlow imaging

Provides clear confirmation of stent apposition  
and lumen size



Quick and easy to use measurement, annotation  
and labelling tools




# Eagle Eye Platinum


# Eagle Eye Platinum ST

## Phased-array, digital IVUS catheter

- Fast, plug-and-play simplicity
- Quick, convenient length estimation
- ChromaFlo imaging
- Unique SyncVision IVUS and iFR co-registration compatibility
- Standard and short-tip variations



Eagle Eye Platinum




Eagle Eye Platinum ST



# Refinity

## Our next generation rotational IVUS catheter

- Low profile, 5Fr guide compatible
- Suitable for radial access
- GlyDx hydrophilic coating
- Soft and flexible tip with short distance transducer
- Smooth proximal to mid-shaft transition




### 45 MHz, high resolution imaging




Refinity ST insertion and recovery forces are significantly lower ( $p < 0.001$ ) than OptiCross\*

\*Data on file at Philips. Bench testing was conducted with three rotational IVUS catheters: Refinity ST, Revolution, OptiCross through a tortuous model







# SpinVision

## Rotational IVUS pullback device

### Accurate

- Automate pullback system - facilitate accurate vessel and lesion measurements
- Maximum pullback accuracy - facilitates clinical research and trials

### Reliable

- Robust and single board design - provides consistent communication with Core system

### Cost effective

- Reusable sled design - for efficient and cost effective use



### SpinVision pullback device

Pullback length	15 cm
Completely reusable pullback device	Yes
Pullback speed	0.5 mm/s 1 mm/s
LED position display	Yes


Ordering information

SpinVision




# Coronary IVUS catheters

## Features and specifications



**Eagle Eye Platinum**



**Eagle Eye Platinum ST**

<b>Frequency/type</b>	20 MHz digital 64 element	20 MHz digital 64 element
<b>Vessel</b>	Coronary, peripheral	Coronary, peripheral
<b>Maximum imaging diameter</b>	20 mm	20 mm
<b>Diameter at transducer</b>	3.5F	3.5F
<b>Tip entry profile</b>	0.019"	0.019"
<b>Tip-to-transducer length</b>	10 mm	2.5 mm
<b>Radiopaque markers</b>	Scanner + 3 marks 10 mm apart	Scanner + 3 marks 10 mm apart
<b>Working length</b>	150 cm	150 cm
<b>Wire lumen length</b>	24 cm	24 cm
<b>Maximum guide wire</b>	0.024"	0.014"
<b>Minimum guide catheter</b>	5F (ID $\geq$ .056")	5F (ID $\geq$ .056")
<b>Minimum sheath size</b>	5F	5F
<b>Delivery platform</b>	Rapid exchange	Rapid exchange
<b>Proprietary imaging</b>	VH IVUS, ChromaFlo	VH IVUS, ChromaFlo
<b>Compatible with SyncVision</b>	Yes	Yes

Eagle Eye Platinum and Eagle Eye Platinum ST are also indicated for peripheral applications. Safety and effectiveness of VH IVUS for use in the characterization of vascular lesions and tissue types has not been established.


### Ordering information

#### Eagle Eye Platinum


85900P

#### Eagle Eye Platinum ST

85900PST



**Revolution**



**Refinity ST**

45 MHz rotational single element	45 MHz rotational single element
Coronary	Coronary
14 mm	14 mm
3.2F	3.0F
0.022"	0.025"
29 mm	20.5 mm
One marker + transducer	One marker + transducer
135 cm	135 cm
2.3 cm	1.55 cm
0.014"	0.014"
6F (ID $\geq$ .064")	5F (ID $\geq$ .056")
N/A	N/A
Rapid exchange	Rapid exchange
N/A	N/A
No	No

## Ordering information

Revolution	Refinity ST
89000	89900




# What is **physiology**

## Two ways to predict ischemia

1

### **iFR** - Instantaneous Free wave Ratio

An instantaneous pressure gradient, across a stenosis during the wave-free period, when resistance is constant and minimised in the cardiac cycle.<sup>1</sup>




$$iFR = \frac{P_{D_{wfp}}}{P_{A_{wfp}}}$$

2

### **FFR** - Fractional Flow Reserve

The ratio of distal mean coronary pressure to mean aortic pressure in the stenotic vessel during maximum hyperaemia.

Represents the very fraction of blood flow that still has been preserved despite the stenosis.<sup>1</sup>



$$FFR = \frac{P_{D_{hyp}}}{P_A}$$



iFR-guided strategy significantly reduces:



Time



Patient discomfort



Procedure cost per patient


1. Davies JE, et al., Use of the Instantaneous Wave-free Ratio or Fractional Flow Reserve in PCI. N Engl J Med. 2017 May 11;376(19):1824-1834.




## iFR vs. FFR:

same wire, same system, fewer steps

iFR workflow



FFR workflow



Single dichotomous cut-point back by data<sup>2,3</sup>



Both DEFINE FLAIR and iFR Swedeheart were designed with the dichotomous cut-point wof iFR in the iFR arm. With comparable MACE rates to FFR, these results mean the 0.89 cut-point for iFR is proven and backed by more than 4500 patients of outcome data.<sup>1</sup>

1. Davies JE, et al., Use of the Instantaneous Wave-free Ratio or Fractional Flow Reserve in PCI. N Engl J Med. 2017 May 11;376(19):1824-1834.



2. Gotberg M, et al., iFR-SWEDEHEART Investigators.. Instantaneous Wave-free Ratio versus Fractional Flow Reserve to Guide PCI. N Engl J Med. 2017 May

3. An iFR cut-point of 0.89 m tches best with an FFR ischemic cut-point of 0.80 with a specificity of 87.8% and sensitivity of 73.0% (From ADVISE II, and iFR Operator's Manual 505-0101.23).




# OmniWire


World's first solid core pressure guide wire<sup>1</sup>.




Traditional pressure wire



OmniWire with solid core



Workhorse guide wire



1. Data on file D000410086\_A, D000485394\_A

## Ordering information

OmniWire 185cm straight


89185

OmniWire 185cm J-tip

89185J

# Verrata Plus

Pressure wire



## The Quick Connect process





1. Angled Approach



2. Drag



3. Drop



## Ordering information

Verrata Plus 185cm straight    Verrata Plus 185cm J-tip

10185 P

10185J P







# AngioSculpt RX PTCA

Your prescription for complex lesions

The AngioSculpt RX PTCA scoring balloon catheter delivers the proven advantages of the AngioSculpt system combined with a new, tapered tip design that enhances deliverability.

The result is an essential tool in the treatment of a wide range of coronary lesions, including in-stent restenosis (ISR) and type C lesions.

1

Large working range (2 atm up to 20 atm) allows physicians to tailor device to vessel size.

2


Nitinol-enhanced balloon deflation for excellent rewrap and recross capabilities.

3


Electropolished, helical scoring element safely scores lesion circumferentially.

4


Rectangular edges provide a predictable dilation resulting in low dissection rates and minimal device slippage.








Edges lock in



~15-25x scoring force



~1x force post-scoring



## Precision

### Proper placement

- Rectangular scoring edges lock the device in place. Minimal device slippage or 'watermelon seeding,' even in ISR.<sup>1</sup>

## Power

### Enhanced mechanical advantage

- The leading edges are designed to drive outward expansion with up to 15-25 times the force of conventional balloons.<sup>2</sup>
- Helical nitinol scoring element creates a large luminal expansion for stent implantation.<sup>3</sup>

## Safety

### Predictable results

- Post-scoring, outward forces are designed to be equivalent to that of a conventional balloon.
- Low dissection rate of 13.6% (majority were non-flow limiting).<sup>1</sup>

- Mooney M, Teirstein P, Moses J, et al. Final results from the U.S. multi-center trial of the Angiosculpt Scoring Balloon Catheter for the treatment of complex coronary artery lesions. Am J Cardio. 2006;98(suppl B):121M
- Angiosculpt Test Plan ST-1197 (2008), on file at AngioScore inc.
- Costa JR, Mintz GS, Carlier SG, et al. Nonrandomized comparison of coronary stenting under intravascular ultrasound guidance of direct stenting without predilation versus conventional predilation with a semi-compliant balloon versus predilation with a new scoring balloon. Am J Cardiol. 2007;100:812-817.

## Ordering information

Angiosculpt 2.0x6mm	Angiosculpt 2.0x10mm	Angiosculpt 2.0x15mm
2200-2006	2200-2010	2200-2015
Angiosculpt 2.5x6mm	Angiosculpt 2.5x10mm	Angiosculpt 2.5x15mm
2200-2506	2200-2510	2200-2515



	Balloon diameter (mm)	Balloon length (mm)	Catheter length	Guide wire compatibility	Guide catheter compatibility
Angiosculpt 2.0x6mm	2.0	6	137	0.014"	6F
Angiosculpt 2.0x10mm	2.0	10	137	0.014"	6F
Angiosculpt 2.0x15mm	2.0	15	137	0.014"	6F
Angiosculpt 2.5x6mm	2.5	6	137	0.014"	6F
Angiosculpt 2.5x10mm	2.5	10	137	0.014"	6F
Angiosculpt 2.5x15mm	2.5	15	137	0.014"	6F
Angiosculpt 3.0x6mm	3.0	6	137	0.014"	6F
Angiosculpt 3.0x10mm	3.0	10	137	0.014"	6F
Angiosculpt 3.0x15mm	3.0	15	137	0.014"	6F
Angiosculpt 3.5x6mm	3.5	6	137	0.014"	6F
Angiosculpt 3.5x10mm	3.5	10	137	0.014"	6F
Angiosculpt 3.5x15mm	3.5	15	137	0.014"	6F


## Ordering information

Angiosculpt 3.0x6mm	Angiosculpt 3.0x10mm	Angiosculpt 3.0x15mm
2200-3006	2200-3010	2200-3015
Angiosculpt 3.5x6mm	Angiosculpt 3.5x10mm	Angiosculpt 3.5x15mm
2200-3506	2200-3510	2200-3515



# Elca

## Coronary laser atherectomy catheter



### Proven technology

- Treating patients for more than 20 years
- Optimally spaced fibers for improved performance
- Adjustable laser energy settings to satisfy many clinical needs
- Automatic shut-off feature for advanced patient safety

### Advanced performance

- Saline infusion improves safety outcomes<sup>1</sup>
- Slow advancement increases luminal gain<sup>2</sup>
- Two-thirds vessel sizing rule for predictable outcome

### Broad range of indications

- Total occlusion traversable by a guidewire
- Occluded SVGs
- Ostial lesions
- Moderately calcified stenoses
- Long lesions (>20mm)
- Lesions which previously failed PTCA
- Restenosis in 316L stainless steel stents prior to brachytherapy




	Elca 0.9 mm X-80	Elca 0.9 mm X-80 OTW*	Elca 1.4 mm	Elca 1.7 mm	Elca 2.0 mm
Guidewire compatibility (in)	0.014	0.014	0.014	0.014	0.014
Guide catheter compatibility (F)	6	6	6/7	7	8
Minimum vessel diameter (mm)	2.0	2.0	2.2	2.5	3.0
Max tip outer diameter (in)	0.038	0.038	0.057	0.069	0.080
Max shaft outer diameter (in)	0.049	0.049	0.062	0.072	0.084
Working length (cm)	130	130	130	130	130
Fluence (mJ/mm <sup>2</sup> )	30-80	30-80	30-60	30-60	30-60
Repetition rate (Hz)	25-80	25-80	25/40	25-40	25-40
Laser on/off time (sec)	10/5	10/5	5/10	5/10	5/10

### Saline infusion recommendations for coronary interventions


- Always perform 10-20cc bolus infusion of saline through the guide catheter after contrast injections.
- During lasting, infuse saline through the guide catheter at a rate of 2-3cc/second.

### Ordering information

Elca 0.9 mm X-80	Elca 1.4 mm
110-004	110-009
Elca 1.7 mm	Elca 2.0 mm
117-016	120-009



Larger luminal glow with slow advancement<sup>2</sup>



1. Tcheng J.E. et al. (1995). Development of a New Technique for Reducing Pressure Pulse Generation During 308-nm Excimer Laser Coronary Angioplasty. Catheterization and cardiovascular Diagnosis. 34, 15-22.
2. Topaz, On, et al, 2001. Optimal Spaced Excimer Laser Coronary Catheters Performance Analysis, Journal of Clinical Laser Medicine and Surgery, Vol 19, issue 1, 9-14.
3. Data on file at Philips..





# Philips Laser System

A broad range of applications to treat more complex conditions

The Philips Laser System is built on over 20 years of proven technology, safely and reliably photoablates a fully spectrum of morphologies.

The Philips laser catheters are indicated in more vessel types than other atherectomy devices and is the only laser system available for electrophysiology.



## Simple: turn key ready in any room

- Ready in less than 30 seconds
- Touchscreen guided workflows
- Minimum training required
- 360° maneuverability and small lab footprint
- Easy positioning within any lab

## Versatile: more vessels, more indications

- Seven indications for coronary vascular
- Only\* device with level 1 evidence for In-Stent Restenosis (ISR)
- First technology proven effective for ISR in peripheral vasculature<sup>1</sup>
- Only laser technology for lead removal
- Variable fluence and rate settings to match patients' needs

## Proven: excimer technology

- Cool, ultraviolet laser with over 20 years of clinical experience
- More than 600,000 procedures performed
- Only atherectomy device with level 1 clinical evidence showing superiority in safety and efficacy of ELA + PTA versus PTA alone for treating femoropopliteal ISR<sup>1</sup>

\* Data of 2021, subject to change.

1. Dippel et al. Randomized Controlled Study of Excimer Laser Atherectomy for Treatment of Femoropopliteal In-stent Restenosis: Initial ISR Results (2015). JACC 8(1): 92-101.\*



#### Philips Laser System


Power requirements	100V-240V, 16 amp, single phase power
Wavelength	308 nm
Class	Class IV laser system
Length	52 in/132 cm
Height	42 in/107 cm
Width	19 in/48 cm
Weight	450 lb/204 kg

# CVX-300

## Excimer Laser System

A broad range of applications to treat more complex conditions

The CVX-300 Excimer laser system's clinical versatility, high clinical success, low adverse events and well-established reimbursement<sup>1-6</sup> helps you to safely treat more complex conditions in vascular intervention and lead management procedures.



### Coronary atherectomy

- In-stent Restenosis (ISR)\*
- Moderately calcified lesions
- Ostial lesions
- Lesions that previously failed PTCA
- CTO traversable by a guidewire
- Occluded SVG
- Long lesions (>20mm)

### Peripheral atherectomy

- Only\*\* device indicated for fem-pop ISR
- No contraindications

### Lead extraction

- Removal of chronically implanted pacing and defibrillator leads
- 30 HRS indications<sup>4-6</sup>
- Has an estimated 5% annual incidence rate<sup>4-6</sup>

\* ISR is limited to BMS (316L SS) and prior to administering brachytherapy.

\*\* Data of 2021, subject to change.

1. Dippel et al. Randomized Controlled Study of Excimer Laser Atherectomy for Treatment of Femoropopliteal In-stent Restenosis: Initial EXCITE ISR Results (2015). JACC 8(1): 92-101.
2. Doshi et al. Comparison of Excimer Laser Atherectomy versus Orbital plus Rotational Atherectomy for Revascularization. SCAI 2017.
3. Gajananana et al. Global Revascularization & Evaluation of Excimer Laser in the Coronaries. SCAI 2017.
4. Wilkoff, B.L., Love, C.J., Byrd, C.L., Bongiorno, M.G., Carrillo, R.G., Crossley, G.H., et al. (2009). Transvenous Lead Extraction: Heart Rhythm Society Expert Consensus on Facilities, Training, Indications, and Patient Management.
5. Philips data on file. D014953-10. June 2017.
6. Wazni, O et. al. Lead Extraction in the Contemporary Setting: The LEXiCon Study: A Multicenter, Observational Retrospective Study of Consecutive Laser Lead Extractions, J Am Coll Cardiol, 55:579-586.





### Philips Laser System

<b>Power requirements</b>	208-230 VAC single phase power
<b>Wavelength</b>	308 nm
<b>Class</b>	Class IV laser system
<b>Length</b>	49 in / 125 cm
<b>Height</b>	35in/89 cm plus 6.9 in/17.5 cm control panel
<b>Width</b>	24 in/61.3 cm
<b>Weight</b>	650 lb/295 kg

## IGT Devices coronary portfolio

# Ordering information

### IVUS and physiology systems

IntraSight	IntraSight 3	IntraSight 5	IntraSight 7
	797403 (series 3)	797403 (series 5)	797403/797406 (series 7)
IntraSight Mobile		IntraSight Mobile 3 797415 (series 3)	IntraSight Mobile 5 797415 (series 5)
Core		Core Integrated 797402	Core Mobile 797414
SyncVision			SyncVision 797406

### IVUS disposables

Eagle Eye Platinum	<b>Eagle Eye Platinum</b> 85900P
Eagle Eye Platinum ST	<b>Eagle Eye Platinum ST</b> 85900PST
Revolution	<b>Revolution</b> 89000
Refinity	<b>Refinity ST</b> 89900
SpinVision	<b>SpinVision</b>

Products subject to country availability. Please contact your local sales representative.

## Physiology

OmniWire	<b>OmniWire 185cm straight</b>	<b>OmniWire 185cm J-tip</b>
	89185	89185J
Verrata Plus	<b>Verrata Plus 185cm straight</b>	<b>Verrata Plus 185cm J-tip</b>
	10185 P	10185J P

## Therapy disposables

AngioSculpt RX PTCA	<b>Angiosculpt 2.0x6mm</b>	<b>Angiosculpt 2.0x10mm</b>	<b>Angiosculpt 2.0x15mm</b>
	2200-2006	2200-2010	2200-2015
	<b>Angiosculpt 2.5x6mm</b>	<b>Angiosculpt 2.5x10mm</b>	<b>Angiosculpt 2.5x15mm</b>
	2200-2506	2200-2510	2200-2515
	<b>Angiosculpt 3.0x6mm</b>	<b>Angiosculpt 3.0x10mm</b>	<b>Angiosculpt 3.0x15mm</b>
	2200-3006	2200-3010	2200-3015
	<b>Angiosculpt 3.5x6mm</b>	<b>Angiosculpt 3.5x10mm</b>	<b>Angiosculpt 3.5x15mm</b>
	2200-3506	2200-3510	2200-3515
	<b>Elca 0.9 mm X-80</b>	<b>Elca 1.4 mm</b>	
	110-004	110-009	
Elca	<b>Elca 1.7 mm</b>	<b>Elca 2.0 mm</b>	
	117-106	120-009	

## Laser systems

Philips Laser System	<b>Philips Laser System</b>
	LAS-100
CVX-300	CVX-300
	CVX-300-P

Products subject to country availability. Please contact your local sales representative.





© 2021 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners. Approved for external distribution. D2021042101

Products subject to country availability. Please check with your local sales representative.

**Philips**  
3721 Valley Centre Drive, Suite 500  
San Diego, CA 92130 USA

**Philips**  
Excelsiorlaan 41  
1930 Zaventem, Belgium

[www.philips.com/IGTDevices](http://www.philips.com/IGTDevices)