PHILIPS

Correct Storage of Chemical/Thermal Recording Paper

Application Note

Images on standard chemical/thermal paper decompose naturally over time. With proper storage and handling, recordings on those papers have been shown to remain legible for approximately 5 years. **Proper storage and handling before and after recording includes:**

Store the paper under cool and dry conditions. Temperatures must not exceed 40°C (104°F) and relative humidity must be below 80%.

Store identical chemical/thermal papers in manila folders or sheet protectors made of polyester, polyimide, polypropylene, or acetate.

Storing the paper as indicated above will minimize image fading and background development (darkening). However, storing the paper as indicated above does not protect against image fading or background development from the sources below.

To avoid image fading or background development, the paper **must not be exposed to, enclosed with, or come in contact with** the following, either before or after recording:

- Solvent based adhesives such as those found in mounting forms, pressure-sensitive tapes or labels, or Scotch or mending tape. The use of adhesives made of starch or water based formulations are acceptable.
- Plastics (containing plasticizers) such as PVC or other vinyl chlorides that are typically found in sheet protectors, separators and envelopes. Other plasticizers also include: polyethylene glycol, diethyl phthalate, diethyl adipate.
- Glossy and some non-glossy papers such as non-chemical/thermal recording papers or product literature that contain tributyl phosphate, dibutyl phthalate or other organic solvents. This includes PageWriter 4700 series paper.
- Liquid or vaporous forms of solvents such as alcohols, ketones, esters, ethers, etc., some of which are found in felt tip pens.
- Carbonless (NCR) and carbon based forms.
- Petroleum solvents such as toluene, benzene and gasoline.
- Bright light or UV sources such as sunlight, fluorescent and related light sources.
- Dissimilar chemical/thermal papers or similar paper from different manufacturers.
- Certain other chemicals including: castor oil and chemicals found in hand and face creams, citric acid as in fruit juices and ammonia.
- EKG, Ultrasound and related pastes, creams and gels (these contain some of the above or related chemicals).

If original recordings are stored, Philips Medical Systems recommends that records be checked annually to determine their integrity. However, where long term storage is desired, the user should consider photocopying or microfilming, or electronic or optical storage or a fade resistant recording paper.

Philips Medical Systems PTP™ brand thermal papers offer improved archivability. See Philips Application Note 5980-1342E for more information concerning Philips's PTP brand recording papers.

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**Philips Medical Systems is part
of Royal Philips Electronics**

Asia

Tel: +852 2821 5888

Europe, Middle East, Africa

Tel: +31 40 27 63005

Latin America

Tel: +55 11 2125 0764

North America

Tel: +1 425 487 7000

Would you like to know more about our
imaginative products? Please do not hesitate to
contact us. We would be glad to hear from you.

On the web

www.medical.philips.com

Via e-mail

medical@philips.com

By fax

+31 40 27 64 887

By postal service

Philips Medical Systems

Global Information Center

I.B.R.S./C.C.R.I. Numéro 11088

5600 VC Eindhoven

Pays-Bas/The Netherlands

(no stamp required)

© Koninklijke Philips Electronics N.V. 2003
All rights are reserved. Reproduction in whole or in
part is prohibited without the prior written consent
of the copyright holder.

Philips Medical Systems Nederland B.V. reserves the
right to make changes in specifications and/or to
discontinue any product at any time without notice
or obligation and will not be liable for any conse-
quences resulting from the use of this publication.

© Copyright 2001 Philips Medical Systems
Printed in U.S.A. 08/2004
AW-5980-1343-2 Rev. D
5980-1343E