PHILIPS

Esophageal/Rectal

Temperature Probe 21090A, M1837A

With 400 Series Thermistor

Application Note

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PHILIPS

ENGLISH

Esophageal/Rectal Temperature Probe 21090A, M1837A With 400 Series Thermistor

Application Note

Description

Hypothermia and hyperthermia are well recognized clinical conditions necessitating temperature monitoring.

Philips Medical Systems' economical and disposable esophageal/rectal temperature probe is a sensitive and accurate temperature transducer; to be used clinically where continuous temperature monitoring is required. The sensor is designed for insertion into esophagus or rectum and consists of an electrically insulated thermistor permanently secured within a PVC tube. The bullet tipped PVC tube providesatraumatic insertion. Leakage current of the temperature probe (thermistor and PVC tube) is well within clinically accepted non-hazardous levels. The temperature probe is disposable, single patient use only, packaged sterile and available in 9F and 12F size.

Range and Accuracy

The accuracy of the thermister is $\pm 0.1^{\circ}\text{C}$ from 25°C to 45°C .

Indications

Esophageal/Rectal temperature probe is intended to be used with Philips Temperature Monitoring Modules via extension cables having either 1/4 inch phone jack or two-pin connectors.

Directions for Use and Precautions

1. If a patient has to be intubated with an E.T. tube, perform intubation prior to placing the temperature sensor into esophagus.
2. Lubricate the temperature sensor prior to insertion, and place the probe in accordance with currently acceptable medical procedures.
3. Verify position of the sensor by acceptable medical procedures.
4. Align the sensor's connector with the monitor cable's connectors and push firmly to assure a full contact. Forced mating of the connectors without proper alignment may cause damage to the connectors and interruption in electrical continuity.

5. During the surgical procedures which employ electro-cautery, use currently acceptable procedures to minimize conditions of the thermistor and lead wires functioning as an alternate path for radio-frequency current to return to ground, causing localized tissue burns. Procedures which may minimize risk of electro-surgical burns are:

Keep both active and ground electrodes of the electro-cautery system in close proximity so that the sensor is outside of the radio-frequency current field.

Keep the temperature monitor with its associated cables separated from electrocautery systems.

Unusual, fast artificial variations in temperature readings may occur with concomitant applications of the electrocautery system.

6. Reuse may cause degradation of physical or electrical properties. Do not reuse on another patient due to risk of cross-infection.

Electromagnetic Compatibility

The electromagnetic compatibility (EMC) validation of the 21090A and M1837A included testing performed according to the international standard for EMC with medical devices. See the Manufacturer's Declaration for details.

Avoiding Electromagnetic Interference

If electromagnetic interference (EMI) is encountered there are a number of things that can be done to mitigate the problem.

- Eliminate the source. Possible sources of EMI can be turned off or moved away to reduce their strength.
- Attenuate the coupling. If coupling is through accessory leads, the interference may be reduced by moving and/or rearranging the leads.
- Add external attenuators. If EMI becomes an unusually difficult problem, external devices may be of help. A Philips Customer Engineer can be of help in determining the need for external devices.

Adverse Reactions

Adverse reactions reported during applications of these sensors include: accidental tracheal or bronchial intubation accompanied with airway obstruction, esophageal or rectal abrasion and/or perforation, pharyngeal abrasion, and tissue burns due to aberrant electro-cautery RF current pathways.

Disposal

Follow local regulations for disposal.

Contraindications

The use of the esophageal/rectal sensor may be contraindicated in neonates and small infants undergoing tracheostomy or internal jugular vein catheterization and laser surgical procedures.

Sterile-Disposable

CE 0123

 Do not reuse	 Do not use if seal is broken
 Caution see documentation	STERILE  Sterilized via irradiation

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Copyright © 2012
Böblingen GmbH