Bridgesheath

Splittable Hemostatic Introducer System

en	Instructions for Use	2
cz	Návod k použití	4
da	Brugervejledning	6
de	Gebrauchsanweisung	8
el	Οδηγίες Χρήσης	10
es	Instrucciones de uso	12
et	Kasutusjuhend	14
fi	Käyttöohjeet	16
fr	Instructions d'utilisation	18
hu	Használati utasítás	20
it	Istruzioni per l'uso	22
jp	取扱説明書	24
ko	사용 지침	26
lt	Naudojimo instrukcijos	28
nl	Gebruiksaanwijzing	30
no	Bruksinstruksjoner	32
po	Instrukcja użycia	34
pt	InSTRUÇÕES de Utilização	36
sc	使用说明	38
sk	Pokyny na používanie	40
sv	Bruksanvisning	42
tu	Kullanım Talimatları	44

BridgeSheath

Splittable Hemostatic Introducer System

Instructions for Use

This device is intended for one time use only.
Read instructions prior to use.

CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

For U.S.-California Only.

Proposition 65, a State of California voter initiative, requires the following notice:

WARNING: This product and its packaging have been sterilized with ethylene oxide. This packaging may expose you to ethylene oxide, a chemical known to the state of California to cause cancer or birth defects or other reproductive harm.

Indications

For the introduction of various types of pacing leads and catheters.

Warnings

- This product is sensitive to light. Do not use if stored outside the protective outer carton.
- Store in a cool, dark, and dry place.
- Infusion through the sideport can be done only after all air is removed from the unit.

Precautions

- Do not alter this device in any way.
- **Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
- Aspiration and saline flushing of the sheath and dilator should be performed to help minimize the potential for air embolism and clot formation.
- Indwelling introducer sheaths should be internally supported by a catheter, pacing lead, or dilator.
- Dilators, catheters, and pacing leads should be removed slowly from the sheath. Rapid removal may damage the valve resulting in flow of blood or air through the valve.
- Never advance or withdraw the guidewire or sheath when resistance is met. Determine the cause by fluoroscopy and take remedial action.
- When injecting or aspirating through the sheath, use the sideport only.

Adverse Events

Adverse events may include, but are not limited to the following:

- | | |
|-----------------|----------------------------|
| • Air embolus | • Pneumothorax |
| • Blood loss | • Hemotorax |
| • Vessel damage | • Pacing lead displacement |
| • Infection | |

Use Sterile Technique

- A suggested procedure:
1. Peel open package and place contents in sterile field.
 2. Prep skin and drape in area of anticipated venipuncture.
 3. Flush the needle, sheath, and dilator prior to use.
 4. Insert the dilator into the sheath until the dilator cap attaches to the sheath hub.
 5. Insert the needle into vessel. Verify the needle position by observing venous blood return.
 6. Aspirate the needle using the syringe.
 7. Remove the syringe and insert soft tip of the guidewire through the needle into the vessel. Advance the guidewire to required depth. Leave an appropriate amount of the guidewire exposed. At no time should the guidewire be advanced or withdrawn when resistance is met. Determine the cause of resistance before proceeding. Fluoroscopic verification of the guidewire's entrance into the superior vena cava and right atrium is suggested.
 8. Hold the guidewire in place and remove needle. Do not withdraw the guidewire back into the needle as this may result in separation of the guidewire.
 9. Thread the dilator/sheath assembly over the guidewire.
 10. Advance the dilator/sheath assembly with a twisting motion over the guidewire and into the vessel. Fluoroscopic observation may be advisable. Attaching a clamp or hemostat to the proximal end of the guidewire will prevent inadvertently advancing the guidewire entirely into the patient.
 11. Once assembly is fully introduced into the venous system, remove the dilator cap from the sheath by rocking the dilator cap off the sheath hub.
 12. Slowly retract the guidewire and dilator, leaving the sheath in position. The sheath will reduce the loss of blood and the inadvertent aspiration of air.
 13. Aspirate all air from the sheath by using a syringe connected to the sideport.
 14. Flush the introducer with saline through the sideport. If the introducer is to remain in place during lead positioning and testing, periodic flushing is advised.
 15. Introduce pacemaker lead or catheter through the sheath and advance it into position.
 16. Flush sheath with saline immediately before peeling sheath away in order to minimize backbleeding.
 17. Split sheath by sharply snapping the sheath hub tabs and peeling sheath tube apart while withdrawing it from the vessel.

TVI Instructions for Use

The Transvalvular Insertion Tool (TVI) is an optional device for use with the BridgeSheath systems

This device is intended for one time use only.
Read instructions prior to use.

CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

Indications

For the introduction of various types of pacing or defibrillator leads and catheters.

Warnings

Improper use of the TVI can cause air embolism and back bleeding.

Precautions

- Do not alter this device in any way.
- **Single Use Devices:** This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
- When using the TVI tool, lead size may not exceed the size of the TVI tool.
- When using the TVI, always keep the exposed proximal end covered to prevent air embolism and back bleeding.
- When retaining a guidewire, the TVI tool must be placed over the guidewire and not along the side of it.

Use Sterile Technique A suggested procedure:

1. A TVI is to be used at the physician's discretion to open the sheath valve for ease of lead placement.
2. CAUTION! When the TVI is inserted into the BridgeSheath valve housing, all hemostasis is lost and the risk of air embolization and back bleeding exists. Always keep the proximal exposed end of the TVI covered with your thumb while the TVI is in use.
3. To use the TVI, insert the distal end of the TVI into the valve housing by gently pushing the TVI into the sheath.
4. Hold your thumb over the proximal exposed opening of the TVI to prevent air embolism or back bleeding.
5. Advance the pacing lead or catheter through the TVI and into the sheath.
6. As soon as the pacing lead or catheter is resting inside the sheath, pull the TVI back out of the sheath valve housing.
7. The TVI may then either be peeled away or it may be temporarily left resting on the shaft of the pacing lead or catheter.
8. Once the TVI has been withdrawn from the sheath, aspirate the sheath through the sideport, until any air which may have entered the sheath during the procedure is removed. Flush with Heparinized saline.

Graphic Symbols for Medical Device Labeling

Single Use Only /
Do Not Reuse

Do not use if package is
opened or damaged

Sterilized Using
Ethylene Oxide

Store in a cool, dark and
dry place

Caution, read instructions
prior to use

Size

Contents (Numeral repre-sents quantity of units inside)

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Uzunluk

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Non pirojenik

This page has been intentionally left blank

This page has been intentionally left blank

Distributed by:

Spectranetics Corporation

9965 Federal Drive, Colorado Springs, CO 80921 USA

Tel: 1-800-231-0978 • Fax: 719-447-2022

www.spectranetics.com

Spectranetics International B.V.

Plesmanstraat 6 • 3833 LA Leusden • The Netherlands

Tel: +31 33 43 47 050 • Fax: +31 33 43 47 051

Merit Medical Systems, Inc.

Bridge™

Prep Kit

DIAGNOSTIC GUIDE WIRES

Description: Bridge guide wires are fabricated from high quality stainless steel utilizing a sophisticated construction process and are available with or without a PTFE coating. Guide wires are supplied sterile, non-pyrogenic, and are intended for single use only.

Bridge guide wires are packaged in a plastic hoop, which is fitted with a luer hub. This packaging is provided to facilitate compliance with the manufacturer recommended guidelines that the wire be flushed with saline or heparinized saline prior to use. (See directions for use -

Indications: Bridge guide wires are used to facilitate the placement of devices during diagnostic and interventional procedures.


Contraindications: Diagnostic guide wires, are contraindicated for use in the coronary and cerebral vasculature.

Precautions: Angiography should be undertaken only by an experienced angiographer. For One Time Use Only. Guide wires will collect blood and other foreign material in their lumens; neither autoclaving nor ultrasonic cleaning will completely remove foreign material, therefore guide wires are recommended for one time use.

Inspect all guide wires prior to use. Do not use any unit if the package is open or damaged.

Employ an aseptic technique during removal from the package and during use.

All guide wires are secured in the hoop dispenser by the locking J-tip straightener. To avoid damaging the guide wire during removal from the flush hoop, grasp the J-tip straightener near the base and slide it forward approximately 5mm or until the J-tip straightener is no longer attached to the flush hoop adapter. Holding both the guide wire and J-tip straightener, continue to dispense guide wire from the hoop.

Do not use excessive force to advance the moveable core while the guide wire is in a vessel. Advancement with excessive force may cause coil penetration and vessel damage.

Avoid withdrawing PTFE coated guide wires back through a metal needle. The sharp edge of the needle may scrape the coating. It is suggested that a catheter or PTFE vessel dilator replace the needle as soon as the guide wire has reached the appropriate position.


During advancement of the catheter and guide wire within the aorta, it is recommended that the guide wire be removed at the appropriate level of the aorta.

Care should be taken when manipulating a catheter during placement and withdrawal to prevent possible intravascular tissue damage. If resistance is felt when removing a guide wire from a catheter, the guide wire and catheter should be removed as a unit to prevent potential damage to the vessel wall.

A guide wire is a delicate instrument and remains the most failible instrument used in a percutaneous procedure. Any time a guide wire is used there is a possibility of thrombus formation/emboli, vessel wall damage, and plaque dislodgement, which could result in myocardial infarction, cardiac arrhythmia, stroke or death. The physician should be familiar with the use of angiography products and the literature concerning the complications of angiography.

Flush Hoop Instruction Guide:

1. Attach flush filled syringe to flush port luer
2. Rotate syringe clockwise (as pictured)
3. Inject saline into hoop
4. Detach syringe from hoop
5. Dispense guide wire

Note: In order to reduce the potential of clot formation, it is recommended that the guide wire be flushed with saline or heparinized saline prior to use. Attach a filled syringe to the luer hub located at the end of the plastic hoop, flush several times. After flushing, remove guide wire from hoop and use as described above.

Directions for use: The following schematic shows a typical procedure for percutaneous entry utilizing the Seldinger technique. Variations in individual patient anatomy may preclude the utilization of this technique.

1. Vessel puncture with a two part needle

2. When using a two part needle remove stylet leaving cannula in place, insert flexible (distal) end of guide wire through cannula and into vessel


3. "J" guide wires are shipped with "J" straightener to aid in the insertion of the wire into the puncture needle. Advance the straightener until 2-3 mm of the tip extends from the tip. Insert wire into hub and through needle. Remove "J" straightener proximally and discard.

4. Remove needle cannula leaving the guide wire within the lumen of the vessel.

5. Pass dilator or catheter over the quide wire directly into the vessel.

6. Carefully remove the guide wire leaving the catheter in place.

Note: Device is sterilized as stated on the package label.

Distributed by:

Spectranetics Corporation
9965 Federal Drive, Colorado Springs, CO 80921 USA
Tel: 1-800-231-0978 • Fax: 719-447-2022 •
www.spectranetics.com

Spectranetics International B.V.
Plesmanstraat 6 • 3833 LA Leusden • The Netherlands
Tel: +31 33 43 47 050 • Fax: +31 33 43 47 051

0086

Manufacturer: Merit Medical Systems, Inc.
South Jordan, Utah 84095 U.S.A.
Authorized Representative:
Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland

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REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

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