QuickCatTM

Extraction Catheter

Instructions for Use

 ***Spectranetics***[®]

Instructions for Use - Sections by Language

Pg	Language		Instructions For Use
3	English	English	Instructions for Use
6	Bulgarian	български език	Указания за употреба
9	Croatian	Hrvatski	Upute za uporabu
12	Czech	Česky	Návod k Použití
15	Danish	Dansk	Brugsanvisning
18	Dutch	Nederlands	Gebbruiksaanwijzing
21	Estonian	Eesti	Kasutusjuhend
24	Finnish	Suomi	Käyttöohje
27	French	Français	Mode d'emploi
30	German	Deutsch	Gebrauchsanweisung
33	Greek	Ελληνικά	Οδηγίες χρήσης
36	Hungarian	Magyar	Használati utasítások
39	Icelandic	Íslenska	Leiðbeiningar um notkun
42	Italian	Italiano	Istruzioni per l'uso
45	Latvian	Latviešu valoda	Lietošanas instrukcija
48	Lithuanian	lietuvių kalba	Naudojimo instrukcija
51	Norwegian	Norsk bokmål	Bruksinstruksjoner
54	Polish	Polski	Instrukcja użytkowania
57	Portuguese	Português	Instruções de Utilização
60	Romanian	Română	Instrucțiuni de utilizare
63	Russian	Русский язык	Инструкции по применению
66	Serbian	Српски језик	Uputstvo za upotrebu
69	Slovak	slovenčina	Návod na použitie
72	Slovenian	Slovenščina	Navodila za Uporabo
75	Spanish	Español	Instrucciones de Uso
78	Swedish	Svenska	Bruksanvisning
81	Turkish	Türkçe	Kullanım Talimatları

Table of Contents

1.	DEVICE DESCRIPTION	3
2.	INDICATIONS AND INTENDED USE	3
3.	CONTRAINDICATIONS	3
4.	WARNINGS	3
5.	PRECAUTIONS	3
6.	POTENTIAL ADVERSE EVENTS	4
7.	HOW SUPPLIED	4
8.	COMPATIBILITY	4
9.	DIRECTIONS FOR USE	4
10.	MANUFACTURER'S LIMITED WARRANTY	5
11.	NONSTANDARD SYMBOLS	5

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to observe warnings and precautions may result in complications. Any recommendations within these instructions are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.

1. DEVICE DESCRIPTION

The QuickCat Extraction Catheter is a single use, disposable dual-lumen catheter designed for removing soft emboli and thrombi from vessels of the arterial system. The distal tip of the catheter is pliable, tapered, and smooth, providing atraumatic passage into the arterial system. The device is supplied sterile and is intended for one time use only.

2. INDICATIONS AND INTENDED USE

The QuickCat Extraction Catheter is indicated for removal of fresh, soft emboli and thrombi from vessels in the arterial system.

Product is intended for single use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheath, angiographic catheters and guidewires may be employed.

3. CONTRAINDICATIONS

- Use in vessels with a diameter < 1.5 mm
- The venous system
- The removal of fibrous, adherent or calcified material (e.g. chronic clot, atherosclerotic plaque)

4. WARNINGS

- Do not use without a guidewire, as vessel injury may result.
- Do not attempt to advance or retract the catheter against resistance until the cause of resistance has been determined by fluoroscopy or other means. Manipulation of the catheter against resistance may result in kinking of the catheter and/or vessel damage.
- If excessive slack or a loop in the guidewire is observed between the guide catheter and the monorail segment of the QuickCat Extraction Catheter during the procedure, the guidewire may become kinked within the vessel during catheter advancement or retraction. Remove the slack or loop in the guidewire before advancing or retracting the QuickCat Catheter to avoid catheter and/or vessel damage.
- If flow into the syringe stops or is restricted, do NOT attempt to flush the extraction lumen of the QuickCat Extraction Catheter while the catheter is inside the patient. Serious injury or death may result.
- Do not use a bent, kinked or damaged catheter as this may lead to vessel injury and/or an inability to advance or withdraw the catheter.
- Do not use for delivery or infusion of diagnostic, embolic, or therapeutic materials into blood vessels.

5. PRECAUTIONS

- Prior to use, check that all of the connections between components are secure and that the system has been completely primed or vacuum may be compromised.
- Do not tighten the hemostatic valve excessively onto the catheter shaft, or catheter damage may result.
- Use caution when crossing or retracting the QuickCat Extraction Catheter across a freshly deployed drug-eluting stent.
- Do not re-sterilize, re-process, or re-use the device.
- DO NOT sterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing.
- Reuse of this single use device could lead to serious patient injury or death and voids manufacturer warranties.
- Do not replace system components with alternate components.

6. POTENTIAL ADVERSE EVENTS

- Access site bleeding/hematoma
- Anastomotic disruption
- Abrupt closure or total occlusion of treated graft or vessel
- Distal embolization of debris resulting in pulmonary compromise and/or limb ischemia
- Local or systemic infection
- Arterial spasm
- Arteriovenous fistula formation
- Drug reactions, adverse reaction to contrast medium
- Acute myocardial infarction
- Vessel dissection, perforation, rupture or injury
- Emergent surgery
- Death

7. HOW SUPPLIED**7.1 STERILIZATION**

Supplied sterilized by radiation within snap-top tray, sealed in peel-open package.

Intended for single use; do not re-sterilize, re-process or re-use.

7.2 INSPECTION PRIOR TO USE

Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether packaging has been compromised. Upon removal from the package, thoroughly inspect the product to ensure that no damage, bends or kinks have occurred during shipment.

8. COMPATIBILITY**DEVICE COMPONENTS**

Catheter:

The QuickCat Extraction Catheter is compatible with 6F guide catheters [I.D. \geq 0.068" (1.73 mm)] and 0.014" (0.36 mm) guidewires. The catheter has a working length of 145 cm and its distal portion is covered with a hydrophilic coating. There is a radiopaque marker approximately 1 mm from the tip.

Vacuum Assembly:

The vacuum assembly consists of a 7.0" (177.8 mm) extension tubing, one-way stopcock and a 30 mL vacuum syringe with adjustable locking plunger. A 40-micron pore filter basket is provided to assist in filtering of blood and thrombotic material for visual or laboratory analysis.

9. DIRECTIONS FOR USE**9.1 STORAGE CONDITIONS**

Store in a cool, dry location. Avoid extended exposure to light.

9.2 DEVICE PREPARATION

1. Using sterile technique, open the pouch and transfer the tray into the sterile field.
2. Remove the catheter from the protective coil and inspect for any bends or kinks.
3. Fill the 30 ml vacuum syringe with 5-10 ml of normal saline.
4. Connect the vacuum syringe to the stopcock, the stopcock to the extension tubing, and the extension tubing to the hub of the QuickCat Extraction Catheter. Ensure that all connections are secure.
5. Flush the system with saline to ensure complete priming.
6. Turn the stopcock to the closed position.
7. Remove the 30 mL vacuum syringe and empty any excess saline solution.
8. Re-attach the empty 30 mL vacuum syringe to the system assembly.

9.3 PROCEDURE

1. Using standard technique, cannulate the target vessel with a guidewire and a guiding catheter (I.D. \geq 0.068" [1.73mm]) with an attached hemostatic valve.
2. Load the QuickCat Extraction Catheter onto the guidewire.

WARNING: If any resistance is felt during manipulation, the cause of the resistance must be determined before the catheter is advanced or withdrawn. Manipulation of the catheter against resistance may result in catheter or vessel damage. If the catheter becomes kinked during use, carefully remove from the patient and continue procedure with a new QuickCat Extraction Catheter.

3. Under fluoroscopy, advance the QuickCat Extraction Catheter to the target site.
4. Tighten the hemostatic valve sufficiently to prevent backflow, but not so tight as to impede catheter movement.

5. With the stopcock in the closed position, pull back the plunger on the 30 mL syringe to the desired amount of extraction volume. Twist the plunger clockwise to lock the syringe in the desired vacuum position.
6. Confirm the correct catheter position via fluoroscopy.
7. Turn the stopcock to the open position to begin the extraction. Slowly advance the QuickCat Extraction Catheter into the target vessel. Blood will enter the syringe until the vacuum is depleted.
8. If blood does not enter syringe within 5 seconds, remove the catheter from the patient. Either flush the catheter outside of the patient or replace with a new catheter.

WARNING: If flow into the syringe stops or is restricted, do NOT attempt to flush the extraction lumen of the QuickCat Extraction Catheter while the catheter is inside the patient. Serious injury or death may result. Remove catheter from the patient and flush the extraction lumen before attempting to reuse or use a new QuickCat Extraction Catheter.

9. Once the extraction process is complete, turn the stopcock to the closed position and remove the catheter from the patient. Extracted blood and thrombus may be filtered for visual and/or laboratory analysis using the 40-micron pore filter basket provided.


9.4 POST-PROCEDURE CARE

After the procedure, hospital standards of care should be followed for removing the sheath and providing hemostasis to prevent bleeding at the vascular access site.

10. MANUFACTURER'S LIMITED WARRANTY

Manufacturer warrants that the QuickCat catheter is free from defects in material and workmanship when used by the stated "Use By" date and when package is unopened and undamaged immediately before use. Manufacturer's liability under this warranty is limited to replacement or refund of the purchase price of any defective QuickCat catheter. Manufacturer will not be liable for any incidental, special, or consequential damages resulting from use of the QuickCat catheter. Damage to the QuickCat catheter caused by misuse, alteration, improper storage or handling, or any other failure to follow these Instructions for Use will void this limited warranty. **THIS LIMITED WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** No person or entity, including any authorized representative or reseller of Manufacturer, has the authority to extend or expand this limited warranty and any purported attempt to do so will not be enforceable against Manufacturer.

11. NONSTANDARD SYMBOLS

Wire Compatibility			
Rapid Exchange(RX)		Vessel Diameter	
Working Length		Guide Catheter Compatibility	
Sheath Compatibility		Quantity	QTY
CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician.			Rx ONLY

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©2019 Spectranetics Corporation
Spectranetics Corporation
Spectranetics International B.V.