Support Catheter

# Instructions for Use

CE  
0086

## Instructions for Use - Sections by Language

### Table of Contents

Pg	Language		Instructions For Use
3	English	English	Instructions for Use
5	Bulgarian	български език	Инструкции за употреба
7	Croatian	Hrvatski	Upute za uporabu
9	Czech	Česky	Návod k Použití
11	Danish	Dansk	Brugsanvisning
13	Dutch	Nederlands	Gebruiksaanwijzing
15	Estonian	Eesti	Kasutusjuhend
17	Finnish	Suomi	Käyttöohjeet
19	French	Français	Mode d'emploi
21	German	Deutsch	Gebrauchsanweisung
23	Greek	Ελληνικά	Οδηγίες χρήσης
25	Hungarian	Magyar	Használati utasítások
27	Icelandic	Íslenska	Leiðbeiningar um notkun
29	Italian	Italiano	Istruzioni per l'uso
31	Latvian	Latviešu valoda	Lietošanas instrukcija
33	Lithuanian	lietuvių kalba	Naudojimo instrukcija
35	Norwegian	Norsk bokmål	Bruksinstruksjoner
37	Polish	Polski	Instrukcja użytkowania
39	Portuguese	Português	Instruções de Utilização
41	Romanian	Română	Instrucțiuni de utilizare
43	Russian	Русский язык	Инструкции по применению
45	Serbian	Srpski jezik	Uputstvo za upotrebu
47	Slovak	Slovenčina	Návod na použitie
49	Slovenian	Slovenščina	Navodila za Uporabo
51	Spanish	Español	Instrucciones de Uso
53	Swedish	Svenska	Brugsanvisning
55	Turkish	Türkçe	Kullanım Talimatları

**Table of Contents**

1. DESCRIPTION	3
2. INDICATIONS FOR USE	3
3. CONTRAINDICATIONS	3
4. WARNINGS	3
5. PRECAUTIONS	3
6. ADVERSE EFFECTS	3
7. HOW SUPPLIED	3
8. COMPATIBILITY	4
9. PROCEDURE SET-UP AND DIRECTIONS FOR USE	4
10. MANUFACTURER'S LIMITED WARRANTY	4
11. NON-STANDARD SYMBOLS	4

**CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician**

**1. DESCRIPTION**

The Spectranetics Quick-Cross Select Support Catheters are intravascular catheters. These catheters are available in a variety of lengths and tip configurations. All models have 3 radiopaque markers spaced equally along the distal shaft to aid in estimating geometry within the vascular system. The distal radiopaque marker is positioned within 3 mm of the distal catheter tip. A standard female luer is placed on the proximal end of each model. The catheter is coated with a lubricious, hydrophilic coating.

**2. INDICATIONS FOR USE**

Quick-Cross Select Support Catheters are intended to guide and support a guidewire during access of the coronary or peripheral vasculature, allow for wire exchanges and provide a conduit for the delivery of saline solutions or diagnostic contrast agents.

**3. CONTRAINDICATIONS**

None known.

**4. WARNINGS**

- Maximum Infusion Pressure: **300 psi** for 0.014" & 0.018" catheters and **500 psi** for 0.035" catheters.
- The catheter is designed and intended for one time use only. **Do not re-sterilize and/or reuse.**
- The catheter should only be used by physicians qualified to perform percutaneous vascular interventions.

**5. PRECAUTIONS**

- DO NOT re-sterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing. Reuse of this single use device could lead to serious patient injury or death and voids manufacturer warranties.
- Store in a cool, dry place. Protect from direct sunlight and high temperature (greater than 55°C or 131°F).
- Do not use if device or packaging is damaged.
- Use the catheter prior to the "Use By" date (Expiration Date) specified on the package.

**6. ADVERSE EFFECTS**

Vascular catheterization and/or vascular intervention may result in complications including but not limited to:

- Vessel dissection, perforation, rupture or total occlusion
- Infection
- Hematoma
- Unstable angina
- Embolism
- Hypo/hypertension
- Acute myocardial infarction
- Arrhythmia, including ventricular fibrillation
- Death

**7. HOW SUPPLIED**

The Spectranetics Quick-Cross Select Support Catheters are supplied **STERILE**. The devices are designated and designed for **SINGLE USE ONLY** and must not be re-sterilized and/or reused.

**7.1 Sterilization**

- The sterility of the product is guaranteed only if the package is unopened and undamaged. Before use, visually inspect the sterile package to ensure that the seals have not been broken. Do not use the catheter if the integrity of the package has been compromised. Do not use catheter if its "Use Before Date," found on package labeling, has been passed.

**7.2 Inspection prior to use**

- Before use, examine carefully for defects, all of the equipment to be used. Do not use any equipment if it is damaged.

## 8. COMPATIBILITY

Product Features & Model Numbers Table

Model (Ref)	Working Length	Tip Shape / Type	Hydrophilic Coating Length	GW Compatibility	Min. Guide I.D.	Min. Sheath I.D.	Maximum Catheter O.D.	Marker Band Spacing	Throw Angle	
518-085	135 cm	Angled	100 cm	0.014"	5F	4F	3.2F	15 mm	45°	
518-087	150 cm	Angled		0.36 mm			0.042"			1.07 mm
518-089	90 cm	Angled	60 cm	0.018"	5F	4F	3.4F	15 mm	45°	
518-091	135 cm	Angled	100 cm	0.46 mm			0.044"			1.12 mm
518-093	150 cm	Angled								
518-077	65 cm	Angled	45 cm	0.035"	N/A	5F	4.5F	50 mm	45°	
518-079	90 cm	Angled	70 cm							
518-081	135 cm	Angled	115 cm							
518-083	150 cm	Angled	130 cm							

## 9. PROCEDURE SET-UP AND DIRECTIONS FOR USE

**Note:** Follow instructions for use for all equipment to be used with the Quick-Cross Select Support Catheter.

- Preparation: Using sterile technique, open the sterile package. Gently remove the protective hoop with the catheter from the pouch. Fill a sterile standard luer-lock syringe with sterile saline. Before removing the catheter from the hoop, connect the syringe to the catheter proximal luer fitting, flush the catheter and allow the saline to fill the hoop. Set catheter in hoop aside until ready for use.
- Insertion: Through a previously inserted, appropriately sized guiding catheter or introducer sheath, introduce the catheter over an appropriately sized guidewire (see specifications) using standard technique.
- Advancement: Use fluoroscopic guidance when advancing the catheter to the desired location within the vasculature.
- Removal: Gently withdraw the catheter using standard technique, being careful to maintain guidewire position if the guidewire is to remain in place.

**Infusion:** To perform infusion, withdraw the guidewire and reference the specifications for maximum infusion pressure.

**Note:** Do not exceed the maximum infusion pressures.

After use, dispose of all equipment in accordance with applicable requirements relating to hospital waste, and potentially bio-hazardous materials.

## 10. MANUFACTURER'S LIMITED WARRANTY

Manufacturer warrants that the Quick-Cross Select is free from defects in material and workmanship when used by the stated "Use By" date and when package is unopened and undamaged immediately before use. Manufacturer's liability under this warranty is limited to replacement or refund of the purchase price of any defective Quick-Cross Select. Manufacturer will not be liable for any incidental, special, or consequential damages resulting from use of the Quick-Cross Select. Damage to the Quick-Cross Select caused by misuse, alteration, improper storage or handling, or any other failure to follow these Instructions for Use will void this limited warranty. **THIS LIMITED WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** No person or entity, including any authorized representative or reseller of Manufacturer, has the authority to extend or expand this limited warranty and any purported attempt to do so will not be enforceable against Manufacturer.

## 11. NON-STANDARD SYMBOLS

Distal Marker Spacing	
Shape	
GW Compatibility	
Maximum OD	
Maximum PSI/kPa	
Sheath Compatibility	
Working Length	
Quantity	QTY
CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.	<b>Rx ONLY</b>

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©2018 Spectranetics Corporation
**Spectranetics Corporation**
**Spectranetics International B.V.**