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English

INSTRUCTIONS FOR USE

AngioSculpt® Percutaneous Transluminal Coronary Angioplasty (PTCA) Scoring Balloon Catheter

Rapid Exchange (RX) Delivery System

NOTE: These instructions apply to all balloon diameters and lengths.

STERILE: Sterilized with ethylene oxide gas. **Do not use if the package is open or damaged.**


CONTENTS: One (1) AngioSculpt® Scoring Balloon Catheter.

STORAGE: Store in a dry, dark, cool place.

I. DEVICE DESCRIPTION

The AngioSculpt catheter comprises a standard balloon catheter that incorporates a nitinol component. The proximal end of the catheter is a common PTCA catheter design comprising a hypotube connected to a plastic hub. The hub-hypotube contact area is supported by a relatively soft strain relief to avoid kinking. The hub is used to inflate the balloon and can be connected to a standard inflation device. The distal part of the catheter consists of a conventional nylon balloon and a laser-cut nitinol scoring element having three or more spiral struts that wrap around the balloon. The struts create focal concentrations of dilating force, which minimizes balloon slippage and assists in the luminal expansion of stenotic coronary arteries. The catheter is compatible with standard 0.014" coronary guide wires and 6F guiding catheters. The catheter length is approximately 137 cm and is offered in diameters of 2.0 – 3.5 mm in 0.5 mm increments. The working length of the scoring element ranges from 6 to 20 mm. The catheter is supplied sterile and is intended for a single use. Figure 1 below depicts the AngioSculpt catheter.

Figure 1: AngioSculpt Catheter

II. INDICATIONS

The AngioSculpt Scoring Balloon Catheter is indicated for the treatment of a hemodynamically significant coronary artery stenosis, including in-stent restenosis.

III. CONTRAINDICATIONS

The AngioSculpt catheter should not be used for the following:

- Coronary artery lesions unsuitable for treatment by percutaneous revascularization.
- Coronary artery spasm in the absence of a significant stenosis.

IV. WARNINGS

- This device is intended for single (one) use only. Do not resterilize and/or reuse, as this can potentially result in compromised device performance and increased risk of inappropriate reesterilization and cross contamination.
- To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.

- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.
- When the catheter is exposed to the vascular system, it should be manipulated while under high quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Balloon pressure should not exceed the rated burst pressure (RBP). [The RBP is based on results of in-vitro testing. At least 99.9% of the balloons (with 95% confidence) will not burst at or below their RBP. Use of a pressure monitoring device is recommended to prevent over-pressurization.]
- PTCA with the AngioSculpt device should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed on site (or at a nearby facility) in the event of a potentially injurious or life-threatening complication.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate the balloon.
- Proceed cautiously when using the AngioSculpt catheter in a freshly deployed bare metal or drug eluting stent. The AngioSculpt catheter has not been tested for post-dilation of stents or lesions distal to freshly deployed stents in clinical studies. Bench testing has shown no additional risk when inserting or withdrawing the AngioSculpt catheter through stents (no interference with stent struts, no retention of or damage to AngioSculpt catheter).
- Use the catheter prior to the "Use Before" (expiration) date specified on the package.

V. PRECAUTIONS

- Prior to angioplasty, the catheter should be examined to verify functionality, device integrity and to ensure that its size and length are suitable for the specific procedure for which it is to be used.
- Only physicians trained in the performance of percutaneous transluminal coronary angioplasty should use the catheter system.
- During and after the procedure, appropriate antiplatelet, anticoagulant and coronary vasodilator therapy consistent with institutional practice for coronary stent procedures should be provided to the patient.
- Do not rotate the catheter shaft in excess of 180 degrees when the tip is constrained.
- Do not rotate the catheter luer hub in excess of five (5) turns during use.
- Do not advance or retract the AngioSculpt catheter over the floppy portion of the guide wire.
- Catheter manipulation, including advancement and retraction, should be performed by grasping the hypotube shaft.
- If unusual resistance is felt when the catheter is being manipulated or if it is suspected that the guide wire has become kinked, carefully remove the entire catheter system (AngioSculpt catheter and steerable guide wire) as a unit.
- If fluoroscopic guidance indicates that the AngioSculpt catheter has advanced beyond the end of the guide wire, withdraw the catheter and reload the wire before advancing again.

VI. COMPLICATIONS

Possible complications include, but are not limited to, those listed in Table 1 below:

Table 1: Possible Complications

- Death
- Heart Attack (acute myocardial infarction)
- Total occlusion of the treated coronary artery
- Coronary artery dissection, perforation, rupture, or injury
- Pericardial tamponade
- No/slow reflow of treated vessel
- Emergency coronary artery bypass (CABG)
- Emergency percutaneous coronary intervention
- CVA/stroke
- Pseudoaneurysm
- Restenosis of the dilated vessel
- Unstable chest pain (angina)
- Thrombo-embolism or retained device components
- Irregular heart rhythm (arrhythmias, including life-threatening ventricular fibrillation)
- Severe low (hypotension)/high (hypertension) blood pressure
- Coronary artery spasm
- Hemorrhage or hematoma
- Need for blood transfusion
- Surgical repair of vascular access site
- Creation of a pathway for blood flow between the artery and the vein in the groin (arteriovenous fistula)
- Drug reactions, allergic reactions to x-ray dye (contrast medium)
- Infection

VII. SUMMARY OF CLINICAL STUDIES

Study Design

A multi-center, non-randomized, single-arm, prospective clinical study was conducted to evaluate the safety and efficacy of the AngioSculpt catheter in a wide range of coronary artery lesions in native vessels and following in-stent restenosis. The study population was adult patients scheduled to undergo clinically indicated percutaneous coronary intervention. Following informed consent, 45 patients meeting the study inclusion criteria underwent treatment with the AngioSculpt catheter.

Patient Selection

The following key inclusion criteria were utilized for patient selection:

- Planning to undergo clinically indicated percutaneous coronary intervention in native coronary arteries including in-stent restenosis
- Target lesion reference vessel size 2.0 – 4.0 mm diameter
- Target lesion length \leq 30mm
- Target lesion severity \geq 60% diameter stenosis

Patients meeting the above criteria were to be excluded for any of the following reasons:

- Target lesion in left main coronary artery
- Target lesion in degenerated saphenous vein graft
- Target lesion in a severe angulation of >60 degrees
- Target lesion distal to freshly deployed stent

- Target lesion demonstrating major dissection prior to deployment of AngioSculpt catheter
- Visible thrombus (by angiography) at target lesion site

Objectives

The primary safety objective was to demonstrate the incidence and severity of device-related complications (MACE: death, Q wave/non-Q wave MI, TLR) at one month follow-up.

The primary performance objective was to demonstrate successful percutaneous revascularization (defined as a reduction in target lesion diameter stenosis to \leq 50% following completion of all interventions and the absence of in-hospital MACE).

Methods

Patients were treated in standard fashion and according to institutional practice for percutaneous coronary interventions including the use of anti-coagulants (e.g. heparin) and antiplatelet agents (ASA, clopidogrel). There was no attempt to alter standard practice for these study patients other than using the investigational device.

Coronary angiography in the views best demonstrating the target lesion was performed prior to device deployment. When technically possible, IVUS of the target lesion was performed prior to device deployment at one of the two clinical sites.

A cineangiogram was taken of the deployed device immediately prior to and during inflation. Coronary angiography of the target lesion following completion of each device treatment (and prior to adjunctive stent placement) was performed in the original views. IVUS of the target lesion was performed following device treatment (and prior to adjunctive stent placement).

Additional clinically indicated interventions (e.g. stent placement) were performed and coronary angiography of the target lesion following completion of all interventions was performed in the original views.

Post-intervention medication was administered according to institutional protocol for percutaneous coronary interventions involving stents. EKG was recorded immediately post-procedure and at 24 hours or pre-discharge (whichever was earlier). CPK (and Troponin-I if CPK is abnormally elevated), at 4 and 12 hours post-procedure and at 24 hours/pre-discharge (whichever was earlier) were recorded.

Patients treated with the AngioSculpt catheter underwent a follow up examination 14-28 days from the interventional procedure which consisted of an office visit or telephone questionnaire for assessment of vital status, post-discharge myocardial infarction, post-discharge CABG surgery/PCI, angina class, and an EKG assessment.

Results

Following informed consent, 45 consecutive patients (age 63 ± 10.7 years, males 69%) referred for percutaneous coronary intervention and meeting the study inclusion criteria underwent treatment with the AngioSculpt catheter. The age and gender of the patients studied reflect the typical distribution in non-selected patients referred for percutaneous coronary intervention.

Patients were treated with standard peri-procedural medications including aspirin, clopidogrel, intravenous heparin and glycoprotein 2b/3a inhibitors at the discretion of the investigator.

All patients were available and underwent clinical follow-up at 14-28 days following treatment with the AngioSculpt catheter.

The AngioSculpt catheter was utilized in 45 patients and a total of 46 lesions. Of these 46 lesions, 32 were in native vessels and 14 were in-stent restenosis (ISR). The AngioSculpt was used as a standalone treatment in 10 lesions (9 ISR and 1 native vessel) and was used in combination with stenting in the remaining 36 lesions. In these cases, the AngioSculpt was intentionally undersized relative to the reference vessel diameter in order to facilitate stent placement (i.e. "pre-dilation").

There were no device-related major adverse clinical events (MACE: death, Q wave or non-Q wave myocardial infarction, ischemia driven target lesion revascularization) encountered during hospitalization or during the follow-up period (24.8 \pm 8.5 days). There were no device-related coronary artery perforations. There were no device malfunctions.

One patient was re-hospitalized during the follow-up period due to an episode of supra-ventricular tachycardia that was documented to be a pre-existing condition and unrelated to the AngioSculpt catheter. A second patient suffered a perforation of a diagonal branch coronary artery that was associated with a non-Q wave MI during treatment with a conventional angioplasty balloon and in a location and artery that were remote from the site treated with the AngioSculpt. This patient did not require surgical intervention and did not experience any additional MACE during follow-up.

The AngioSculpt catheter was successfully deployed in all 46 lesions. In all 46 attempted lesions, the primary performance endpoint of reduction of the lesion diameter stenosis to $\leq 50\%$ at the completion of the interventional procedure was successfully achieved. In all lesions treated, the AngioSculpt catheter demonstrated stable position during deployment without any significant slippage on angiography. The angiographic results are summarized in Table 2.

Table 2: Angiographic Results

	Pre-AS (n=46)	AS Alone (n=10)	AS Pre-Stent (n=36)	Post-Stent (n=36)
RVD (mm)	2.87±0.41	N/A	N/A	N/A
Length (mm)	15.67±6.14	N/A	N/A	N/A
DS (%)	75.27±12.91	17.46±8.15*	38.68±17.19*	3.81±3.75
MLD (mm)	0.75±0.35	2.49±0.43*	1.83±0.59*	2.91±0.47

RVD= reference vessel diameter

DS= diameter stenosis

MLD= minimum luminal diameter

*p<0.001 compared to pre-AngioSculpt

Intra-vascular ultrasound (IVUS) was performed pre and post treatment with the AngioSculpt catheter to evaluate the morphologic effects of the device on the plaque and to further confirm device safety. The IVUS results demonstrated plaque scoring and luminal expansion following treatment with the AngioSculpt catheter. There were no perforations or other evidence of unanticipated vessel injury on IVUS evaluation. The IVUS results are summarized in Table 3.

Table 3: IVUS Results

	Pre-AS (n=30)	ISR Post-AS (n=11)	De Novo Post-AS (n=19)	De Novo Post-Stent (n=19)
MLA (mm ²)	2.01±0.71	4.55±2.2*	2.65±0.9*	6.28±2.02

MLA= minimum luminal area

*p<0.001 compared to pre-AngioSculpt

The AngioSculpt catheter was successfully deployed in all 46 lesions. In four lesions, the stenosis was so severe that it could not be crossed initially by either the IVUS catheter or the AngioSculpt catheter and was therefore pre-dilated with a small balloon catheter (1.5/2.0 mm) following which the AngioSculpt catheter was successfully deployed.

There were no instances of device malfunction, retained device components or embolization. Each device was carefully inspected following completion of the procedure. There was no evidence of any damage or deterioration in any of the devices.

VIII. MATERIALS REQUIRED FOR USE WITH THE ANGIOSCULPT CATHETER:

WARNING - Use single use items only. Do not resterilize or reuse.

- Femoral, brachial or radial guiding catheter ($\geq 6F$)
- Hemostatic valve

- Contrast medium diluted 1:1 with normal saline
- Sterile heparinized normal saline
- 10cc and 20cc syringes for flushing and balloon prep
- Inflation device (indeflator)
- 0.014" coronary guide wire
- Guide wire introducer
- Guide wire torque device
- Radiographic contrast
- Manifold (for pressure monitoring and contrast injection), extension pressure tubing

IX. INSTRUCTIONS FOR USE

Prior to use of the AngioSculpt, examine carefully for damage and device integrity. Do not use if the catheter has bends, kinks, missing components or other damage.

1. Pre-medicate patients with ASA, Clopidogrel/Ticlopidine, intravenous anti-coagulants, coronary vasodilators and GP2b/3a blockers according to institutional protocol for percutaneous coronary interventions involving stents.
2. Perform coronary angiogram in the view best demonstrating the target lesion prior to device deployment.
3. Position 0.014" coronary guide wire of choice beyond the target lesion.
4. Using sterile technique, remove an appropriately sized (≤ 1.0 x reference vessel diameter (RVD)) AngioSculpt device from the sterile package and place on the sterile field.
5. Inspect the device to ensure that all components are intact.
6. Flush the guide wire lumen with saline by carefully inserting the distal catheter tip into the distal end of a 10 cc syringe and injecting saline until droplets emerge from the proximal guide wire lumen.
7. Attach 20 cc syringe filled with 2-3 cc of radiographic contrast to the catheter balloon inflation port.
8. Aspirate/remove air from the catheter balloon lumen using the 20 cc syringe filled with 2-3 cc of radiographic contrast and leave on vacuum for 30 seconds.
9. Gently release vacuum from the 20 cc syringe and remove it from the balloon inflation port.
10. Attach inflation device (indeflator), filled with 50:50 mixture of radiographic contrast and normal saline, to the balloon inflation port by creating a meniscus. Avoid introducing air bubbles into the catheter balloon lumen.
11. Aspirate using the inflation device, locking in vacuum.

NOTE: All air must be removed from the balloon and displaced with contrast medium prior to inserting into the body (repeat steps 9-11, if necessary).

12. Advance the AngioSculpt device over the coronary guide wire to the target lesion (may predilate with a 2.0 – 2.5 mm diameter balloon catheter if needed to cross the lesion).
- NOTE: When backloading the catheter onto the guide wire, the catheter should be supported, ensuring that the guide wire does not come in contact with the balloon.**
13. Inflate the AngioSculpt balloon per the following recommended protocol:
 - 2 atmospheres
 - Increase the inflation pressure by 2 atmospheres every 10-15 seconds until full device inflation is achieved
 - May inflate to a maximum pressure that is \leq RBP at the discretion of the physician (bearing in mind the estimated inflated diameter of the device at a given pressure)
14. Perform coronary angiogram (in the same view(s) as step 2) of the target lesion following completion of device treatment (and prior to adjunctive stent placement).

15. To remove the AngioSculpt device, apply negative pressure to the inflation device and confirm that the balloon is fully deflated. The catheter should be retracted only by grasping the hypotube shaft.
16. Inspect all components to ensure that the device is intact. Follow institutional procedures for disposal of biohazards. If device malfunction occurs or any defects are noted on inspection, flush the guide wire lumen and clean the outer surface of the device with saline, store the device in a sealed plastic bag contact Spectranetics® or the EC Authorized Representative for further instructions.
17. Complete any additional interventions as clinically indicated (e.g. stent placement).
18. Remove the coronary guide wire and perform coronary angiography (in the same view(s) as step 2) of the target lesion following completion of all interventions.
19. Remove all catheters and manage the arterial access site according to institutional protocol.
20. Continue treatment with ASA, Clopidogrel/Ticlopidine, and GP2b/3a blockers according to institutional protocol for percutaneous coronary interventions involving stents.

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Spectranetics International B.V.
Spectranetics Corporation