Contents of unopened, undamaged package are STERILE and NONPYROGENIC

Disposable: This device is intended for single-use only. Do not clean or resterilize. Store at room temperature.

Shelf life: The recommended shelf life is printed on each package. Storage beyond the expiration date may result in product deterioration.

CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician with appropriate training. Read all instructions prior to use

Indications for Use:

- The Quick-Access Needle Holders are intended to facilitate the placement of guidewires into the vascular system.

Contraindications:

- None.

Warnings:

- The devices are intended for one-time use only. Do not resterilize or reuse. Reuse and/or resterilization of the device will result in cross-patient blood contamination and compromised sterilization.
- Do not advance the guidewire or the Quick-Access Needle Holder if resistance is met.
- Use only appropriate guidewires of size 0.035" and below.

Precautions:

- Procedures should be conducted under fluoroscopic guidance with appropriate x-ray equipment.
- The sealed needle container should be inspected prior to opening. If the seal is broken, or the container has been damaged, sterility cannot be assured.
- Careful attention must be paid to maintain tight connection between the Quick-Access Needle Holder and needle.
- DO NOT resterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing.
- Reuse of this single use device could lead to serious patient injury or death.

Potential Complications or Adverse Events:

- Potential complications related to angioplasty include, but are not limited to clot formation and embolism, nerve damage, vascular perforation requiring surgical repair, damage to the vascular intima, or death.

How Supplied:

- Supplied sterilized and sealed in peel-open package. Intended for single use; do not re-sterilize, re-process or re-use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether packaging has been compromised.
- Upon removal from the package, thoroughly inspect the product to ensure that no damage has occurred during shipment.
- Do not use Quick-Access Needle Holder if its "Use By" date, found on the package labeling, has passed.

Compatibility:

- Guidewire: compatible with guidewire sizes up to 0.035"

Directions for Use:

- Remove the Quick-Access Needle Holder from its packaging and inspect the Quick-Access Needle Holder for damage before insertion.
- Purge the Quick-Access Needle Holder through lumen, observing for leaks and obstructions.
- Attach a needle to the Quick-Access Needle Holder.
- Insert a guidewire of size 0.035" and below until the guidewire tip is 3cm from the base of the Quick-Access Needle Holder's proximal tip.
- Advance the shut-off slider so that the guidewire is locked in position.
- Insert the needle into the target vessel and advance using the standard Seldinger technique.
- Once gaining confirmation of location in the artery via blood back flow, release the slider and insert the guidewire.
- Remove the needle and Quick-Access Needle Holder.
- Proceed with further percutaneous intervention such as balloon dilation and/or stent placement.
- Apply pressure to the insertion site according to standard practice or hospital protocol after removing all catheters.

Warning, Warranty, and Limitations:

- Needle Holder may fail to function for a variety of reasons, including medical complications and/or needle holder failure due to breakage.
- Needle Holder and accessories may be easily damaged before, during, or after insertion despite all efforts in design, manufacturing, and testing.
- No representation or warranty is made that the catheter will not fail or malfunction.
- Catheters and accessories are sold in an "as is" condition. All steps in the Instructions for Use are to aid the user in assessing the quality and performance of the catheter. The entire risk as to the quality and performance of the catheter is borne by the buyer.
- Upstream Peripheral Technologies, Ltd. disclaims all warranties, expressed or implied, with respect to the catheters and accessories, including but not limited to any implied warranty of merchantability or fitness for a particular purpose. Upstream Peripheral Technologies shall not be liable to any person for any medical expenses or any direct or consequential damages resulting from the use of any catheter or accessory or caused by any defect, failure, or malfunction of any catheter or accessory, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Upstream Peripheral Technologies to any representation or warranty with respect to catheters and accessories.

References:

- Gruntzig A.; Kumpe, D.; Technique of PTA with the Gruntzig Balloon Catheter, Amer Journ of Roent 132:547-552 (1979)
- Zeitler, E.; Gruntzig A.; Schoop, W.; (Eds), Percutaneous Vascular Recanalization , Springer-Verlag, Heidelberg/New York (1978)

Description of Graphical Symbols

		 0482	Use by 			
Manufacturer: Upstream Peripheral Technologies c/o ARAN Building 43 HaEshel Street P.O. Box 3067 Caesaria, Israel 38900 Tel: +972-4-623-9014 Email: drottenberg@aran-rd.com	Authorized European Representative: Mednet GmbH Borkstrasse 10 48163 Munster, Germany		 Caution, consult accompanying labeling and instructions	Batch code	Sterilize with ethylene oxide   Single use. Do not reuse or resterilize	REF Catalog number RX sale or use on physician order

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