English / English

Contents of unopened, undamaged package are STERILE and NONPYROGENIC

Disposable: This device is intended for single-use only. Do not clean or resterilize. Store at room temperature in a dry area away from heat and chemical fluids.

Shelf life: The recommended shelf life is printed on each package. Storage beyond the expiration date may result in balloon and catheter deterioration.

CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician with appropriate training. Read all instructions prior to use.

Indications for Use:

- The Quick-Cross Capture Guidewire Retriever is intended to be used in conjunction with a steerable guidewire to access discrete regions of the vasculature and for guidewire exchange. The Quick-Cross Capture Guidewire Retriever is not intended for use in the coronary, cerebral or carotid vasculature.

Contraindications:

- None.

Description

- The Quick-Cross Capture is a coaxial lumen design catheter with a balloon mounted on its distal tip. The outer lumen is for balloon inflation and is accessed through the side female luer at the proximal end of the Y-connector. The inner through lumen permits the use of guidewires to facilitate the advancement of the catheter through and around the stenosis. The wire lumen is accessed through the second port of the Y-connector. The Quick-Cross Capture Catheter is indicated for non-coronary and non-cerebral arteries (e.g., femoral, iliac, tibial, popliteal, and renal arteries). A removable balloon protector is provided. Two radiopaque bands identify the catheter.

Warnings:

- The devices are intended for one-time use only. Do not resterilize or reuse. Reuse and/or resterilization of the device will result in cross-patient blood contamination and compromised sterilization.
- Do not advance the guidewire or the Quick-Cross Capture if resistance is met.
- CAUTION: Pressure in excess of 5 ATM can cause balloon rupture.
- CAUTION: Use an inflation device with a pressure gauge for balloon inflation/deflation.
- Use only appropriate balloon inflation medium. Do not use air or gaseous medium to inflate the balloon.

Precautions:

- Procedures should be conducted under fluoroscopic guidance with appropriate x-ray equipment.
- The sealed catheter container should be inspected prior to opening. If the seal is broken, or the container has been damaged, sterility cannot be assured.
- Careful attention must be paid to maintain tight catheter connections and aspiration before proceeding to avoid air introduction into the system.
- Before removing the catheter from the sheath, it is important that the balloon is completely deflated.
- CAUTION: The catheter system should not be advanced against resistance. The cause of the resistance should be identified and corrective action should be taken. If resistance is felt upon catheter removal from an artery, then the system, guidewire, and the sheath should be removed together as a unit, particularly if balloon rupture or leakage is suspected.
- Damage may result from kinking, stretching, or forcefully wiping the catheter. Care should be used when handling.
- DO NOT resterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing.
- Reuse of this single use device could lead to serious patient injury or death.

Potential Complications or Adverse Events:

- Potential complications related to angioplasty include, but are not limited to clot formation and embolism, nerve damage, vascular perforation requiring surgical repair, damage to the vascular intima, or death.
- Potential complications related to introducing the catheter into the body include, but are not limited to, infection, hematoma formation, and air emboli.
- Potential balloon separation following balloon rupture or abuse and the subsequent need to use a snare or other medical interventional technique to retrieve the pieces.

How Supplied:

- Supplied sterilized and sealed in peel-open package. Intended for single use; do not re-sterilize, re-process or re-use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether packaging has been compromised.
- Upon removal from the package, thoroughly inspect the product to ensure that no damage, bends or kinks have occurred during shipment.
- Do not use catheter if its "Use By" date, found on the package labeling, has passed.

Compatibility:

- Sheath: Compatible with 6 French Introducer Sheath.
- Guidewires: Compatible with 0.014", 0.018", or 0.035" guidewires.

Directions for Use:

- Remove the Quick-Cross Capture from its packaging and inspect the catheter for damage before insertion.
- Advance the blue balloon protector of the catheter so that it covers the funnel at the catheter's distal tip to facilitate introduction to introducer. Purge the catheter's through lumen, observing for leaks and obstructions.
- When inflating the Balloon, use either 50/50 or 75/25 solution of saline and contrast medium, respectively.
- Attach an inflation device with the pressure gauge half filled with contrast solution to the balloon port of the catheter.
- Point inflation device with pressure gauge nozzle downward, aspirate until all air is removed from the balloon and bubbles no longer appear in the contrast solution.
- Prior to catheter insertion, confirm the appropriate guidewire for the size/length Quick-Cross Capture being used.
- Enter the vessel percutaneously over a 0.014", 0.018", or 0.035" guidewire, using the standard Seldinger technique. If an insertion sheath is used, be sure that it is a 6 French compatible sheath. After inserting the balloon protector into the sheath, advance the catheter, and then pull the protector to the proximal end of the catheter.
- In the true lumen, advance the Quick-Cross Capture using accepted angioplasty techniques (see references) until gaining confirmation fluoroscopically that the funnel, identified by the distal marker, is proximal to the occlusion.
- Pull the catheter backwards about 1cm, letting the funnel regain its original shape, to allow easier guidewire exchange.
- Inflate the balloon with contrast medium until 1.0 ATM is reached for the 6mm balloon or 3.0 ATM for the 8 mm balloon.
- Advance the guidewire across the lesion.
- Exchange guidewires as necessary, using the two radiopaque markers to confirm guidewire location and trajectory.
- Deflate the catheter's balloon by drawing a vacuum via an inflation device with a pressure gauge. The greater the vacuum applied and held during the withdrawal, the lower the deflated balloon profile.
- Gently withdraw the catheter, leaving the guidewire in place to complete the interventional procedure. As the balloon exits the vessel, use a smooth, gentle, steady motion. If resistance is felt upon removal, then the catheter should be removed under fluoroscopic guidance, particularly if the balloon ruptures or leakage is suspected.
- Then proceed with further percutaneous intervention such as balloon dilation and/or stent placement.
- Apply pressure to the insertion site according to standard practice or hospital protocol after removing all catheter

English / English

Upstream Peripheral Technologies Balloon Sizing Chart for 6mm Funnel

Applied Pressure	Diameter at applied pressure
1 ATM (nominal pressure)	6.5 mm

Upstream Peripheral Technologies Balloon Sizing Chart for 8mm Funnel

Applied Pressure	Diameter at applied pressure
3 ATM (nominal pressure)	8 mm

References:

- a) Gruntzig A.; Kumpe, D.; Technique of PTA with the Gruntzig Balloon Catheter, *Amer Journ of Roent* 132:547-552 (1979)
- b) Zeitler, E.; Gruntzig A.; Schoop, W.; (Eds), *Percutaneous Vascular Recanalization*, Springer-Verlag, Heidelberg/New York (1978)

Description of Graphical Symbols

	EC REP	CE 0482	Use by 	LOT	STERILE EO	
Manufacturer: Upstream Peripheral Technologies c/o ARAN Building 43 HaEshel Street P.O. Box 3067 Caesaria, Israel 38900 Tel: +972-4-623-9014 Email: drottenberg@aran-rd.com	Authorized European Representative: Mednet GmbH Borkstrasse 10 48163 Munster, Germany		 Caution, consult accompanying labeling and instructions	Batch code	Sterilize with ethylene oxide  	REF Catalog number RX sale or use on physician order

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