**Quick-Cross<sup>®</sup>**  
Support Catheter

# Instructions for Use

**Spectranetics<sup>®</sup>**

CE  
2797

## Instructions for Use - Sections by Language

<b>Pg</b>	<b>Language</b>		<b>Instructions for Use</b>
3	English	English	Instructions for Use
5	Bulgarian	български език	Указания за употреба
7	Croatian	Hrvatski	Upute za uporabu
9	Czech	Česky	Návod k Použití
11	Danish	Dansk	Brugsanvisning
13	Dutch	Nederlands	Gebruiksaanwijzing
15	Estonian	Eesti	Kasutusjuhend
17	Finnish	Suomi	Käyttöohje
19	French	Français	Mode d'emploi
21	German	Deutsch	Gebrauchsanweisung
23	Greek	Ελληνικά	Οδηγίες χρήσης
25	Hungarian	Magyar	Használati utasítások
27	Icelandic	Íslenska	Leiðbeiningar um notkun
29	Italian	Italiano	Istruzioni per l'uso
31	Latvian	Latviešu valoda	Lietošanas instrukcija
33	Lithuanian	lietuvių kalba	Naudojimo instrukcija
35	Norwegian	Norsk bokmål	Bruksinstruksjoner
37	Polish	Polski	Instrukcja użytkowania
39	Portuguese	Português	Instruções de Utilização
41	Romanian	Română	Instrucțiuni de utilizare
43	Russian	Русский язык	Инструкции по применению
45	Serbian	Српски језик	Uputstvo za upotrebu
47	Slovak	slovenčina	Návod na použitie
49	Slovenian	Slovenščina	Navodila za Uporabo
51	Spanish	Español	Instrucciones de Uso
53	Swedish	Svenska	Bruksanvisning
55	Turkish	Türkçe	Kullanım Talimatları

**Table of Contents**

<b>1.</b>	<b>DESCRIPTION.....</b>	<b>3</b>
<b>2.</b>	<b>INDICATIONS FOR USE.....</b>	<b>3</b>
<b>3.</b>	<b>SPECIFICATIONS.....</b>	<b>3</b>
<b>4.</b>	<b>CONTRAINDICATIONS.....</b>	<b>3</b>
<b>5.</b>	<b>WARNINGS.....</b>	<b>3</b>
<b>6.</b>	<b>PRECAUTIONS.....</b>	<b>3</b>
<b>7.</b>	<b>ADVERSE EVENTS.....</b>	<b>4</b>
<b>8.</b>	<b>HOW SUPPLIED.....</b>	<b>4</b>
<b>9.</b>	<b>COMPATIBILITY.....</b>	<b>4</b>
<b>10.</b>	<b>DIRECTIONS FOR USE.....</b>	<b>4</b>
<b>11.</b>	<b>WARRANTY INFORMATION.....</b>	<b>4</b>
<b>12.</b>	<b>NON-STANDARD SYMBOLS.....</b>	<b>4</b>


**1. DESCRIPTION**

The Spectranetics Quick-Cross Support Catheters are intravascular catheters, available in 9 models. All models have 3 radiopaque markers spaced equally along the distal shaft to aid in estimating geometry within the vascular system. The distal radiopaque marker is positioned within 3 mm of the distal catheter tip. A standard female luer is placed on the proximal end of each model. The distal 40 cm of each catheter model is coated with a lubricious, hydrophilic coating.

Model number 518-032 and 518-065 has a shaft of varying stiffness with a proximal shaft diameter of 3.0 Fr. tapering to a distal shaft diameter of 2.0 Fr and is compatible with a 0.014 inch or smaller guidewire.

Model numbers 518-033, 518-034, and 518-035 have a shaft of varying stiffness with a proximal shaft diameter of 3.4 Fr. tapering to a distal shaft diameter of 2.3 Fr and are compatible with a 0.018 inch or smaller guidewire.

Model numbers 518-036, 518-037, 518-038, and 518-066 have a shaft of varying stiffness with a proximal shaft diameter of 4.8 Fr. tapering to a distal shaft diameter of 3.8 Fr and are compatible with a 0.035 inch or smaller guidewire.

**2. INDICATIONS FOR USE**

The Spectranetics Quick-Cross Support Catheters are guide wire exchange and infusion devices designed for use in the vascular system. The catheters are intended to support a guidewire during access of vasculature, allow for exchange of guidewires, and provide a conduit for the delivery of saline solutions or diagnostic contrast agents.

**3. SPECIFICATIONS**

	518-032	518-065	518-033	518-034	518-035	518-066	518-036	518-037	518-038
<b>GW Compatibility, inch</b>	0.014	0.014	0.018	0.018	0.018	0.035	0.035	0.035	0.035
<b>Working Length, cm</b>	135	150	90	135	150	65	90	135	150
<b>Minimum guidewire length, cm</b>	180	180	150	180	180	150	150	180	180
<b>Distal marker spacing, mm</b>	15	15	15	15	15	50	50	50	50
<b>Outer Diameter (O.D.), inch</b>	0.039	0.039	0.044	0.044	0.044	0.063	0.063	0.063	0.063
<b>Distal Shaft diameter, inch</b>	0.026	0.026	0.030	0.030	0.030	0.050	0.050	0.050	0.050
<b>Tip Profile, inch</b>	0.020	0.020	0.023	0.023	0.023	0.041	0.041	0.041	0.041
<b>Minimum Guide Catheter, Fr.</b>	5	5	5	5	5	6	6	6	6
<b>Sheath Compatibility, Fr.</b>	4	4	4	4	4	5	5	5	5

**4. CONTRAINDICATIONS**

No known contraindications

**5. WARNINGS**

- This catheter should only be used by physicians qualified to perform percutaneous, vascular interventions.
- Catheter manipulation should only occur under fluoroscopy.
- The catheter should not be advanced through an area of resistance unless the source of resistance is identified by fluoroscopy and appropriate steps are taken to reduce or remove the obstruction.
- The catheter should not be advanced into a vessel having a diameter smaller than the catheter outer diameter.

**6. PRECAUTIONS**

- DO NOT resterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing. Reuse of this single use device could lead to serious patient injury or death and voids manufacturer warranties.
- Maximum infusion pressure is 300 psi.
- The catheter is designed and intended for intravascular use only.
- This catheter is designed and intended for one time use only.
- Careful inspection before use should verify that the catheter has not been damaged in shipment and that its condition is suitable for the procedure.
- Only use guidewires of the recommended diameter and length.
- If the catheter is used for infusion, reference the table of flow rates and ensure infusion pressure does not exceed the recommendations.
- Avoid introducing air or any other gas through the catheter into the vascular system.

**7. ADVERSE EVENTS**

Vascular catheterization and/or vascular interventions may result in complications including but not limited to:

- Vessel dissection, perforation, rupture or total occlusion
- Unstable angina
- Embolism
- Hypo/hypertension
- Acute myocardial infarction
- Arrhythmia, including ventricular fibrillation
- Death

**8. HOW SUPPLIED**

 The Spectranetics Quick-Cross Support Catheters are supplied **STERILE**. The devices are designated and designed for **SINGLE USE ONLY** and must not be resterilized and/or reused.

**8.1 Sterilization**

- The sterility of the product is guaranteed only if the package is unopened and undamaged. Before use, visually inspect the sterile package to ensure that the seals have not been broken. Do not use the catheter if the integrity of the package has been compromised. Do not use catheter if its "Use Before Date," found on package labeling, has been passed.

**8.2 Inspection Prior to Use**

- Before use, examine carefully for defects, all of the equipment to be used. Do not use any equipment if it is damaged.
- Store in a cool, dry place. Protect from direct sunlight and high temperature (greater than 60°C or 140°F).

**9. COMPATIBILITY**

- See the table in Specifications.
- After use, dispose of all equipment in accordance with applicable specific requirements relating to hospital waste, and potentially bio-hazardous materials.

**10. DIRECTIONS FOR USE**
**Procedure Set-Up**

- Follow the "Directions for Use" section below.

**Note: Follow instructions for use for all equipment to be used with the Quick-Cross® Support catheters. For example, guiding catheters, introducer sheaths, and guidewires.**

1. Preparation: Using sterile technique, open the sterile package. Gently remove the protective hoop with the catheter from the pouch. Fill a sterile standard luer-lock syringe with sterile saline. Before removing the catheter from the hoop, connect the syringe to the catheter proximal luer fitting, flush the catheter and allow the saline to fill the hoop. Set catheter in hoop aside until ready for use.
2. Insertion: Through a previously inserted, appropriately sized guiding catheter or introducer sheath, introduce the catheter over an appropriate sized guidewire (see specifications) using standard technique.
3. Advancement: Use fluoroscopic guidance when advancing the catheter to the desired location within the vasculature.
4. Removal: Gently withdraw the catheter using standard technique, being careful to maintain guidewire position if the guidewire is to remain in place.
5. Infusion: To perform infusion, withdraw the guidewire and reference the chart below. Note: Do not exceed 300 psi inlet infusion pressure.

**Quick-Cross Infusion Flow Rates (ml/second) at 150 and 300 psi Injection Pressures for Saline and Contrast Solutions**

Model	Size	Length	Sterile Saline		Contrast*	
			150 psi	300 psi	150 psi	300 psi
518-032	0.014	135	1.1	1.6	0.4	1.0
518-065	0.014	150	1.0	1.5	0.4	0.7
518-033	0.018	90	2.0	2.9	0.8	1.6
518-034	0.018	135	1.8	2.5	0.7	1.2
518-035	0.018	150	1.7	2.4	0.6	1.2
518-066	0.035	65	8.7	12.4	5.8	10.4
518-036	0.035	90	6.8	10.0	4.2	7.2
518-037	0.035	135	4.7	7.8	3.4	6.1
518-038	0.035	150	5.4	8.0	3.2	5.5

\* 75/25 Optiray 320 contrast / Sterile Saline mix

**11. WARRANTY INFORMATION**
**Manufacturer's Limited Warranty**

Manufacturer warrants that the Quick-Cross catheter is free from defects in material and workmanship when used by the stated "Use By" date and when package is unopened and undamaged immediately before use. Manufacturer's liability under this warranty is limited to replacement or refund of the purchase price of any defective Quick-Cross catheter. Manufacturer will not be liable for any incidental, special, or consequential damages resulting from use of the Quick-Cross catheter. Damage to the Quick-Cross catheter caused by misuse, alteration, improper storage or handling, or any other failure to follow these Instructions for Use will void this limited warranty. **THIS LIMITED WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** No person or entity, including any authorized representative or reseller of Manufacturer, has the authority to extend or expand this limited warranty and any purported attempt to do so will not be enforceable against Manufacturer.

**12. NON-STANDARD SYMBOLS**

Importer		Tip Profile	
GW Compatibility		Sheath Compatibility	
Distal Marker Spacing		Working Length	
Outer Diameter (O.D.)		Manufacturer	
Authorized Representative in the European Community			

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©2021 Spectranetics Corporation
**Spectranetics Corporation**