SightRailTM

Manual Dilator Sheath

Instructions for Use

 ***Spectranetics***[®]

CE
2797

Instructions for Use - Sections by Language

Pg	Language		Instructions for Use
3	English	English	Instructions for Use
7	Bulgarian	Български Език	Указания за Употреба
11	Croatian	Hrvatski	Upute za Uporabu
15	Czech	Česky	Návod k Použití
19	Danish	Dansk	Brugsanvisninger
23	Dutch	Nederlands	Gebruiksaanwijzing
27	Estonian	Eesti	Kasutusjuhis
31	Finnish	Suomi	Käyttöohjeet
35	French	Français	Mode d'emploi
39	German	Deutsch	Gebrauchsanweisung
43	Greek	Ελληνικά	Οδηγίες Χρήσης
47	Hungarian	Magyar	Használati Utasítások
51	Icelandic	Íslenska	Leiðbeiningar um Notkun
55	Italian	Italiano	Istruzioni per l'uso
59	Latvian	Latviešu Valoda	Lietošanas Pamācība
63	Lithuanian	Lietuvių Kalba	Naudojimo Instrukcija
67	Norwegian	Norsk Bokmål	Bruksinstruksjoner
71	Polish	Polski	Instrukcja Użytkowania
75	Portuguese	Português	Instruções de Utilização
79	Romanian	Română	Instrucțiuni de Utilizare
83	Russian	Русский Язык	Инструкции по Применению
87	Serbian	Српски Језик	Uputstvo za Upotrebu
91	Slovak	Slovenčina	Návod na Použitie
95	Slovenian	Slovenščina	Navodila za Uporabo
99	Spanish	Español	Instrucciones de Uso
103	Swedish	Svenska	Bruksanvisning
107	Turkish	Türkçe	Kullanım Talimatları

Table of Contents

1. DEVICE DESCRIPTION.....	3
2. INDICATIONS FOR USE.....	3
3. CONTRAINDICATIONS.....	3
4. WARNINGS.....	3
5. PRECAUTIONS.....	3
6. POTENTIAL ADVERSE EVENTS.....	4
7. INDIVIDUALIZATION OF TREATMENT.....	4
8. HOW SUPPLIED.....	4
9. COMPATIBILITY.....	4
10. DIRECTIONS FOR USE.....	4
11. MANUFACTURER'S LIMITED WARRANTY.....	5
12. NON-STANDARD SYMBOLS.....	6

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

1. DEVICE DESCRIPTION

The SightRail Manual Dilator Sheath set consists of an inner and outer polymer sheath that may be used individually or as a set in a telescoping fashion. The inner sheath terminates with beveled tips. The outer sheath terminates with one beveled tip and one blunt tip. Reference Figure 1.

Figure 1: SightRail Manual Dilator Sheath Tip Configuration

Exterior stripe markings aligned with the beveled tips permit visual identification of sheath orientation. Arrow markings on the inner sheath indicate a region in which the inner and outer sheath tips are aligned. Reference Figure 2.


Figure 2: SightRail Manual Dilator Sheath Marking Features

The sheath set size is also marked on the inner sheath. There are multiple diameter and length options available, which are referenced in Table 1.

Use with other devices

The SightRail Manual Dilator Sheath may be used in conjunction with the Spectranetics Lead Locking Device (LLD™). Follow the "Instructions for Use" for other devices used.

Table 1 provides SightRail Manual Dilator Sheath models and sizing specifications.

Table 1: SightRail Manual Dilator Sheath Models and Sizes

Model	Size	Inner / Outer Sheath Length (cm)
550-008	8.5 F	43 / 33
550-010	10 F	43 / 33
550-011	11.5 F	43 / 33
550-013	13 F	43 / 33
555-508	8.5 F Long	51 / 41
555-510	10 F Long	51 / 41
555-511	11.5 F Long	51 / 41
555-513	13 F Long	51 / 41

2. INDICATIONS FOR USE

The SightRail Manual Dilator Sheath is intended for use in patients requiring the percutaneous dilation of tissue to facilitate the removal of cardiac leads.

3. CONTRAINDICATIONS

None known.

4. WARNINGS

- The SightRail Manual Dilator Sheaths should be used at institutions with cardiothoracic surgical capabilities by physicians knowledgeable in the techniques and devices for lead removal. Complication prevention and management protocols should be in place and routinely practiced. The recommendations for lead management of the Heart Rhythm Society¹ (HRS) and European Heart Rhythm Association² (EHRA) are highly recommended for best results.
- When using the SightRail Manual Dilator Sheaths, DO NOT insert sheaths over more than one lead at the same time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.
- Maintain appropriate traction on the lead being extracted during advancement of the inner or outer SightRail Manual Dilator Sheaths.
- DO NOT maintain a stationary position with SightRail Manual Dilator Sheath tips at the Superior Vena Cava (SVC)-Right Atrial (RA) junction as it may result in damage to this delicate area during subsequent lead extraction and reinsertion procedures (e.g., manipulating the dilator sheath or implanting a new lead).

5. PRECAUTIONS

- Do not resterilize or re-use this device, as these actions can compromise device performance and/or increase the risk of cross-contamination due to inappropriate reprocessing.
- Re-use of this single use device could lead to serious patient injury or death

and voids manufacturer warranties.

- Do not alter the sheath from its original state prior to use.
- When the SightRail Manual Dilator Sheath is in the body, it should only be manipulated under fluoroscopic observation with radiographic equipment that provides high quality images.
- Prior to the procedure, evaluate the physical dimensions of the lead in relation to the specifications of the dilator sheath to determine compatibility.
- If selectively removing leads with the intent to leave one or more chronic leads implanted intact, the non-targeted leads must be subsequently tested to ensure that they were not damaged or dislodged during the procedure.
- When advancing dilator sheaths, use proper sheath technique. Maintain adequate tension and coaxial alignment on the lead to minimize the risk of vessel wall or cardiac structure damage.
- If excessive scar tissue or calcification prevents safe advancement of dilator sheaths, consider an alternate approach.
- Excessive force with dilator sheaths used intravascularly may result in damage to the vascular system requiring emergency surgical repair.
- If the lead breaks, evaluate fragment for retrieval.
- If hypotension develops, rapidly evaluate; treat as appropriate.
- Due to rapidly evolving lead technology, this device may not be suitable for dilation of tissue around all types of leads. If there are questions or concerns regarding compatibility of this device with particular leads, contact the lead manufacturer.
- Do not pull on the lead because it may stretch, distort, or break, making subsequent removal more difficult. Damage to a lead may prevent passage of a lead locking device through the lumen and/or make dilation of scar tissue more difficult.
- When removing a chronic pacing lead, be aware that if it is freed spontaneously during the extraction procedure, the lead tip may become trapped in the upper vasculature. Dilator sheaths, advanced at least to the innominate vein, are often necessary to extract the lead tip through the scar tissue at the site of venous entry, and to avoid a venotomy.
- If the dilator sheath fails to progress after initial success, or if advancing the sheath was difficult, remove the sheaths one at a time to inspect the tips. If the inner sheath tip is distorted or frayed, the other end may be used. A new sheath set may also be used to continue treatment.
- When advancing a sheath around a bend, be cognizant of the sheath's beveled tip orientation.

6. POTENTIAL ADVERSE EVENTS

Potential adverse events related to the procedure of intravascular removal of leads include (listed generally in order of increasing potential effect):

- Dislodging or damaging non-targeted lead
- Chest wall hematoma
- Thrombosis
- Arrhythmias
- Bacteremia
- Hypotension
- Pneumothorax
- Migrating fragment from lead
- Migration of vegetation from lead
- Pulmonary embolism
- Laceration or tearing of vascular structures or the myocardium
- Hemopericardium
- Cardiac tamponade
- Hemothorax
- Stroke
- Death

7. INDIVIDUALIZATION OF TREATMENT

Weigh the relative risks and benefits of intravascular lead removal procedures in cases when:

- Dual coil ICD leads are being removed;
- The lead to be removed has a sharp bend or evidence of fracture
- The lead shows evidence of insulation disintegration raising the concern of pulmonary embolism
- Vegetations are attached directly to the lead body

8. HOW SUPPLIED

8.1 Sterilization

- Single use only; not for resterilization or reprocessing
- Ethylene oxide sterilized
- Non-pyrogenic
- Sterility guaranteed if package is unopened and undamaged
- Store device in a dry cool place (below 60°C / 140°F) until use.

8.2 Inspection Prior to Use

- Prior to use, visually inspect the sterile package to ensure that seals have not been broken. The sheaths should be examined carefully for defects or damage. Do not use if the sheaths have apparent defects or are damaged or the Use-by Date has been exceeded.

9. COMPATIBILITY

Information for determining SightRail Manual Dilator Sheath dimensional compatibility is shown in Table 1.

10. DIRECTIONS FOR USE

10.1 Procedure Setup

- Obtain a thorough patient history, including blood type. Appropriate blood products should be readily available.
- Determine the manufacturer, model number, physical dimensions, and implant date of the target lead. Perform radiographic evaluation of position, type and condition of target lead.
- Use a procedure room that has high quality fluoroscopy, pacing equipment, defibrillator, and thoracotomy and pericardiocentesis trays.

- Prep and drape the patient's chest for possible thoracotomy.
- Arrange for immediate surgical back-up.
- Establish back-up pacing as needed.
- Using sterile technique, open the sterile package.
- Make available, and open as needed, any other adjunct devices such as lead locking devices or related accessories.

10.2 Clinical Technique

1. Surgically expose the proximal end of the target lead(s) and remove the lead(s) from its/their connections (if connected).
2. Remove all suture and tie-down materials.
3. Cut off all proximal fittings, if present, using clippers or other cutters. It is important to cut the lead very close to the connector (but past any crimp joints) leaving as long a portion of the target lead to work with as possible. Avoid closing off the interior lumen (or coil) of the lead when cutting it.
4. Advance a lead locking device down the lead lumen as distal as possible and deploy the locking mechanism.
5. Unless the lead exterior is damaged, degraded or too thin, tie a suture at the proximal end of the lead to use as a traction element (the suture may be attached to the proximal portion of a lead locking device).
6. For an active fixation lead, attempt to unscrew the lead fixation mechanism via counterclockwise rotation by rotating the lead (and lead locking device, if used) counterclockwise.
7. Gently apply traction to the lead to determine if it is still engaged in tissue. If the lead is sufficiently free of binding tissue, gently pull on the lead locking device (if used) and lead to remove it.
8. If the lead is not removed from the vessel with gentle pulling, dilator sheaths (or other retrieval devices) may help separate the lead from any tissue encapsulation by advancing the dilator sheath over the target lead and any accompanying lead locking device.
9. The SightRail Manual Dilator Sheath set is assembled by inserting the inner sheath into the proximal end of the outer sheath until it appears beyond the distal end of the outer sheath. The SightRail Manual Dilator Sheaths may be advanced in alternating, telescoping fashion to progressively dilate tissue along the length of the lead. Alternatively, the inner or outer sheath can be used independently.

Exterior markings on the SightRail Manual Dilator Sheath may be used to aid in identification of the approximate beveled tip orientation and relative position of the inner and outer sheaths. Always use fluoroscopic monitoring when manipulating the sheaths within the vascular system.

10. Always maintain adequate tension and coaxial alignment on the lead to support the maneuvering of dilator sheaths to properly guide them within the patient anatomy.

With too little tension, the sheaths may damage the vein. Too much tension may cause a myocardial avulsion. The size of the dilator sheath should be large enough that the sheath can be advanced over the lead without causing the lead to buckle or its exterior jacket to distort, but the sheath should not be overly loose.

Rotating the sheath(s) during advancement may facilitate progress through resistant scar tissue.

11. For cardiac leads, if the lead has not been freed by the time the dilator sheath nears the myocardial lead fixation point, position the sheath end against the myocardium. Countertraction may be applied to dilate the remaining tissue at the lead tip. This is accomplished by holding the dilator sheath about one centimeter from the myocardium while steadily and gently applying tension to the lead. Rotation of the sheath may help dilate remaining tissue at the lead tip. The blunt end of the outer sheath may be used as an alternative countertraction platform.
12. When all binding tissue surrounding the target lead has been successfully dilated, the lead may slide freely out of the body with applied traction.
13. Once a lead is removed from the patient, the sheath may then be used as a conduit for a guidewire to facilitate the implantation of a new lead.
14. At the completion of the extraction procedure, withdraw the SightRail Manual Dilator Sheath(s) from the patient and dispose of the product according to local biological handling and disposal procedures.

References

1. Wilkoff B.L., et al. Transvenous Lead Extraction: Heart Rhythm Society Expert Consensus on Facilities, Training, Indications, and Patient Management. Heart Rhythm. July 2009.
2. Deharo J.C., et al. Pathways for training and accreditation for transvenous lead extraction: a European Heart Rhythm Association position paper. Europace (2012) 14, 124-134
3. Verma A, Wilkoff BL. Intravascular pacemaker and defibrillator lead extraction: A state-of-the-art review. Heart Rhythm, Vol. 1, No. 6, December 2004; 739-745.
4. Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
5. Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.
6. Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991;101:989-997.
7. Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. PACE 1990; 13:1871-1875.

11. MANUFACTURER'S LIMITED WARRANTY

Manufacturer warrants that SightRail Manual Dilator Sheaths are free from defects in material and workmanship when used by the stated "Use By" date and when package is unopened and undamaged immediately before use. Manufacturer's liability under this warranty is limited to replacement or refund of the purchase price of any defective SightRail Manual Dilator Sheath. Manufacturer will not be liable for any incidental, special, or consequential damages resulting from use of the SightRail Manual Dilator Sheath. Damage

Instructions for Use

English / English

to the SightRail Manual Dilator Sheath caused by misuse, alteration, improper storage or handling, or any other failure to follow these Instructions for Use will void this limited warranty. **THIS LIMITED WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** No person or entity, including any authorized representative or reseller of Manufacturer, has the authority to extend or expand this limited warranty and any purported attempt to do so will not be enforceable against Manufacturer.

12. NON-STANDARD SYMBOLS

Quantity	QTY
Inner Sheath Minimum Inner Diameter	
Inner Sheath Maximum Outer Diameter	
Inner Sheath Length	
Outer Sheath Minimum Inner Diameter	
Outer Sheath Maximum Outer Diameter	
Outer Sheath Length	
Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.	
Importer	
Size	
Long	

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Références

1. Wilkoff B.L., et al. Transvenous Lead Extraction: Heart Rhythm Society Expert Consensus on Facilities, Training, Indications, and Patient Management. Heart Rhythm. July 2009.
2. Deharo J.C., et al. Pathways for training and accreditation for transvenous lead extraction: a European Heart Rhythm Association position paper. Europace (2012) 14, 124-134

3. Verma A, Wilkoff BL. Intravascular pacemaker and defibrillator lead extraction: A state-of-the-art review. Heart Rhythm, Vol. 1, No. 6, December 2004; 739-745.
4. Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
5. Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.
6. Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991;101:989-997.
7. Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. PACE 1990; 13:1871-1875.

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- Wilkoff B.L., et al. Transvenous Lead Extraction: Heart Rhythm Society Expert Consensus on Facilities, Training, Indications, and Patient Management. Heart Rhythm. July 2009.
- Deharo J.C., et al. Pathways for training and accreditation for transvenous lead extraction: a European Heart Rhythm Association position paper. Europace (2012) 14, 124-134
- Verma A, Wilkoff BL. Intravascular pacemaker and defibrillator lead extraction: A state-of-the-art review. Heart Rhythm, Vol. 1, No. 6, December 2004; 739-745.
- Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
- Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.
- Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991;101:989-997.
- Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. PACE 1990; 13:1871-1875.

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
Atsauces

- Wilkoff B.L., et al. Transvenous Lead Extraction: Heart Rhythm Society Expert Consensus on Facilities, Training, Indications, and Patient Management. Heart Rhythm. July 2009.
- Deharo J.C., et al. Pathways for training and accreditation for transvenous lead extraction: a European Heart Rhythm Association position paper. Europace (2012) 14, 124-134
- Verma A, Wilkoff BL. Intravascular pacemaker and defibrillator lead extraction: A state-of-the-art review. Heart Rhythm, Vol. 1, No. 6, December 2004; 739-745.
- Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
- Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992: 10:735-748.
- Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991;101:989-997.
- Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. PACE 1990: 13:1871-1875.

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Referințe

1. Wilkoff B.L., et al. Transvenous Lead Extraction: Heart Rhythm Society Expert Consensus on Facilities, Training, Indications, and Patient Management. Heart Rhythm. July 2009.
2. Deharo J.C., et al. Pathways for training and accreditation for transvenous lead extraction: a European Heart Rhythm Association position paper. Europace (2012) 14, 124-134
3. Verma A, Wilkoff BL. Intravascular pacemaker and defibrillator lead extraction: A state-of-the-art review. Heart Rhythm, Vol. 1, No. 6, December 2004; 739-745.
4. Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.

5. Byrd CL, et al. Lead extraction: indications and techniques. *Cardiol Clin* 1992; 10:735-748.
6. Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. *J Thorac Cardiovasc Surg* 1991;101:989-997.
7. Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. *PACE* 1990; 13:1871-1875.

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©2020 Spectranetics Corporation
Spectranetics Corporation
Spectranetics International B.V.