## Instructions for Use

### Instructions for Use - Sections by Language

Pg	Language		Instructions for Use
3	English	English	Instructions for Use
8	Chinese	中文	使用说明


**Table of Contents**

1.	DESCRIPTION .....	3
2.	INDICATIONS FOR USE .....	3
3.	CONTRAINDICATIONS .....	3
4.	WARNINGS .....	3
5.	PRECAUTIONS .....	3
6.	ADVERSE EVENTS .....	4
7.	HOW SUPPLIED .....	4
8.	COMPATIBILITY .....	4
9.	DIRECTIONS FOR USE .....	5
10.	MANUFACTURER'S LIMITED WARRANTY .....	6
11.	NON-STANDARD SYMBOLS .....	6

**US Only CAUTION: Federal law restricts this device to sale by or on the order of a physician.**

**1. DESCRIPTION**

The VisiSheath Dilator Sheath is a single lumen sheath used independently or as a support for an inner sheath to facilitate tissue dilation. One end is terminated with a 45° angle cut, while the other end is blunt, see Figure 1. Both ends contain a radiopaque marker band to enable fluoroscopic identification of tip location and orientation. An additional exterior mark aligned with the tip of the 45° angle cut permits visual identification of sheath orientation. There are multiple diameter and length options available, see Table 1.

**Figure 1: VisiSheath Dilator Sheath Tip Configuration**

**2. INDICATIONS FOR USE**

The VisiSheath Dilator Sheath is intended for use in patients requiring the percutaneous dilation of tissue surrounding cardiac leads, indwelling catheters and foreign objects. The device is also intended for use in the introduction and support of intravascular catheters.

**3. CONTRAINDICATIONS**

None known.

**4. WARNINGS**

- Dilator sheaths should be used only at institutions with thoracic surgical capabilities by physicians knowledgeable in the techniques and devices for lead or catheter removal.
- When using dilator sheaths, do not insert sheaths over more than one lead or catheter at a time. Severe vessel damage, including venous wall laceration requiring surgical repair, may occur.
- Do not maintain a stationary position with the VisiSheath tip at the Superior Vena Cava (SVC)-right atrial (RA) junction as it may result in damage to this delicate area during subsequent lead extraction and reinsertion procedures (e.g., manipulating the dilator sheath or implanting a new lead).
- Weigh the relative risks and benefits of intravascular lead/catheter dilation procedures especially in cases when:
  - the object to be dilated away from adherent tissue is of a dangerous shape or configuration,
  - the likelihood of lead/catheter disintegration may result in increased risk of fragment embolization,
  - vegetations are attached directly to the lead/catheter body.

**5. PRECAUTIONS**

**For single use only.** The VisiSheath Dilator Sheath must not be resterilized and/or reused.

- DO NOT resterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing. Reuse of this single use device could lead to serious patient injury or death and voids manufacturer warranties.
- Do not alter the sheath from its original state prior to use.
- When the VisiSheath Dilator Sheath is in the body, it should be manipulated only under fluoroscopic observation with radiographic equipment that provides high quality images.
- Prior to the procedure, evaluate the physical dimensions of the lead, catheter, or inner sheath in relation to the specifications of the dilator sheath to determine possible incompatibility.
- If selectively removing leads/catheters with the intent to leave one or more chronic leads/catheters implanted intact, the non-targeted leads/catheters must be subsequently tested to ensure that they were not damaged or dislodged during the procedure.
- When advancing dilator sheaths, use proper sheath technique. Maintain adequate tension and coaxial alignment on the lead/catheter to minimize the risk of vessel wall damage.

## English / English

- If excessive scar tissue or calcification prevents safe advancement of dilator sheaths, consider an alternate approach.
- Excessive force with dilator sheaths used intravascularly may result in damage to the vascular system requiring emergency surgical repair.
- If the lead/catheter breaks, evaluate fragment for retrieval.
- If hypotension develops, rapidly evaluate; treat as appropriate.
- Due to rapidly evolving lead/catheter technology, this device may not be suitable for dilation of tissue around all types of leads/catheters. If there are questions or concerns regarding compatibility of this device with particular leads/catheters, contact the lead/catheter manufacturer.
- Do not pull on the lead/catheter because it may stretch, distort, or break, making subsequent removal more difficult. Damage to a lead may prevent passage of a lead locking device through the lumen and/or make dilation of scar tissue more difficult.
- If a lead locking device has not been used, be aware that damage to the lead caused by pulling on it may prevent subsequent passage of a lead locking device through the lumen and/or make dilation of scar tissue more difficult.
- If removing a chronic pacing lead, be aware that if it is freed spontaneously during the extraction procedure, the lead tip may become trapped in the upper vasculature. Dilator sheaths, advanced at least to the innominate vein, are often necessary to extract the lead tip through the scar tissue at the site of venous entry, and to avoid a venotomy.
- If the dilator sheath fails to progress after initial success, or if advancing the sheath was difficult, remove the sheath to inspect the tip. If the tip is distorted or frayed, exchange the damaged sheath for a new sheath before continuing treatment.
- When advancing a sheath around a bend, keep the point of the sheath's beveled tip oriented toward the inside of the bend.

## 6. ADVERSE EVENTS

### Potential Adverse Events

Potential adverse events related to the procedure of intravascular removal of leads/catheters include (listed generally in order of increasing potential effect):

- dislodging or damaging non-targeted lead/catheter
- chest wall hematoma
- thrombosis
- arrhythmias
- bacteremia
- hypotension
- pneumothorax
- migrating fragment from lead/catheter
- migration of vegetation from lead/catheter
- pulmonary embolism
- laceration or tearing of vascular structures or the myocardium
- hemopericardium
- cardiac tamponade
- hemothorax
- stroke
- death

## 7. HOW SUPPLIED

### 7.1 Sterilization

For single use only. Do not re-sterilize and/or reuse. The VisiSheath Dilator Sheath is supplied sterile and non-pyrogenic. Sterility is guaranteed only if the packaged is unopened and undamaged.

### 7.2 Device / Packaging Care

Store devices in a dry cool place (below 60° C / 140° F) until use.

### 7.3 Inspection Prior to Use

Prior to use, visually inspect the sterile package to ensure that seals have not been broken. The sheath should be examined carefully for defects or damage. Do not use if it has apparent defects or is damaged.

## 8. COMPATIBILITY

### Use with other devices

The VisiSheath Dilator Sheath may be used as a support sheath for compatibly sized inner sheaths, including the Spectranetics Laser Sheath (SLS™II) used for removal of cardiac leads. It may also be used in conjunction with a Spectranetics Lead Locking Device (LLD®).

Be sure to closely follow the "Instructions for Use" for each device used.

Table 1 provides physical dimensions and specifications for use in determining VisiSheath compatibility with other devices.

English / English

Table 1: VisiSheath Dilator Sheath Size Compatibility Specifications

Model (Ref)	Length (cm)	Sheath Size (Diameter)	Minimum Inner Diameter (mm/in/F)	Maximum Outer Diameter (mm/in/F)	SLS™II Compatibility
501-012	43	S	4.2 /	5.5 /	12 F
501-112	33		0.168 /	0.215 /	
501-212	23		12.8	16.4	
501-014	43	M	5.0 /	6.5 /	14 F
501-114	33		0.198 /	0.253 /	
501-214	23		15.0	19.3	
501-016	43	L	5.9 /	7.5 /	16 F
501-116	33		0.236 /	0.293 /	
501-216	23		17.9	22.4	

**9. DIRECTIONS FOR USE**

**9.1 Procedure Setup**

- Obtain a thorough patient history, including blood type. Appropriate blood products should be readily available.
- Determine the manufacturer, model number, physical dimensions, and implant date of the target lead/catheter. Perform radiographic evaluation of position, type and condition of target lead/catheter.
- Use a procedure room that has high quality fluoroscopy, pacing equipment, defibrillator, and thoracotomy and pericardiocentesis trays.
- Prep and drape the patient’s chest for possible thoracotomy.
- Arrange for immediate surgical back-up.
- Establish back-up pacing as needed.
- Using sterile technique, open the sterile package.
- Make available, and open as needed, any other adjunct devices such as inner sheaths (laser or non-laser), lead locking devices, or related accessories.

**9.2 Clinical Technique**

1. Surgically expose the proximal end of the target lead/catheter(s) and remove the lead/catheter from its connections (if connected).
2. Remove all suture and tie-down materials.
3. Cut off all proximal fittings, if present, using clippers or other cutters. It is important to cut the lead/catheter very close to the connector (but past any crimp joints) leaving as long a portion of the target catheter/lead(s) to work with as possible. Avoid closing off the interior lumen (or coil) of the lead/catheter when cutting it.

**PRECAUTION:** Do not pull on the lead/catheter because it may stretch, distort, or break, making subsequent removal more difficult. Damage to a lead may prevent passage of a lead locking device through the lumen and/or make dilation of scar tissue more difficult.

4. Advance a lead locking device down the lead lumen as distal as possible and deploy the locking mechanism.
5. Unless the lead/catheter exterior is damaged, degraded or too thin, tie a suture at the proximal end of the lead/catheter to use as a traction element (the suture may be attached to the proximal portion of a lead locking device).
6. For an active fixation lead/catheter, attempt to unscrew the lead/catheter fixation mechanism via counterclockwise rotation by rotating the catheter/lead (and lead locking device if used) counterclockwise.
7. Gently apply traction to the lead/catheter to determine if it is still engaged in tissue. If the lead is sufficiently free of binding tissue, gently pull on the lead locking device (if used) and catheter/lead to remove it.

**PRECAUTION:** If a lead locking device has not been used, be aware that damage to the lead caused by pulling on it may prevent subsequent passage of a lead locking device through the lumen and/or make dilation of scar tissue more difficult.

**PRECAUTION:** If removing a chronic pacing lead, be aware that if it is freed spontaneously during the extraction procedure, the lead tip may become trapped in the upper vasculature. Dilator sheaths, advanced at least to the innominate vein, are often necessary to extract the lead tip through the scar tissue at the site of venous entry, and to avoid a venotomy.

8. If the lead/catheter is not removed from the vessel with gentle pulling, dilator sheaths (or other retrieval devices) may help separate the lead/catheter from any tissue encapsulation by advancing the dilator sheath over the target lead/catheter and any accompanying lead locking device.
9. A size-compatible inner sheath, including a laser sheath, may be used with a VisiSheath Dilator Sheath acting as a support sheath. The VisiSheath Dilator Sheath is loaded over the inner sheath by inserting the distal end of inner sheath through the proximal end of the dilator sheath until it appears beyond the proximal end of the

## English / English

dilator sheath. The inner sheath and VisiSheath Dilator Sheath may be advanced in alternating fashion to progressively dilate tissue along the length of the lead/catheter.

- Always maintain adequate tension and coaxial alignment on the lead/catheter to support the maneuvering of the dilator sheaths and any accompanying inner sheaths to properly guide them within the patient anatomy.

With too little tension, the sheaths may damage the vein. Too much tension may cause a myocardial avulsion. The size of the dilator sheath should be large enough that the sheath can be advanced over the lead/catheter without causing the lead/catheter to buckle or its exterior jacket to distort, but the sheaths should not be overly loose.

Rotating the sheaths during advancement may facilitate progress through resistant scar tissue.

Always use fluoroscopic monitoring when manipulating the dilator sheath within the vascular system. Never use sheaths on more than one lead/catheter at a time.

**PRECAUTION:** If the dilator sheath fails to progress after initial success, or if advancing the sheath was difficult, remove the sheath to inspect the tip. If the tip is distorted or frayed, exchange the damaged sheath for a new sheath before continuing treatment.



**PRECAUTION:** When advancing a sheath around a bend, keep the point of the sheath's beveled tip oriented toward the inside of the bend.

- For cardiac leads, if the lead has not been freed by the time the dilator sheath nears the myocardial lead fixation point, position the sheath end against the myocardium, with any inner sheath retracted several centimeters into the dilator sheath. Countertraction may be applied to dilate the remaining tissue at the lead tip. This is accomplished by holding the dilator sheath about one centimeter from the myocardium while steadily and gently applying tension to the lead. Rotation of the sheath may help dilate remaining tissue at the lead tip.
- When all binding tissue surrounding the target lead/catheter has been successfully dilated, the lead/catheter may slide freely out of the body with applied traction.

### 10. MANUFACTURER'S LIMITED WARRANTY

Manufacturer warrants that the VisiSheath is free from defects in material and workmanship when used by the stated "Use By" date and when package is unopened and undamaged immediately before use.

### 11. NON-STANDARD SYMBOLS

<b>Working Length</b>		<b>Non-Pyrogenic</b>	
<b>Inner Diameter</b>		<b>Outer Diameter</b>	
<b>Sterilized Using Ethylene Oxide</b>		<b>Do not use if package is damaged</b>	
<b>Upper Limit of Temperature</b>		<b>Consult Instructions for Use</b>	
<b>Single Use</b>		<b>Keep Dry</b>	
<b>Quantity</b>	<b>QTY</b>	<b>Catalog Number</b>	
<b>CAUTION: Federal law restricts this device to sale by or on the order of a physician.</b>			

### References

- Verma A, Wilkoff BL. Intravascular pacemaker and defibrillator lead extraction: A state-of-the-art review. Heart Rhythm, Vol. 1, No. 6, December 2004; 739-745.
- Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
- Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.
- Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991;101:989-997.
- Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. PACE 1990; 13:1871-1875.

English / English

**【Registrant】**

Name of the Registrant: Spectranetics Corporation  
Residence: 9965 Federal Drive, Colorado Springs, Colorado 80921, USA  
Contact Information: 1-800-231-0978

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Name of the Manufacturer: Spectranetics Corporation  
Residence: 9965 Federal Drive, Colorado Springs, Colorado 80921, USA  
Production Address: 9965 Federal Drive, Colorado Springs, Colorado 80921, USA  
Contact Information: 1-800-231-0978

**【Agent】**

Name of Agent: Philips (China) Investment Ltd. Co.  
Residence: Building A1, #718 Lingshi Road, JingAn District, Shanghai  
Contact Information: 8008100038

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**【After Sales Service Unit】 Philips (China) Investment Ltd. Co.**

Residence: Building A1, #718 Lingshi Road, JingAn District, Shanghai  
Contact Information: 8008100038

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1. Verma A, Wilkoff BL. Intravascular pacemaker and defibrillator lead extraction: A state-of-the-art review. Heart Rhythm, Vol. 1, No. 6, December 2004; 739-745.
2. Smith HJ, et al. Five-years experience with intravascular lead extraction. PACE 1994; 17:2016-2020.
3. Byrd CL, et al. Lead extraction: indications and techniques. Cardiol Clin 1992; 10:735-748.
4. Byrd CL, et al. Intravascular techniques for extraction of permanent pacemaker leads. J Thorac Cardiovasc Surg 1991;101:989-997.
5. Byrd CL, et al. Intravascular lead extraction using Locking Stylets and sheaths. PACE 1990; 13:1871-1875.

This Page Intentionally Left Blank

This Page Intentionally Left Blank

This Page Intentionally Left Blank

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住所：9965 Federal Drive, Colorado Springs, Colorado 80921, USA

生产地址：9965 Federal Drive, Colorado Springs, Colorado 80921, USA

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©2019 Spectranetics Corporation