AngioSculpt[®]X

Drug-Coated PTCA Scoring Balloon Catheter

**AngioSculpt[®]X Drug-Coated Percutaneous Transluminal Coronary Angioplasty (PTCA) Scoring Balloon Catheter
Rapid Exchange (RX) Delivery System**

THIS PAGE INTENTIONALLY LEFT BLANK

TABLE OF CONTENTS

Pg	Lauguage		Instructions for Use
4	English	English	Instuuctions for Use
9	German	Deutsch	Gebrauchsanweisung
14	French	Français	Mode d'emploi
20	Italian	Italiano	Istruzioni per l'uso
25	Spanish	Español	Instrucciones de Uso
30	Danish	Dansk	Brugsanvisninger
35	Swedish	Svenska	Bruksanvisning
40	Norwegian	Norsk Bokmål	Bruksinstruksjoner
45	Dutch	Nederlands	Gebruiksaanwijzing
50	Portuguese	Português	Instruções de Utilização
56	Greek	Ελληνικά	Οδηγίες Χρήσης
62	Finnish	Suomi	Käyttöohjeet
67	Polish	Polski	Instrukcja Użytkowania

AngioSculpt[®]X

Drug-Coated PTCA
Scoring Balloon Catheter

AngioSculpt[®]X Drug-Coated Percutaneous Transluminal Coronary
Angioplasty (PTCA) Scoring Balloon Catheter
Rapid Exchange (RX) Delivery System

English / English

INSTRUCTIONS FOR USE

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. FAILURE TO OBSERVE ALL WARNINGS AND PRECAUTIONS MAY RESULT IN COMPLICATIONS.

NOTE: These instructions apply to all balloon diameters and lengths.

STERILE: Sterilized with ethylene oxide gas. **Non-pyrogenic. Do not use if the package is open or damaged.**

CONTENTS: One (1) AngioSculpt[®]X Drug-Coated PTCA Scoring Balloon Catheter.

STORAGE: Store in a dry, dark, cool place.

I. DEVICE DESCRIPTION


1. PTCA Catheter Description

AngioSculpt[®]X Drug-Coated PTCA Scoring Balloon Catheter is a standard PTCA catheter with a scoring balloon near the distal tip. The distal end of the catheter has a conventional nylon-blend balloon and a nitinol scoring element with three spiral struts that wrap around the balloon. The struts create focal concentrations of dilating force, which minimize balloon slippage and assist in the luminal expansion of stenotic arteries. The scoring balloon is coated with a specialized formulation that includes the anti-proliferative drug, paclitaxel. The drug-coated scoring balloon is designed to expand to a specified diameter and length at a specified pressure. Conventional radiopaque markers aid in positioning the balloon in the stenosis.

Figure 1 shows the distal section of the catheter with the scoring balloon. The proximal end of the balloon is connected to a balloon inflation channel.

The product is offered on a rapid exchange (RX) delivery platform, and is available in balloon diameters of 2.0 – 3.5 mm in 0.5 mm increments, and in scoring balloon lengths of 10, 15 and 20 mm. The catheter length is approximately 137 cm and is compatible with 0.014-inch guide wires and 6F guide catheters. The catheter is supplied sterile and is intended for a single use.

Figure 1: Distal Section of AngioSculpt[®]X Drug-Coated PTCA Scoring Balloon Catheter

2. Drug Coating Description

The drug coating is a non-polymer based formulation, consisting of paclitaxel as the active pharmaceutical ingredient and the excipient nordihydroguaiaric acid (NDGA). The coating covers the working balloon surface, part of the balloon cones, and the scoring element of the AngioSculpt[®]X catheter at an average surface concentration of 3µg/mm². The key functional characteristic of the formulation is to allow for release of paclitaxel to the tissue of the vascular wall during inflation.

II. INDICATIONS

The AngioSculpt[®]X Drug-Coated PTCA Scoring Balloon is indicated for the treatment of hemodynamically significant coronary artery stenosis, including in-stent restenosis, for the purpose of improving myocardial perfusion.

III. CONTRAINDICATIONS

The AngioSculpt[®]X catheter should not be used for the following:

- Coronary artery lesions unsuitable for treatment by percutaneous revascularization.
- Coronary artery spasm in the absence of a significant stenosis.
- Patients with known hypersensitivity to paclitaxel or paclitaxel-related compounds.
- Patients who cannot receive recommended anti-platelet and/or anticoagulant therapy
- Women who are breastfeeding, pregnant, or are intending to become pregnant or men intending to father children.

IV. WARNINGS

- To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.
- When the catheter is exposed to the vascular system, it should be manipulated while under high quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Do not exceed the rated burst pressure (RBP) during balloon inflation. The RBP is based on results of in-vitro testing. At least 99.9% of the balloons (with 95% confidence) will not burst at or below their RBP. Use of a pressure monitoring device is recommended to prevent over-pressurization.
- PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of a potential cardiovascular injury or life-threatening complication.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate the balloon.
- Proceed cautiously when using the AngioSculpt[®]X catheter in a freshly deployed stent. The AngioSculpt[®]X catheter has not been tested for post-dilatation of stents or lesions distal to freshly deployed stents in clinical studies.
- Use the device prior to the expiration date specified on the package.

V. PRECAUTIONS

- Do not immerse the AngioSculpt[®]X catheter in a saline bath. Contact with fluids prior to insertion may compromise therapeutic drug delivery. In order to ensure therapeutic drug delivery, replace any device where the balloon has come into contact with fluids prior to use.
- The AngioSculpt[®]X catheter should be handled with dry sterile gloves whenever possible prior to use. Care should be taken to minimize contact with the coated balloon portion of the AngioSculpt[®]X catheter during preparation and insertion.
- To ensure therapeutic drug delivery:
 - Never inflate the AngioSculpt[®]X catheter prior to reaching the target lesion.
 - The AngioSculpt[®]X catheter should be advanced to the target site in an efficient manner (i.e. ≤ 3 minutes) and immediately inflated.
- Maintain balloon inflation for a minimum of 30 seconds.
- Predilatation of the lesion with an uncoated standard PTCA balloon should be performed.
- Always advance and retrieve the AngioSculpt[®]X catheter under negative pressure. Whenever possible, the AngioSculpt[®]X catheter should be the final treatment of the vessel.
- The AngioSculpt[®]X catheter is intended to treat a single lesion in a single patient; do not re-use the AngioSculpt[®]X catheter to dilate additional lesions.
- Prior to angioplasty, examine the catheter to verify functionality, catheter integrity and to ensure that its size and length are suitable for the specific lesion for which it is to be used.
- Only physicians trained in the performance of percutaneous transluminal coronary angioplasty should use the AngioSculpt[®]X catheter.

- Appropriate dual antiplatelet, anticoagulant and coronary vasodilator therapy should be administered before, during, and after treatment with the AngioSculpt®X catheter. Antiplatelet therapy for less than 3 months following treatment with the AngioSculpt®X catheter has not been studied and patients in the first-in-human study were administered dual antiplatelet therapy consisting of aspirin plus either clopidogrel or ticlopidine for a minimum of 3 months following treatment with the AngioSculpt®X catheter.
- Do not rotate the catheter shaft in excess of 180 degrees when the tip is constrained.
- Do not rotate the catheter luer hub in excess of five (5) turns during use.
- Do not advance or retract the AngioSculpt®X catheter over the floppy portion of the guide wire.
- Catheter manipulation, including advancement and retraction, should be performed by grasping the hypotube shaft.
- If unusual resistance is felt when the catheter is being manipulated or if it is suspected that the guide wire has become kinked, carefully remove the entire catheter system (AngioSculpt®X catheter and steerable guide wire) as a unit.
- If fluoroscopic guidance indicates that the AngioSculpt®X catheter has advanced beyond the end of the guide wire, withdraw the catheter and reload the wire before advancing again.
- It is not recommended that the AngioSculpt®X catheter be used in conjunction with other drug coated balloons or drug eluting stents to treat the same lesion in the same procedure or within 90 days. The safety of combinations of different drug-device products has not been assessed.
- DO NOT sterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing.
- Reuse of this single use device could lead to serious patient injury or death and voids manufacturer warranties.

VI. DRUG INFORMATION

Mechanism of Action

The AngioSculpt®X Drug-Coated PTCA Scoring Balloon Catheter coating contains paclitaxel, an anti-proliferative pharmaceutical agent that specifically binds to and stabilizes microtubules. By blocking microtubule depolymerization, paclitaxel affects inhibition of smooth muscle cell and fibroblast proliferation and migration as well as secretion of extracellular matrix. The combination of these effects results in the inhibition of neointimal hyperplasia and therefore prevents restenosis.

Pharmacokinetics

The maximum paclitaxel plasma level measured after 10 minutes in the porcine model was 4.6 ng/mL, significantly lower than the myelosuppression level of ≥ 85 ng/mL. After 24 hours, the paclitaxel plasma level dropped to less than 0.4 ng/mL. The excipient NDGA could not be detected in plasma (LOQ: 2 ng/mL). The observed peak levels and total serum exposure of the AngioSculpt®X catheter are therefore roughly 50 to 200-fold lower than that reported for pharmaceutical paclitaxel (3 h. infusion, Taxol™ package insert: T1/2 = 12.1 to 20.2 hours, Cmax = 2170 to 3650 ng/mL, and AUC = 7952 to 15007 ng•h/mL).

After 10 minutes, paclitaxel tissue levels at the target lesion reached 7.1% of the nominal drug dose, or 142 ng/mg. Between 1 and 28 days, the tissue level stabilized at 1.1% of the nominal drug dose, or 15 ng/mg. The mean concentration of paclitaxel in coronary arteries was approximately 10 μ M which is within the range of paclitaxel concentrations (0.1 to 10 μ M/L) known to inhibit the proliferation of human arterial smooth muscle cells (Axel et al., 1997). The excipient NDGA could not be detected in the artery wall.

Drug Interactions

No peer-reviewed drug interaction studies have evaluated the AngioSculpt®X catheter. Therefore, the instructions for use of all concurrent drug therapies should be consulted for interactions with paclitaxel. When using the AngioSculpt®X catheter on a patient who is taking a drug with known paclitaxel interactions or when initiating a drug therapy in a patient who has recently undergone treatment with the AngioSculpt®X catheter, the potential for both local and systemic drug interactions should be considered. Paclitaxel is metabolized by cytochrome P450 isoenzymes CYP2C8 and CYP3A4, and is a substrate of P-glycoprotein. Agents that compete with or inhibit these isoenzymes

may lead to increased levels of paclitaxel. Caution should be exercised when administering paclitaxel with known substrates or inhibitors of the cytochrome P450 isoenzymes CYP2C8 and CYP3A4, especially when no reliable clinical data exists to describe the effects of the drug interactions.

Carcinogenicity, Genotoxicity, and Reproductive Toxicology

No long-term animal studies have been published in which the carcinogenic potential of paclitaxel was assessed. Similarly, there are no adequate and well-controlled studies on the effects of paclitaxel in pregnant women or in men intending to father children. Paclitaxel's mechanism of action involves interfering with cell proliferation through microtubule stabilization, which may result in the loss of chromosomes during cell division. Although paclitaxel was not mutagenic in the Ames test or in the CHO/HGPRT and salmonella assays, this indirect action has been shown to cause DNA fragmentation within in vitro and in vivo micronucleus genotoxicity assays. Paclitaxel has also been reported to cause chromosomal aberrations in primary human lymphocytes.

The treating physician should weigh the potential medical benefits of the AngioSculpt®X catheter against these potential genotoxic and reproductive risks.

VII. ADVERSE EFFECTS

Possible adverse effects include, but are not limited to, the following:

<ul style="list-style-type: none"> • Death • Heart Attack (acute myocardial infarction) • Total occlusion of the treated coronary artery • Coronary artery dissection, perforation, rupture, or injury • Pericardial tamponade • No/slow reflow of treated vessel • Emergency coronary artery bypass (CABG) • Emergency percutaneous coronary intervention • CVA/stroke • Pseudoaneurysm • Restenosis of the dilated vessel • Unstable angina • Thromboembolism or retained device components 	<ul style="list-style-type: none"> • Irregular heart rhythm (arrhythmias, including life-threatening ventricular arrhythmias) • Severe low (hypotension)/high (hypertension) blood pressure • Coronary artery spasm • Hemorrhage or hematoma • Need for blood transfusion • Surgical repair of vascular access site • Creation of a pathway for blood flow between the artery and the vein in the groin (arteriovenous fistula) • Drug reactions, allergic reactions to x-ray dye (contrast medium) • Infection • Allergic reaction to the drug coating
--	---

There may be other potential adverse events that are unforeseen at this time.

VIII. SUMMARY OF CLINICAL STUDY

Multi-Center First-In-Human Clinical Trial STUDY DESIGN

PATENT-C was a first-in-human, multi-center, controlled, randomized, single-blind, prospective clinical investigation comparing novel paclitaxel-coated AngioSculpt scoring balloon with otherwise identical commercially available bare AngioSculpt balloons in patients with significant ($\geq 70\%$) coronary artery bare metal in-stent restenosis. The study evaluated the efficacy, safety and acute tolerance of paclitaxel-coated AngioSculpt balloons in inhibiting bare metal stent (BMS) restenosis.

A total of 61 patients were enrolled at five clinical sites: 4 sites in Germany and 1 site in Brazil. Thirty-three patients were randomized to paclitaxel-coated balloon and 28 patients were randomized to the bare control group. Protocol required follow-up included quantitative coronary angiography (QCA) performed at 6 months following the index procedure and clinical follow-up at 30 days and 6, 12, and 24 months after the index procedure.

A blinded, independent core lab performed QCA analysis of all the films for the interventional procedure, all 6 month follow-up angiographies, target lesion revascularization up to 6 month follow-up, and other unscheduled angiographies or reinterventions up to 6 months. Additionally, a blinded, independent Clinical Events Committee (CEC) adjudicated all deaths; myocardial infarctions; target lesion and target vessel revascularizations; and device failures.

The primary efficacy objective was to demonstrate reduced in-segment late lumen loss (LLL) in the coated-balloon group compared to the bare group at 6 months per QCA analysis. Secondary endpoints included procedural success, major adverse cardiac events (MACE) through 6 months, and the individual clinical endpoints of death, stent thrombosis, myocardial infarction and coronary revascularization through 6 months. Additionally, patients were followed for the same individual clinical endpoints of death, stent thrombosis, myocardial infarction and coronary revascularization through 12 months for the purpose of assessing long-term safety.

Results were analyzed on an intent-to-treat (ITT) and per-protocol (PP) basis. The ITT population was comprised of patient grouping according to randomization assignment regardless of actual treatment received. The PP population was comprised of the patient grouping according to the actual treatment received. All 28 patients (and 30 lesions) randomized to the bare group were successfully treated with a bare AngioSculpt device, and were eligible for inclusion in the ITT and PP analyses. As 3 patients in the coated group did not receive treatment with a coated AngioSculpt, 3 of the 33 patients and lesions in the coated ITT population were excluded from the PP population, leaving 30 patients/lesions in the PP population.

PRIMARY EFFICACY ENDPOINT

Late Lumen Loss at 6 Months

The PATENT-C Trial met the primary endpoint of the study of in-segment LLL (calculated as difference between the angiographic in-segment minimum lumen diameter [MLD] post procedure and at 6 month follow-up) evaluated by QCA. LLL at 6 months in the bare balloon was 0.48±0.51 for both the ITT and PP analysis compared to 0.17±0.40 (p = 0.01) and 0.12±0.26 (p = 0.009) in the coated ITT and PP groups respectively (Table 1).

Table 1: LLL at 6 Month Follow-up by QCA (ITT and PP Population)

Population		Lesions available for analysis	Bare Balloon	Paclitaxel-coated Balloon	p-value
In-segment	ITT	27 bare/	0.48±	0.17±	0.01
LLL, mm		27 coated	0.51*	0.40*	
In-segment	PP	27 bare/	0.48±	0.12±	0.009
LLL, mm		24 coated	0.51*	0.26*	

* Values are mean±SD

SECONDARY ENDPOINTS AND ONGOING SAFETY SURVEILLANCE

Procedural Success

In accordance with the protocol, which defined procedural success as < 50% diameter stenosis (by Core Lab analysis) and absence of in-hospital MACE, the three patients in whom an AngioSculpt could not be crossed/advanced to the lesion were considered device failures, but not procedural failures. Furthermore, as the protocol stated that procedural success rate is calculated as the number of patients achieving procedural success divided by the total number of patients treated with the drug-coated AngioSculpt device, the 3 patients who were not treated with an AngioSculpt device were removed from the calculation of procedural success rate. Success rate was therefore 100% in the drug-coated group and 92.6% in the bare balloon group. The source of failure in the bare group included one patient who experienced an in hospital MACE and another patient with > 50% final residual diameter stenosis by QCA. The success rates are shown in Table 2.

Table 2: Procedural Success Rate

Procedural Success Rate	Bare Balloon	Paclitaxel-coated Balloon	p-value
Success rate, success/total patients treated with AngioSculpt device (%)	25/27† (92.6%)*	30/30 (100%)*	0.2

Values are n (%), * Angiogram unavailable for one patient in the bare group, so cannot be certain of procedural success; therefore N = 27 patients.

MACE up to 6 Months

Adjudicated MACE was defined as cardiac death, target vessel MI, or clinically driven target lesion revascularization (TLR) through 6 months. MACE rates were calculated by counting one event type per patient: once a patient met an endpoint, subsequent events of the same endpoint were not counted for that patient. Within 6 months after the index procedure the cumulative MACE rate was 32.1% in the bare group and

6.1% in the coated group. The rate of target vessel MI (TV-MI) was 7.1% in the bare group and zero in the coated group. The rate of clinically driven TLR was 32.1% in the bare group compared to 3% in the coated group. No patients died of cardiac causes in the bare group and 1 patient (3%) in the coated group experienced cardiac death. The cardiac death in the coated patient was determined to be secondary to non-target vessel MI resulting from non-TVR performed the day prior to death. Table 3 and Table 4 summarize the MACE rates up to 6 months post procedure in the ITT and PP populations.

Table 3: Secondary endpoint of MACE at 6 months of follow-up, ITT Population

Secondary Endpoint Occurrence	Bare Balloon, N= 28 patients	Paclitaxel-Coated Balloon, N = 33 patients	p-value
Any MACE	9 (32.1%)	2 (6.1%)	0.016
Cardiac Death	0	1 (3.0%)	1.0
Target Vessel MI	2* (7.1%)	0	0.21
Clinically driven TLR	9* (32.1%)	1 (3.0%)	0.004

Values are n (%), * One patient in the bare group experienced two TV-MIs and one TLR, but is counted as having reached the endpoints once per usual convention.

Table 4: Secondary endpoint of MACE at 6 months of follow-up, PP population

Secondary Endpoint Occurrence	Bare Balloon, N= 28 patients	Paclitaxel-Coated Balloon, N = 30 patients	p-value
Any MACE	9 (32.1%)	1 (3.3%)	0.005
Cardiac Death	0	1 (3.3%)	1.0
Target Vessel MI	2* (7.1%)	0	0.23
Clinically driven TLR	9* (32.1%)	0	0.0006

Values are n (%), * One patient in the bare group experienced two TV-MIs and one TLR, but is counted as having reached the endpoints once per usual convention.

Individual Clinical Endpoints up to 6 Months

All patients were followed for the individual clinical endpoints of stent thrombosis (occurrence and timing), death (cardiac/non-cardiac), MI (target vessel or any infarction), and repeat revascularization (clinically driven TLR, clinically driven TVR, any revascularization) through 6 months. Individual clinical endpoint rates were calculated by counting one event of each type per patient; once. There were no cases of stent thrombosis in either group. In addition to the cardiac death noted in the preceding section (MACE through 6 months), one patient in the bare group experienced non-cardiac death (ITT and PP analyses). The rates of MI were similar in both groups. The rate of clinically driven TVR was 35.7% in the bare group, including 9 patients (32.1%) that also met the endpoint of clinically driven TLR (ITT and PP). In the bare group, one patient experienced clinically driven TLR (and therefore also TVR) in the ITT analysis. As this patient was not included in the PP population for reason of not receiving a drug-coated AngioSculpt, no patients in the coated group PP population experienced clinically driven TVR or TLR. The individual clinical endpoint rates through 6 months are shown in tables 5 and 6, which summarize events up to 12 months.

Any Death, Any Myocardial Infarction, Coronary Revascularization, and Stent Thrombosis up to 24 months

Although not part of the formal secondary endpoint analysis, all patients were also followed for death, MI, coronary revascularization and stent thrombosis through 24 months as part of ongoing safety surveillance. Between 6 and 12 months of follow up one patient in the coated group died, secondary to TV-MI and possible very late stent thrombosis (ST). The cardiac death in this patient was classified as not related to the study device or procedure per investigator assessment and this was confirmed by the CEC. The CEC made the additional determination of possible very late ST and TV-MI as a conservative measure on the basis that there was not sufficient source documentation to exclude these events. As the possible stent thrombosis occurred > 1 year after the index procedure (1 year and 3 days after), it was classified as very late per the ARC consensus document; however, as it was within the 1 year follow-up window of 12±2 months, it was included in the calculation of 12 month event rates. Table 5 (ITT) and 6 (PP) summarize the event rates up to 24 months of follow up.

Table 5: Adjudicated Death, Myocardial Infarction, Coronary Revascularization, and Stent Thrombosis up to 24 months, ITT population

Clinical Endpoint Event	Bare Balloon, N = 28 patients	(n)*	Paclitaxel-Coated Balloon, N = 33 patients	(n)*	p-value
Events Through 6 Months					
Any Stent Thrombosis	0	28	0	33	-
Any Death	1 (3.6%)	28	1 (3.0%)	33	1.0
Cardiac	0	28	1 (3.0%)	33	1.0
Non-Cardiac	1 (3.6%)	28	0	33	0.46
Any MI	2 (7.1%)	28	1 (3.0%)	33	0.6
Target Vessel MI	2 (7.1%)	28	0	33	0.21
Non-target vessel MI	0	28	1 (3.0%)	33	1.0
Clinically driven TLR	9 (32.1%)	28	1 (3.0%)	33	0.004
Clinically driven TVR†	10 (35.7%)	28	1 (3.0%)	33	0.0016
Any revascularization‡	12 (42.9%)	28	6 (18.2%)	33	0.05
Events Through 12 Months					
Any Stent Thrombosis	0	27	1 (3.2%)	31	1.0
Any Death	1 (3.6%)	28	2 (6.3%)	32	1.0
Cardiac	0	27	2 (6.3%)	32	0.5
Non-Cardiac	1 (3.6%)	28	0	31	0.47
Any MI	2 (7.4%)	27	2 (6.3%)	32	1.0
Target Vessel MI	2 (7.4%)	27	1 (3.2%)	31	0.6
Non-target Vessel MI	0	27	1 (3.1%)	32	1.0
Clinically driven TLR	9 (33.3%)	27	1 (3.2%)	31	0.004
Clinically driven TVR†	10 (37.0%)	27	1 (3.2%)	31	0.0016
Any revascularization‡	12 (44.4%)	27	7 (21.9%)	32	0.094
Events Through 24 Months					
Any Stent Thrombosis	0	27	1 (3.2%)	31	1.0
Any Death	1 (3.6%)	28	3 (9.4%)	32	0.6
Cardiac	0	27	2 (6.3%)	32	0.5
Non-Cardiac	1 (3.6%)	28	1 (3.3%)	30	1.0
Any MI	3 (11.1%)	27	2 (6.3%)	32	0.7
Target Vessel MI	3 (11.1%)	27	1 (3.2%)	31	0.3
Non-target Vessel MI	0	27	1 (3.2%)	31	1.0
Clinically driven TLR	9 (33.3%)	27	1 (3.3%)	30	0.004
Clinically driven TVR†	10 (37.0%)	27	2 (6.7%)	30	0.008
Any revascularization‡	12 (44.4%)	27	7 (22.6%)	31	0.097

Values are n(%).
 * Number of patients with available data/in denominator for calculation. For example, 12 month rates in the bare group, one patient died (cardiac) prior to 6 month follow-up, therefore n = 27 except for any/non-cardiac death, which is n = 28. For the 12 month rates in the coated group, one patient with no 12 month follow-up and another patient with cardiac death, non-TV-MI, and any revascularization; therefore, n = 32 for those parameters and n = 31 for the remaining parameters.
 †TVR includes TVR/non-TLR and TLR
 ‡Any revascularization included patients with TLR, TVR, and non-TV

Table 6: Adjudicated Death, Myocardial Infarction, Coronary Revascularization, and Stent Thrombosis up to 24 months, PP population

Clinical Endpoint Event	Bare Balloon, N = 28 patients	(n)*	Paclitaxel-Coated Balloon, N = 30 patients	(n)*	p-value
Events Through 6 Months					
Any Stent Thrombosis	0	28	0	30	-
Any Death	1 (3.6%)	28	1 (3.3%)	30	1.0
Cardiac	0	28	1 (3.3%)	30	1.0
Non-Cardiac	1 (3.6%)	28	0	30	0.48
Any MI	2 (7.1%)	28	1 (3.3%)	30	0.61
Target Vessel MI	2 (7.1%)	28	0	30	0.23
Non-target vessel MI	0	28	1 (3.3%)	30	1.0
Clinically driven TLR	9 (32.1%)	28	0	30	0.0006
Clinically driven TVR†	10 (35.7%)	28	0	30	0.0003
Any revascularization‡	12 (42.9%)	28	5 (16.7%)	30	0.04
Events Through 12 Months					
Any Stent Thrombosis	0	27	1 (3.6%)	28	1.0
Any Death	1 (3.6%)	28	2 (6.9%)	29	1.0
Cardiac	0	27	2 (6.9%)	29	0.5
Non-Cardiac	1 (3.6%)	28	0	28	1.0
Any MI	2 (7.4%)	27	2 (6.9%)	29	1.0
Target Vessel MI	2 (7.4%)	27	1 (3.6%)	28	0.6
Non-target Vessel MI	0	27	1 (3.5%)	29	1.0
Clinically driven TLR	9 (33.3%)	27	0	28	0.0007
Clinically driven TVR†	10 (37.0%)	27	0	28	0.0003
Any revascularization‡	12 (44.4%)	27	6 (20.7%)	29	0.09
Events Through 24 Months					
Any Stent Thrombosis	0	27	1 (3.6%)	28	1.0
Any Death	1 (3.6%)	28	3 (10.3%)	29	0.6
Cardiac	0	27	2 (6.9%)	29	0.5
Non-Cardiac	1 (3.6%)	28	1 (3.7%)	27	1.0
Any MI	3 (11.1%)	27	2 (6.9%)	29	0.7
Target Vessel MI	3 (11.1%)	27	1 (3.6%)	28	0.4
Non-target Vessel MI	0	27	1 (3.6%)	28	1.0
Clinically driven TLR	9 (33.3%)	27	0	27	0.002
Clinically driven TVR†	10 (37.0%)	27	1 (3.7%)	27	0.005
Any revascularization‡	12 (44.4%)	27	6 (21.4%)	28	0.09

Values are n(%).
 * Number of patients with available data/in denominator for calculation. For example, 12 month rates in the bare group, one patient died (cardiac) prior to 6 month follow-up, therefore n = 27 except for any/non-cardiac death, which is n = 28. For the 12 month rates in the coated group, one patient with no 12 month follow-up and another patient with cardiac death, non-TV-MI, and any revascularization; therefore, n = 29 for those parameters and n = 28 for the remaining parameters.
 †TVR includes TVR/non-TLR and TLR
 ‡Any revascularization included patients with TLR, TVR, and non-TV

MATERIALS REQUIRED FOR USE WITH THE ANGIOSCULPT®X CATHETER

WARNING - Use single use items only. Do not resterilize or reuse.

- Femoral, brachial or radial guiding catheter (≥ 6F)
- Hemostatic valve
- Contrast medium diluted 1:1 with normal saline
- Sterile heparinized normal saline
- 10cc and 20cc syringes for flushing and balloon prep
- Inflation device (indeflator)
- 0.014" coronary guide wire
- Guide wire introducer
- Guide wire torque device
- Radiographic contrast
- Manifold (for pressure monitoring and contrast injection), extension pressure tubing

IX. DIRECTIONS FOR USE

Preparation and Use of the AngioSculpt®X Catheter

Prior to use of the AngioSculpt®X catheter, examine carefully for damage and catheter integrity. Do not use if the catheter has bends, kinks, missing components or other damage. Do not use if inner package is open or damaged.

1. Pre-medicate patient with appropriate dual antiplatelet therapy, anticoagulant, and coronary vasodilator therapy consistent with institutional protocol for percutaneous coronary interventions.
2. Perform coronary angiogram in the view best demonstrating the target lesion prior to device deployment.
3. Position 0.014" coronary guide wire of choice beyond the target lesion.
4. Pre-dilate the lesion with an uncoated standard PTCA balloon sized smaller than the reference vessel.
5. If multiple AngioSculpt®X catheters are required to complete treatment, the sequentially used AngioSculpt®X catheters should be minimally sized and positioned so that the balloons overlap as necessary to cover the lesion and margins of the predilatation segment. The AngioSculpt®X catheter should extend approximately 2 mm proximally and distally of the pre-dilatation segment. Care should be taken not to extend the entire dilatation segment unnecessarily. An example is shown in Figure 2.


Figure 2. Balloons are sized to overlap only as much as needed to treat lesion appropriately

6. Using sterile technique, remove an appropriately sized (≤ 1.0 x reference vessel diameter (RVD)) AngioSculpt®X catheter from the sterile package and place on the sterile field.

Caution: The AngioSculpt®X balloon should be handled with dry sterile gloves whenever possible prior to use. Care should be taken to minimize unnecessary contact with the coated balloon portion of the AngioSculpt®X catheter during preparation and insertion.

7. Remove the stylet from the guidewire lumen and the protective tubing from the balloon.
8. Inspect the catheter to ensure that all components are intact.
9. Flush the guide wire lumen with saline, taking care to keep the balloon dry.
Caution: Do not immerse the AngioSculpt®X catheter in a saline bath. Replace any device if its balloon has come in contact with fluids prior to use.
10. Attach 20 cc syringe filled with 2-3 cc of radiographic contrast to the catheter balloon inflation port.
11. Aspirate/remove air from the catheter balloon lumen using the 20 cc syringe filled with 2-3 cc of radiographic contrast and leave on vacuum for 30 seconds.
12. Gently release vacuum from the 20 cc syringe and remove it from the balloon inflation port.
13. Attach inflation device (indeflator), filled with 50:50 mixture of radiographic contrast and normal saline, to the balloon inflation port by creating a meniscus. Avoid introducing air bubbles into the catheter balloon lumen.
14. Aspirate using the inflation device, locking in vacuum.

NOTE: All air must be removed from the balloon and displaced with contrast medium prior to inserting into the body (repeat steps 11-14 if necessary).

Caution: To ensure therapeutic delivery, never inflate the AngioSculpt®X catheter prior to reaching the target lesion.

15. Advance the AngioSculpt®X catheter over the coronary guide wire to the target lesion.

NOTE: When back loading the catheter onto the guide wire, the catheter should be supported, ensuring that the guide wire does not come in contact with the balloon. Do not advance or retract the AngioSculpt®X catheter over the floppy portion of the guide wire. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.

Caution: To ensure appropriate drug delivery, the AngioSculpt®X catheter should be advanced to the target site in an efficient manner and immediately inflated. Maintain balloon inflation for a minimum of 30 seconds.

16. Position the balloon relative to the lesion, ensuring coverage at least 2 mm proximally and distally beyond the margins of the predilatation lesion segment, and inflate the balloon to the appropriate pressure (reference Compliance Chart included in the product packaging).
17. Inflate the AngioSculpt®X balloon per the following recommended protocol:
 - 2 atmospheres
 - Increase the inflation pressure by 2 atmospheres every 10-15 seconds until full device inflation is achieved
 - May inflate to a maximum pressure that is $<$ RBP at the discretion of the physician (bearing in mind the estimated inflated diameter of the device at a given pressure)
18. Perform coronary angiogram (in the same view(s) as step 2) of the target lesion following completion of device treatment.
19. To remove the AngioSculpt®X catheter, apply negative pressure to the inflation device and confirm that the balloon is fully deflated. The catheter should be retracted only by grasping the hypotube shaft.
NOTE: Do not rotate the catheter shaft in excess of 180 degrees when the tip is constrained. Do not rotate the catheter luer hub in excess of five (5) turns during use. Catheter manipulation, including advancement and retraction, should be performed by grasping the catheter shaft.
20. Inspect all components to ensure that the catheter is intact. Follow institutional procedures for disposal of biohazards. If device malfunction occurs or any defects are noted on inspection, flush the guide wire lumen and clean the outer surface of the catheter with saline, store the catheter in a sealed biohazard bag, and contact Spectranetics or the EC Authorized Representative for further instructions.
NOTE: Whenever possible, the AngioSculpt®X catheter should be the final treatment of the vessel.
21. Remove the coronary guide wire and perform coronary angiography (in the same view(s) as step 2) of the target lesion following completion of all interventions.
22. Remove all catheters and manage the arterial access site according to institutional protocol.

X. MANUFACTURER'S LIMITED WARRANTY

Manufacturer warrants that the AngioSculpt®X Drug-Coated PTCA Scoring Balloon Catheter is free from defects in material and workmanship when used by the stated "Use By" date and when package is unopened and undamaged immediately before use. Manufacturer's liability under this warranty is limited to replacement or refund of the purchase price of any defective AngioSculpt®X Drug-Coated PTCA. Manufacturer will not be liable for any incidental, special, or consequential damages resulting from use of the AngioSculpt®X Drug-Coated PTCA. Damage to the AngioSculpt®X Drug-Coated PTCA caused by misuse, alteration, improper storage or handling, or any other failure to follow these Instructions for Use will void this limited warranty. **THIS LIMITED WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** No person or entity, including any authorized representative or reseller of Manufacturer, has the authority to extend or expand this limited warranty and any purported attempt to do so will not be enforceable against Manufacturer.

Medikamentenbeschichteter AngioSculpt®X Perkutane transluminale koronare Angioplastie (PTCA) Scoring-Ballonkatheter Rapid Exchange (RX) Zuführsystem

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No long-term animal studies have been published in which the carcinogenic potential of paclitaxel was assessed. Similarly, there are no adequate and well-controlled studies on the effects of paclitaxel in pregnant women or in men intending to father children. Paclitaxel's mechanism of action involves interfering with cell proliferation through microtubule stabilization, which may result in the loss of chromosomes during cell division. Although paclitaxel was not mutagenic in the Ames test or in the CHO/HGPRT and salmonella assays, this indirect action has been shown to cause DNA fragmentation within *in vitro* and *in vivo* micronucleus genotoxicity assays. Paclitaxel has also been reported to cause chromosomal aberrations in primary human lymphocytes. The treating physician should weigh the potential medical benefits of the AngioSculpt®X catheter against these potential genotoxic and reproductive risks.

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© 2019 Spectranetics Corporation
Spectranetics International B.V.
The Spectranetics Corporation