The logo features a thick, black, curved line above the text. The text "Quick-Cross" is in a bold, italicized sans-serif font, with a trademark symbol (TM) to the right. Below it, the word "SELECT" is in a larger, bold, italicized sans-serif font. The "S" in "SELECT" is connected to the "C" in "Quick-Cross" by a thin, curved line.

Quick-CrossTM
SELECT

Support Catheter

Instructions for Use

The logo consists of a stylized diamond shape with a white circle inside, followed by the word "Spectranetics" in a bold, italicized sans-serif font, with a registered trademark symbol (®) to the right.

 **Spectranetics[®]**

CE
2797

Instructions for Use - Sections by Language

Table of Contents

| Pg | Language | | Instructions For Use |
|----|------------|-----------------|---------------------------|
| 3 | English | English | Instructions for Use |
| 5 | Bulgarian | български език | Инструкции за употреба |
| 7 | Croatian | Hrvatski | Upute za uporabu |
| 9 | Czech | Česky | Návod k Použití |
| 11 | Danish | Dansk | Brugsanvisning |
| 13 | Dutch | Nederlands | Gebruiksaanwijzing |
| 15 | Estonian | Eesti | Kasutusjuhend |
| 17 | Finnish | Suomi | Käyttöohjeet |
| 19 | French | Français | Mode d'emploi |
| 21 | German | Deutsch | Gebrauchsanweisung |
| 23 | Greek | Ελληνικά | Οδηγίες χρήσης |
| 25 | Hungarian | Magyar | Használati utasítások |
| 27 | Icelandic | Íslenska | Leiðbeiningar um notkun |
| 29 | Italian | Italiano | Istruzioni per l'uso |
| 31 | Latvian | Latviešu valoda | Lietošanas instrukcija |
| 33 | Lithuanian | lietuvių kalba | Naudojimo instrukcija |
| 35 | Norwegian | Norsk bokmål | Bruksinstruksjoner |
| 37 | Polish | Polski | Instrukcja użytkowania |
| 39 | Portuguese | Português | Instruções de Utilização |
| 41 | Romanian | Română | Instrucțiuni de utilizare |
| 43 | Russian | Русский язык | Инструкции по применению |
| 45 | Serbian | Srpski jezik | Uputstvo za upotrebu |
| 47 | Slovak | Slovenčina | Návod na použitie |
| 49 | Slovenian | Slovenščina | Navodila za Uporabo |
| 51 | Spanish | Español | Instrucciones de Uso |
| 53 | Swedish | Svenska | Brugsanvisning |
| 55 | Turkish | Türkçe | Kullanım Talimatları |

Table of Contents

| | |
|--------------------------------------------|---|
| 1. DESCRIPTION | 3 |
| 2. INDICATIONS FOR USE | 3 |
| 3. CONTRAINDICATIONS | 3 |
| 4. WARNINGS | 3 |
| 5. PRECAUTIONS | 3 |
| 6. ADVERSE EFFECTS | 3 |
| 7. HOW SUPPLIED | 3 |
| 8. COMPATIBILITY | 4 |
| 9. PROCEDURE SET-UP AND DIRECTIONS FOR USE | 4 |
| 10. MANUFACTURER'S LIMITED WARRANTY | 4 |
| 11. NON-STANDARD SYMBOLS | 4 |

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician

1. DESCRIPTION

The Spectranetics Quick-Cross Select Support Catheters are intravascular catheters. These catheters are available in a variety of lengths and tip configurations. All models have 3 radiopaque markers spaced equally along the distal shaft to aid in estimating geometry within the vascular system. The distal radiopaque marker is positioned within 3 mm of the distal catheter tip. A standard female luer is placed on the proximal end of each model. The catheter is coated with a lubricious, hydrophilic coating.

2. INDICATIONS FOR USE

Quick-Cross Select Support Catheters are intended to guide and support a guidewire during access of the coronary or peripheral vasculature, allow for wire exchanges and provide a conduit for the delivery of saline solutions or diagnostic contrast agents.

3. CONTRAINDICATIONS

None known.

4. WARNINGS

- Maximum Infusion Pressure: **300 psi** for 0.014" & 0.018" catheters and **500 psi** for 0.035" catheters.
- The catheter is designed and intended for one time use only. **Do not re-sterilize and/or reuse.**
- The catheter should only be used by physicians qualified to perform percutaneous vascular interventions.

5. PRECAUTIONS

- DO NOT re-sterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing. Reuse of this single use device could lead to serious patient injury or death and voids manufacturer warranties.
- Store in a cool, dry place. Protect from direct sunlight and high temperature (greater than 55°C or 131°F).
- Do not use if device or packaging is damaged.
- Use the catheter prior to the "Use By" date (Expiration Date) specified on the package.

6. ADVERSE EFFECTS

Vascular catheterization and/or vascular intervention may result in complications including but not limited to:

- Vessel dissection, perforation, rupture or total occlusion
- Infection
- Hematoma
- Unstable angina
- Embolism
- Hypo/hypertension
- Acute myocardial infarction
- Arrhythmia, including ventricular fibrillation
- Death

7. HOW SUPPLIED

The Spectranetics Quick-Cross Select Support Catheters are supplied **STERILE**. The devices are designated and designed for **SINGLE USE ONLY** and must not be re-sterilized and/or reused.

7.1 Sterilization

- The sterility of the product is guaranteed only if the package is unopened and undamaged. Before use, visually inspect the sterile package to ensure that the seals have not been broken. Do not use the catheter if the integrity of the package has been compromised. Do not use catheter if its "Use Before Date," found on package labeling, has been passed.

7.2 Inspection prior to use

- Before use, examine carefully for defects, all of the equipment to be used. Do not use any equipment if it is damaged.

8. COMPATIBILITY

Product Features & Model Numbers Table

| Model (Ref) | Working Length | Tip Shape / Type | Hydrophilic Coating Length | GW Compatibility | Min. Guide I.D. | Min. Sheath I.D. | Maximum Catheter O.D. | Marker Band Spacing | Throw Angle | |
|-------------|----------------|------------------|----------------------------|------------------|-----------------|------------------|-----------------------|---------------------|-------------|---------|
| 518-085 | 135 cm | Angled | 100 cm | 0.014" | 5F | 4F | 3.2F | 15 mm | 45° | |
| 518-087 | 150 cm | Angled | | 0.36 mm | | | 0.042" | | | 1.07 mm |
| 518-089 | 90 cm | Angled | 60 cm | 0.018" | 5F | 4F | 3.4F | 15 mm | 45° | |
| 518-091 | 135 cm | Angled | 100 cm | 0.46 mm | | | 0.044" | | | 1.12 mm |
| 518-093 | 150 cm | Angled | | | | | | | | |
| 518-077 | 65 cm | Angled | 45 cm | 0.035" | N/A | 5F | 4.5F | 50 mm | 45° | |
| 518-079 | 90 cm | Angled | 70 cm | | | | | | | |
| 518-081 | 135 cm | Angled | 115 cm | | | | | | | |
| 518-083 | 150 cm | Angled | 130 cm | | | | | | | |
| | | | | | | | 0.89 mm | | | 1.50 mm |

9. PROCEDURE SET-UP AND DIRECTIONS FOR USE

Note: Follow instructions for use for all equipment to be used with the Quick-Cross Select Support Catheter.

- Preparation: Using sterile technique, open the sterile package. Gently remove the protective hoop with the catheter from the pouch. Fill a sterile standard luer-lock syringe with sterile saline. Before removing the catheter from the hoop, connect the syringe to the catheter proximal luer fitting, flush the catheter and allow the saline to fill the hoop. Set catheter in hoop aside until ready for use.
- Insertion: Through a previously inserted, appropriately sized guiding catheter or introducer sheath, introduce the catheter over an appropriately sized guidewire (see specifications) using standard technique.
- Advancement: Use fluoroscopic guidance when advancing the catheter to the desired location within the vasculature.
- Removal: Gently withdraw the catheter using standard technique, being careful to maintain guidewire position if the guidewire is to remain in place.

Infusion: To perform infusion, withdraw the guidewire and reference the specifications for maximum infusion pressure.

Note: Do not exceed the maximum infusion pressures.

After use, dispose of all equipment in accordance with applicable requirements relating to hospital waste, and potentially bio-hazardous materials.

10. MANUFACTURER'S LIMITED WARRANTY

Manufacturer warrants that the Quick-Cross Select is free from defects in material and workmanship when used by the stated "Use By" date and when package is unopened and undamaged immediately before use. Manufacturer's liability under this warranty is limited to replacement or refund of the purchase price of any defective Quick-Cross Select. Manufacturer will not be liable for any incidental, special, or consequential damages resulting from use of the Quick-Cross Select. Damage to the Quick-Cross Select caused by misuse, alteration, improper storage or handling, or any other failure to follow these Instructions for Use will void this limited warranty. **THIS LIMITED WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** No person or entity, including any authorized representative or reseller of Manufacturer, has the authority to extend or expand this limited warranty and any purported attempt to do so will not be enforceable against Manufacturer.

11. NON-STANDARD SYMBOLS

| | |
|---------------------------------------------------------------------------------------------|----------------|
| Distal Marker Spacing | |
| Shape | |
| Importer | |
| GW Compatibility | |
| Maximum OD | |
| Maximum PSI/kPa | |
| Sheath Compatibility | |
| Working Length | |
| Quantity | QTY |
| CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. | Rx ONLY |

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©2020 Spectranetics Corporation
Spectranetics Corporation
Spectranetics International B.V.