# stellarex™

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
### DEVICE DESCRIPTION

#### PTA Catheter

The Stellarex™ 0.014" OTW drug-coated angioplasty balloon (Stellarex balloon) consists of an over-the-wire (OTW) coaxial lumen catheter with a distally mounted semi-compliant inflatable balloon and an atraumatic tapered tip. The balloon is coated with a proprietary coating containing the drug paclitaxel.

The catheter is compatible with a 0.014" (0.36 mm) guidewire. Each device has a protective sheath over the drug-coated balloon portion of the catheter. A compliance chart is included on the product label for each device.

The balloon has two radiopaque markers for positioning the balloon relative to the treatment area. In addition, the 150mm length balloon has a second pair of markers to define the middle of the balloon length. The radiopaque marker bands indicate the working length of the balloon and facilitate fluoroscopic visualization during delivery and placement. The drug-coating on the balloon is a formulation consisting of paclitaxel as the active pharmaceutical agent. The paclitaxel coating covers the working length of the balloon body.

#### Drug Coating

The drug coating consists of the active pharmaceutical ingredient paclitaxel and excipients. The drug coating covers the working length of the balloon component of the catheter. The drug coating is evenly distributed across the balloon surface at a concentration of 2µg/mm<sup>2</sup>. The key functional characteristic of the drug coating is to allow for release of paclitaxel to the tissue of the vascular wall during balloon inflation.

### INDICATIONS FOR USE

The Stellarex 0.014" OTW drug-coated angioplasty balloon is indicated for the treatment of de-novo or re-stenotic lesions up to 270mm in length (single-vessel) or 320mm in length (multi-vessel) in native popliteal and infra-popliteal arteries to establish blood flow and to maintain vessel patency.

### CONTRAINDICATIONS

The Stellarex balloon is contraindicated for use in:

- Patients with known hypersensitivity to paclitaxel or structurally related compounds.
- Patients who cannot receive recommended anti-platelet and/or anti-coagulation therapy.
- Women who are breastfeeding, pregnant or are intending to become pregnant or men intending to father children.

### WARNINGS

- **A signal for increased risk of late mortality has been identified by a meta-analysis of randomized controlled trials published in December 2018 by Katsanos et. al. following the use of paclitaxel-coated balloons and paclitaxel-eluting stents for femoropopliteal arterial disease beginning approximately 2-3 years post-treatment compared with the use of non-drug coated devices. There is uncertainty regarding the magnitude and mechanism for the increased late mortality risk, including the impact of repeat paclitaxel-coated device exposure. Physicians should discuss this late mortality signal and the benefits and risks of available treatment options with their patients.**
- The safety of utilizing a combination of Stellarex balloons with a total drug dose exceeding 9.4mg paclitaxel has not been evaluated and should not be used.
- The Stellarex balloon is supplied STERILE for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risk of patient infection and risk of compromised device performance.
- The Stellarex balloon should be used only by physicians who are experienced and knowledgeable of the clinical and technical aspects of percutaneous transluminal angioplasty.
- Prior to use of the Stellarex balloon, physicians should read and understand the instructions for use. Failure to follow the indications, contraindications, restrictions, warnings and precautions may result in complications.
- Do not use after the "Use By" date.
- The Stellarex balloon contains paclitaxel, a known genotoxin. Do not use the Stellarex balloon in women who are breastfeeding, pregnant or intending to become pregnant, or in men intending to father children.
- Never use air or any gaseous medium to inflate the Stellarex balloon.
- When the Stellarex balloon is exposed to the vascular system, it should be manipulated under high quality fluoroscopic observation.
- Do not manipulate the Stellarex balloon in an inflated state.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Resistance may cause damage to device or lumen. Carefully withdraw the catheter.
- Precautions to prevent or reduce clotting should be considered. Physician experience and discretion will determine the appropriate anticoagulation therapy for each patient.
- The Stellarex balloon should be used with caution for procedures involving calcified lesions.
- Allergic reactions to contrast medium should be identified before PTA angioplasty.
- The outer foil pouch is not a sterile barrier. The inner Tyvek pouch is the product sterile barrier. Do not allow the Tyvek inner pouch to contact the sterile field.

### PRECAUTIONS

- The Stellarex balloon should not be inflated in excess of the rated burst pressure (RBP).
- Do not use contrast media that is contraindicated for intravascular use with this device.
- Carefully inspect the Stellarex balloon and package prior to use. Do not use the catheter if it is damaged or if the size, shape or condition is unsuitable for the intended procedure.
- Do not immerse or wipe the balloon section of the Stellarex balloon with any fluid as the integrity of the drug coating may be damaged or compromised. Replace any Stellarex balloon where the balloon has come into contact with fluids prior to use.
- Use sterile gloves to handle the Stellarex balloon prior to use. Care should be taken to minimize contact with the coated balloon portion of the device.
- Pre-dilatation of the target lesion is recommended for highly stenosed and difficult to cross lesions. Pre-dilate using a PTA catheter undersized by at least 0.5 mm with respect to the reference vessel diameter. If positioning of the Stellarex balloon is difficult while attempting to cross the lesion, remove the catheter and attempt a second pre-dilatation.
- Avoid saline solution contact with the Stellarex balloon coating when flushing the wire lumen.
- Never inflate the Stellarex balloon outside the body or prior to reaching the target lesion as it may disrupt the coating integrity.

- Do not attempt to pass the Stellarex balloon through a smaller French size guide catheter or introducer sheath than indicated on the label. Refer to package label for guide catheter and introducer sheath compatibility.
- For proper drug delivery to the target lesion, maintain inflation of the Stellarex balloon for a minimum of 60 seconds. In order to optimize lesion dilatation, longer inflation times may be performed at the discretion of the operator.
- Use of the Stellarex balloon has not been studied in conjunction with other interventional techniques. If provisional (bail out) stenting is required, a bare metal stent indicated for treatment of the lower extremities should be used.
- The optimal duration of antiplatelet therapy for each patient is at the discretion of the physician.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

### LATE MORTALITY SIGNAL FOR PACLITAXEL COATED DEVICES

A meta-analysis of randomized controlled trials published in December 2018 by Katsanos et. al. identified an increased risk of late mortality at 2 years and beyond for paclitaxel-coated balloons and paclitaxel-eluting stents used to treat femoropopliteal arterial disease. In response to these data, FDA performed a patient-level meta-analysis of long-term follow-up data from the pivotal premarket randomized trials of paclitaxel-coated devices used to treat femoropopliteal disease using available clinical data through May 2019. The meta-analysis also showed a late mortality signal in study subjects treated with paclitaxel-coated devices compared to patients treated with uncoated devices. Specifically, in the 3 randomized trials with a total of 1090 patients and available 5-year data, the crude mortality rate was 19.8% (range 15.9% - 23.4%) in patients treated with paclitaxel-coated devices compared to 12.7% (range 11.2% - 14.0%) in subjects treated with uncoated devices. The relative risk for increased mortality at 5 years was 1.57 (95% confidence interval 1.16 - 2.13), which corresponds to a 57% relative increase in mortality in patients treated with paclitaxel-coated devices. As presented at the June 2019 FDA Advisory Committee Meeting, an independent meta-analysis of similar patient-level data provided by VIVA Physicians, a vascular medicine organization, reported similar findings with a hazard ratio of 1.38 (95% confidence interval 1.06 - 1.80). Additional analyses have been conducted and are underway that are specifically designed to assess the relationship of mortality to paclitaxel-coated devices.

The presence and magnitude of the late mortality risk should be interpreted with caution because of multiple limitations in the available data, including wide confidence intervals due to a small sample size, pooling of studies of different paclitaxel-coated devices that were not intended to be combined, substantial amounts of missing study data, no clear evidence of a paclitaxel dose effect on mortality, and no identified pathophysiologic mechanism for the late deaths.

Paclitaxel-coated balloons and stents improve blood flow to the legs and decrease the likelihood of repeat procedures to reopen blocked blood vessels compared to uncoated devices. The benefits of paclitaxel-coated devices (e.g., reduced reinterventions) should be considered in individual patients along with potential risks (e.g., late mortality).

### USE IN SPECIAL POPULATIONS

The safety and effectiveness of the Stellarex balloon has not been established in pediatric patients (< 18 years of age).

Use of the Stellarex balloon in patients ≥ 18 years of age and older is at the discretion of the physician.

### DRUG INFORMATION

#### Mechanism of Action

The Stellarex balloon catheter coating contains paclitaxel, an anti-proliferative pharmaceutical that specifically binds to and stabilizes microtubules, and the excipient polyethylene glycol. Paclitaxel affects inhibition of smooth muscle cell and fibroblast proliferation/migration as well as secretion of extracellular matrix by blocking microtubule proliferation. The combination of these effects results in the inhibition of neointimal hyperplasia and therefore restenosis.

#### Drug Interactions

Formal drug interaction studies have not been conducted for the Stellarex balloon. The respective instructions for use for all drugs used in conjunction with the Stellarex balloon should be consulted for interactions with paclitaxel. Consideration should be given to the potential for systemic and local drug interactions in the vessel wall in a patient who is taking a drug with known interactions to paclitaxel or when deciding to initiate drug therapy in a patient who has been treated with the Stellarex balloon.

The metabolism of paclitaxel is catalyzed by cytochrome P450 isoenzymes CYP2C8 and CYP3A4 and it is a substrate of P-glycoprotein. Potential drug interactions may occur with any drug that affects these isoenzymes. In the absence of formal drug interaction studies, caution should be exercised when administering paclitaxel.

#### Carcinogenicity, Genotoxicity and Reproductive Toxicology

No long-term studies have been performed to evaluate the carcinogenic potential of the Stellarex balloon.

**WARNING: The Stellarex balloon contains paclitaxel, a known genotoxin. Do not use the Stellarex balloon in women who are breastfeeding, pregnant or intending to become pregnant, or in men intending to father children.**

### POTENTIAL COMPLICATIONS / ADVERSE EVENTS

#### Adverse Events

Potential complications which may be associated with a peripheral balloon dilation procedure include, but may not be limited to, the following:

- |   |   |
|---|---|
| • Abnormal heart rhythms  | • Infection or pain at insertion site               |
| • Allergic reaction to contrast medium, antiplatelet therapy, or catheter system components (drug, excipients, and materials) | • Inflammation                                      |
| • Amputation/loss of limb   | • Ischemia  |
| • Aneurysm  | • Occlusion   |
| • Arrhythmias   | • Pain or tenderness                                |
| • Arterio-venous fistula (AVF)  | • Pseudoaneurysm                                    |
| • Bleeding  | • Renal failure                                     |
| • Death   | • Restenosis  |
| • Embolism/device embolism  | • Sepsis/infection                                  |
| • Fever   | • Shock   |
| • Hematoma  | • Stroke/CVA  |
| • Hemorrhage  | • Thrombosis  |
| • Hypertension/hypotension  | • Vessel dissection, perforation, rupture, or spasm |

Potential complications which may be associated with the addition of paclitaxel to a PTA balloon catheter include, but may not be limited to, the following:

- Allergic/immunologic reaction to paclitaxel
- Alopecia
- Anemia
- Edema (non-pulmonary)
- Gastrointestinal symptoms (Diarrhea, nausea, pain, vomiting)
- Hemolysis
- Hematologic dyscrasia (including neutropenia, leucopenia, thrombocytopenia)
- Histologic changes in vessel wall including inflammation, cellular damage, or necrosis
- Myalgia/arthralgia

- Peripheral neuropathy
- Rash
- Transfusion

## PATIENT COUNSELING INFORMATION

Physicians should advise patients on the following:

- Risks associated with a PTA procedure
- Risks associated with a paclitaxel coated PTA catheter
- Pre- and post-procedure care including antiplatelet therapy

## HOW SUPPLIED

The Stellarex balloon is supplied STERILE for single use only (ethylene oxide sterilization). The Stellarex balloon is contained within an inner Tyvek pouch within an outer foil pouch. The pouches are contained within a single unit box.

**WARNING: The outer foil pouch is not a sterile barrier. The inner Tyvek pouch is the product sterile barrier. Do not allow the Tyvek inner pouch to contact the sterile field.**

**WARNING: The Stellarex balloon is supplied STERILE for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risk of patient infection and risk of compromised device performance.**

## STORAGE

The Stellarex balloon should be stored at room temperature in a dry location in its original packaging. The device should be used prior to the "Use By" date printed on the device packaging.

## RECOMMENDED ITEMS

Prepare the following items using sterile technique:

- 10 cc syringe filled with sterile heparinized saline
- Three-way stopcock
- Contrast media - the standard inflation medium is a 1:1 mixture of contrast medium and sterile saline.

**CAUTION: Do not use contrast media that is contraindicated for intravascular use.**

- Appropriately sized exchange guidewire (refer to product labeling)
- Appropriately sized hemostatic introducer sheath (refer to product labeling)
- Inflation device with manometer

## INSPECTION PROCEDURES

Inspect the Stellarex balloon and packaging. Do not use if packaging or product damage is evident. This product may be a potential biohazard, handle and dispose of it in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Inspect the Stellarex balloon "Use By" date. Use before the "Use By" date.


**CAUTION: Carefully inspect the Stellarex balloon prior to use. Do not use the catheter if it is damaged or if the size, shape or condition is unsuitable for the intended procedure.**

## USE OF MULTIPLE STELLAREX BALLOONS

**WARNING: The safety of utilizing a combination of Stellarex balloons with a total drug dose exceeding 9.4mg paclitaxel has not been evaluated and should not be used.**

Balloon Diameter (mm)	Total Nominal Dose per Balloon Size(mg)		
	Balloon Length (mm)		
	40	80	150
2.0	0.5	1.0	2.0
2.5	0.7	1.3	2.5
3.0	0.8	1.6	3.0
3.5	0.9	1.8	3.4
4.0	1.0	2.1	3.9

If multiple Stellarex balloons are required to treat a lesion, the sequentially used Stellarex balloons should be angiographically positioned so that the marker bands of consecutively placed balloons overlap a minimum of 10 mm and the most proximal and most distal balloons extend 10 mm beyond the predilated segment. The use of an arterial land marking system (eg, radiopaque ruler) must be used to ensure appropriate placement of the Stellarex balloons.

## PRE-DILATATION

**CAUTION: Pre-dilatation of the target lesion is recommended for highly stenosed and difficult to cross lesions. Pre-dilate using a PTA catheter undersized by at least 0.5 mm with respect to the reference vessel diameter. If positioning of the Stellarex balloon is difficult while attempting to cross the lesion, remove the catheter and attempt a second pre-dilatation.**

Target lesion pre-dilatation using a PTA balloon is recommended for tight stenoses or difficult to cross lesions. Pre-dilate using a PTA catheter undersized by at least 0.5 mm with respect to the reference vessel diameter. Pre-dilatation assists in subsequent tracking and ensures complete vessel contact of the Stellarex balloon during inflation. Limit the pre-dilatation balloon length to avoid vessel injury outside the vessel area intended for treatment with the Stellarex balloon.

## CATHETER PREPARATION

**WARNING: The outer foil pouch is not a sterile barrier. The inner Tyvek pouch is the product sterile barrier. Do not allow the Tyvek inner pouch to contact the sterile field.**

**CAUTION: Use sterile gloves to handle the Stellarex balloon prior to use. Care should be taken to minimize contact with the coated balloon portion of the device.**

1. Select the appropriate size balloon for the procedure.
2. Remove inner Tyvek pouch from the outer foil pouch and carton outside of the sterile field.
3. Remove the catheter hoop from the Tyvek inner pouch.
4. Carefully remove the catheter from the hoop.
5. Apply negative pressure to the balloon prior to removing the balloon protector sheath and lumen stylet.
6. Grasp the proximal end of the balloon protector with one hand. Using the opposite hand gently slide the thumb and forefinger toward the distal end of the balloon protector to straighten the balloon.
7. Grasp the balloon protector at roughly the midpoint of the balloon protector and pull it away from the balloon catheter. The balloon protector and wire lumen stylet should be removed together.
8. Flush the guidewire lumen with heparinized saline solution through the wire lumen marked "THRU."  
**CAUTION: Avoid saline solution contact with the Stellarex balloon coating when flushing the wire lumen.**
9. Fill a 10 cc syringe with approximately 4 cc of equal volume (1:1) of contrast media and saline.
10. Evacuate air from the balloon and balloon lumen:
  - a. Attach the syringe to the balloon lumen, marked "BALLOON"
  - b. Apply negative pressure and aspirate for 15 seconds. Slowly release the pressure to neutral, allowing contrast media to fill the shaft of the catheter.
  - c. Disconnect the syringe from the "BALLOON" port of the catheter.
  - d. Remove all air from the syringe. Reconnect the syringe to the "BALLOON" port.
  - e. Apply negative pressure on the balloon until air no longer returns to the device.
  - f. Slowly release the device pressure to neutral.
  - g. Repeat as necessary to remove all air from the balloon and lumen.
11. Replace the syringe with an inflation device with manometer, taking care not to introduce air into the catheter.

**CAUTION: Do not immerse or wipe the balloon section of the Stellarex balloon with any fluid as the integrity of the drug coating may be damaged or compromised. Replace any Stellarex balloon where the balloon has come into contact with fluids prior to use.**

## CATHETER INSERTION AND DILATATION

The Stellarex balloon can be introduced percutaneously through an appropriate sized introducer sheath.

**CAUTION: Do not attempt to pass the Stellarex balloon through a smaller French size guide catheter or introducer sheath than indicated on the label. Refer to package label for guide catheter and introducer sheath compatibility.**

1. Apply negative pressure to the balloon.
2. Place the prepared catheter over a pre-positioned guidewire, which has been placed through the lesion, and introduce the catheter percutaneously. Negative pressure should be maintained during advancement over the guidewire.
3. Advance the catheter tip to the treatment location. A suitable length guidewire should be used at all times to maintain control and position of the guidewire.

**CAUTION: Use fluoroscopic guidance to manipulate the Stellarex balloon catheter during the procedure.**

**WARNING: If resistance is encountered at any time during the insertion procedure, do not force passage. Resistance may cause damage to the device or lumen. Carefully withdraw the catheter.**

4. Position the catheter in the treatment location with the working area of the balloon overlapping the end of the treatment location by at least 10mm. The radiopaque marker bands indicate the working length of the balloon. The position of the balloon catheter may only be changed with the guidewire in place.
5. Inflate the balloon to dilate the target area according to the compliance chart printed on the device packaging. Do not exceed rated burst pressure.  
**CAUTION: Assure treatment of the entire target lesion area with Stellarex balloon(s). Always manipulate the Stellarex balloon under fluoroscopic observation when in the body.**  
**CAUTION: For proper drug delivery to the target lesion, maintain inflation of the Stellarex balloon for a minimum of 60 seconds. In order to optimize lesion dilatation, longer inflation times may be performed at the discretion of the operator.**
6. Deflate the balloon and apply negative pressure.
7. With the guidewire in place and with negative pressure in the balloon, withdraw the catheter. Do not retract the catheter unless the balloon is free and fully deflated.
8. If more than one Stellarex balloon is required to treat a single lesion, the working length of the balloons must overlap by at least 10 mm. A new, unused balloon must be used for each deployment.
9. Results should be verified by angiography.
10. If a Stellarex balloon has entered the vasculature and cannot be deployed, the balloon CANNOT be re-inserted for deployment.

## POST-TREATMENT DILATATION OR STENTING

If required, post-treatment balloon dilatation is allowed.

**CAUTION: If provisional (bail out) stenting is required, a bare metal stent indicated for treatment of the lower extremities should be used.**

## DISPOSAL

**CAUTION: After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.**

## WARRANTY

Although this product has been manufactured under carefully controlled conditions, Spectranetics® has no control over the conditions under which this product is used. Spectranetics® therefore disclaims all warranties, both express and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Spectranetics® shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind Spectranetics® to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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©2020 Spectranetics
**Spectranetics Corporation**
**Spectranetics International B.V.**