QuickCat[™]

Extraction Catheter

Instructions for Use

Spectranetics[®]

Instructions for Use - Sections by Language

Pg	Language		Instructions For Use
3	English	English	Instructions for Use
7	Chinese	中文	使用说明

Table of Contents

1.	DEVICE DESCRIPTION	3
2.	DEVICE COMPONENTS	3
3.	INDICATIONS AND INTENDED USE	3
4.	CONTRAINDICATIONS	3
5.	WARNINGS	3
6.	PRECAUTIONS	4
7.	POTENTIAL ADVERSE REACTIONS	4
8.	DEVICE PREPARATION	4
9.	PROCEDURE.....	4
10.	POST-OPERATIVE CARE.....	5
11.	HOW SUPPLIED	5
12.	STORAGE CONDITIONS	5
13.	WARRANTY INFORMATION	5

Read all instructions carefully prior to use. Observe all warnings and precautions, and keep in mind these instructions that run throughout this Instructions for Use. Failure to observe these warnings or precautions may lead to complications. All the suggestions made in these instructions are designed to serve general guidelines. They are not intended to replace an agreement in system or professional clinical judgment in terms of patient care.

WARNING: U.S. Federal Law restricts this device to sale to only internists or similar physicians.

1. DEVICE DESCRIPTION

QuickCat™ thrombus aspiration catheter is a disposable double-inner-cavity catheter that is designed to remove soft emboli and thrombosis in the vessels of arterial system. The end of the catheter is flexible, conical, and smooth, providing a channel into arterial system to prevent damages. The device is supplied sterile and for a single use.

2. DEVICE COMPONENTS

Catheter: The QuickCat™ thrombus aspiration catheter matches a 6F guiding catheter [I.D. \geq 0.068 (1.73 mm)] and a 0.014" (0.36 mm) wire. The working length of the catheter is 145 cm, and its end is covered with a hydrophilic coating. There is a radiopaque mark at 1 mm of the tip.

Vacuum Assembly: Vacuum assembly comprises a 7.0 (177.8 mm) extension tube, a one-way hemostatic valve, and a 30 ml vacuum syringe with an adjustable locking piston. A 40 μ m filter basket is used to assist in filtering blood and thrombotic substances for a visual or laboratory analysis.

3. INDICATIONS AND INTENDED USE

QuickCat™ thrombus aspiration catheter is used to remove the fresh and soft emboli and thrombosis in the vessels of arterial system.

The product can be used for one time only and should be operated by physicians trained and experienced in diagnostics and interventional procedures.

Standard techniques can be used to arrange vascular access sheathing canals, angiographic catheters and guidewires.

4. CONTRAINDICATIONS

- Use in blood vessel with a diameter < 1.5 mm
- Meridian system
- Removal of fibrous, adherent or calcified materials, such as chronic clots and atherosclerotic plaques

5. WARNINGS

- Do not use it if no guidewire is available, which may lead to a vascular injury.
- Do not try to advance or retract catheter in a resisting reverse direction until the origin of the resistance is determined by fluoroscopy or any other methods. An operation in a resisting reverse direction of the catheter may twist the catheter or cause a vascular injury.
- During the process of handling, if it is too loose or a guidewire loop is observed between guiding catheter and monorail portion of the QuickCat™ thrombus aspiration catheter, the guidewire may twist in the blood vessels when the catheter is advanced or retracted. Before the QuickCat™ catheter is advanced or retracted, remove any loose or looped guidewire to avoid a catheter damage and/or vascular injury.
- In the case of a flow into syringe or restriction, when the catheter is still in patient's body, do not try to flush the aspiration inner cavity of the QuickCat™ thrombus aspiration catheter. By doing so, it may cause a severe damage or death.

- Do not use curved, twisted or damaged catheter, because it may lead to a vascular injury and/or inability to advance or retract the catheter.
- Do not transfer or infuse any diagnostic, inserted, or therapeutic material into blood vessels.

6. PRECAUTIONS

- Prior to use, check if all connections between components are safe, and confirm that the system is completely prepared and whether the vacuum state is appropriate.
- Do not excessively tighten the hemostatic valve on the catheter shaft, or it may cause damage to catheter.
- Caution should be exercised when the QuickCat™ thrombus aspiration catheter is traversed or retracted to pass through the newly deployed drug coating site.
- Do not re-disinfect, reprocess or reuse the device.
- Do not re-disinfect or re-use this device, which may be harmful to the performance of the device, or increase the risk of cross-infection due to improper reprocessing.
Reusing this disposable device may cause a serious patient injury or even death and void the manufacturer's warranty.
- Do not replace any system components with alternative components.

7. POTENTIAL ADVERSE REACTIONS

- Bleeding or hematoma of surgical site
- Damage of anastomosis
- Sudden closure or total occlusion of the treated site or blood vessels
- Pulmonary injury or limb ischemia due to peripheral embolism of fragments
- Local or systemic infection
- Arteriospasm
- Formation of arteriovenous fistula
- Drug reactions, adverse reactions to contrast agents
- Acute myocardial infarction
- Vascular anatomy, perforation, rupture or damage
- Emergency operation
- Death

8. DEVICE PREPARATION

1. Using an aseptic technique, open the bag and transfer the tray to a sterile area.
2. Remove the catheter from the protective ring, and check whether there is any bending or twisting.
3. Fill the 30 ml vacuum syringe with 5 ml – 10 ml of normal saline.
4. Connect the vacuum syringe to piston, connect the piston to extension tube, and connect the extension tube to the holder of the QuickCat™ thrombus aspiration catheter. Make sure all connections are secure.
5. Flush the system with normal saline to ensure complete filling.
6. Rotate the hemostatic valve to the Closed position.
7. Remove the 30 ml vacuum syringe and fill it with saturated saline solution.
8. Fit the 30 ml vacuum syringe onto the system assembly again.

9. PROCEDURE

1. Using a standard technique, insert the guidewire and the guiding catheter (I.D. ≥ 0.068 [1.73 mm]) with a hemostatic valve into the target vessel.
2. Load the QuickCat™ thrombus aspiration catheter on the guidewire.

WARNING: If any resistance is felt during the operation, the cause for the resistance must be determined before the catheter is advanced or retracted. Operating the catheter in a direction opposite to the resistance may lead to a vascular injury. If the catheter is twisted during use, carefully remove it from the patient body, and continue to use a new QuickCat™ thrombus aspiration catheter.

3. During fluoroscopy, push the QuickCat™ thrombus aspiration catheter to the target position.
4. Adequately tighten the hemostatic valve to prevent backflow, but do not tighten it too much to prevent the movement of catheter.
5. The piston is in the Closed position, and pull back the piston of the 30 ml syringe to a desired position for extraction. Turn the piston clockwise and lock the syringe in the desired vacuum position.
6. Determine the correct catheter location through fluoroscopy.
7. Turn the piston to the Open position to start aspiration. Slowly advance the QuickCat™ thrombus aspiration catheter to the target blood vessel. Blood will flow into the syringe until the space is filled up.

- If no blood enters the syringe within 5 seconds, remove the catheter from the patient's body. Outside the patient's body, flush the catheter or replace it with a new catheter.

WARNING: In the case of a flow into the syringe or restriction, when the catheter is still in the patient's body, do not try to flush the aspiration inner cavity of the QuickCat™ thrombus aspiration catheter. By doing so, it may cause a severe damage or death. Remove the catheter from the patient's body and flush the aspiration inner cavity before reuse, or use a new QuickCat™ thrombus aspiration catheter.

- Once the extraction process is completed, turn the piston to the Closed position, and remove the catheter from the patient's body. Use a filter basket with 40 µm tiny holes to filter the extracted blood and blood clots for a visual and/or laboratory analysis.

10. POST-OPERATIVE CARE

After operation, the sheath should be removed by following the care standards of hospital, and provide hemostasis to prevent major hemorrhage of the vascular access position.

11. HOW SUPPLIED

It is sterilized through radiation in the top plate, with a sealed tear-type opened package. It is for a single use; do not re-disinfect, reprocess or reuse it. If the package is unopened or intact, it will be sterile. If the package is suspected of having been opened, do not use the product again. After removing the package, thoroughly check the product to ensure that no damage, bending or twisting occurred during the course of loading.

12. STORAGE CONDITIONS

Store it in a cool and dry place. Avoid excessive exposure to sunlight.

Table 1: Non-standard graphical symbols for medical device labels

Wire Compatibility		Working Length	
Fast Exchange		Vascular Diameter	
Sheath Compatibility		Guiding Catheter Compatibility	
		Non-Pyrogenic	NON-PYROGENIC

13. WARRANTY INFORMATION

Manufacturer limited warranty

If the product is used before the specified "Use before this date" and the package is unopened and intact before use, the manufacturer will authorize the QuickCat™ catheter to be exempted from responsibilities for material defects and processes. Under this warranty, the manufacturer's responsibility is limited to replacing or partially refunding any defective QuickCat™ catheter based on the purchase price. The manufacturer is not responsible for any incidental, special or consequential damages caused by use of QuickCat™ catheter. Any QuickCat™ catheter damage due to misuse, modification, improper storage or operation, or any other failures to comply with these use rules is invalid under this limited warranty. **This limited warranty explicitly supersedes all other warranties, express or implied, including the default warranties of marketability or suitability for a particular purpose.** No individual or entity, including any authorized representative or reseller of the manufacturer, has the right to extend or develop this limited warranty and any purported attempt to do; to the manufacturer, these are non-enforceable.

Attachment to the Instructions for Use

[Product model] 60090-01

[Composition]

		Composition
1	Aspiration catheter	Holder
		Pressure release band
		Proximal tube
		Middle tube
		Distal double-cavity tube
		Hydrophilic coating
		Radioactive label
2	Vacuum tube	
3	Filter mesh	
4	Vacuum syringe	
5	Rotating hemostatic valve	

[Sterilization method] Electron beam sterilization

[Length of validity] 3 years

[Manufacturer] Spectranetics Corporation

[Registered address] 9965 Federal Drive Colorado Springs Colorado 80921 USA

[Production address] 9965 Federal Drive Colorado Springs Colorado 80921 USA

[After-sale service institution] Philips (China) Investment Ltd. Co.

[Address] Building A1, #718 Lingshi Road, JingAn District, Shanghai

[Telephone No.] 8008100038

[Product standard number] 20163772416

[Medical device registration certificate number] 20163772416

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©2019 Spectranetics Corporation
Manufactured by Spectranetics Corporation