AngioSculpt®

PTA Scoring Balloon Catheter

with HydroCross™ Coating

English

INSTRUCTIONS FOR USE

AngioSculpt® Percutaneous Transluminal Angioplasty (PTA) Scoring Balloon Catheter with HydroCross™ Coating

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE: FAILURE TO OBSERVE ALL WARNINGS AND PRECAUTIONS MAY RESULT IN COMPLICATIONS.

STERILE: Sterilized with ethylene oxide gas. Non-pyrogenic. **Do not use if the package is open or damaged.**

CONTENTS: One (1) AngioSculpt PTA Scoring Balloon Catheter with HydroCross™ Coating

STORAGE: Store in a dry, dark, cool place.

I. DEVICE NAME

The device name is AngioSculpt PTA Scoring Balloon Catheter with HydroCross™ Coating; the generic device name is PTA Balloon Dilatation Catheter.

II. DEVICE DESCRIPTION

The AngioSculpt Scoring Balloon Catheter is a standard balloon dilatation catheter with a scoring balloon near the distal tip. One lumen is used for inflation of the balloon with contrast medium; the other lumen permits the use of a guide wire to facilitate advancement of the catheter to and through the stenosis to be dilated. The product is offered on an over-the-wire (OTW) delivery platform.

The distal end of the catheter has a conventional nylon-blend balloon and a nitinol scoring element with spiral struts that wrap around the balloon. There are four struts on balloon diameters 2.0 mm - 3.5 mm and five struts on 4.0, 5.0, and 6.0 mm diameter balloons. The struts create focal concentrations of dilating force, which minimize balloon slippage and assists in the luminal expansion of stenotic arteries.

The balloon has radiopaque markers to aid in positioning the balloon in the stenosis, and is designed to provide an expandable segment of known diameter and length at a specific pressure.

As shown in Figure 1 below, the catheter has a segment which is coated with a hydrophilic coating (HydroCross™ Coating).


Figure 1

Product specification information including balloon diameters, balloon lengths, guidewire compatibility, sheath

compatibility and catheter lengths can be found in Table 1.

Table 1: Product Specifications

Balloon Diameter (mm)	Balloon Length (mm)	Guide Wire Compatibility	Sheath Compatibility (F)	Catheter Length (cm)
2.0	100	0.014"	6	137
2.0	100	0.014"	6	155
2.5	100	0.014"	6	137
2.5	100	0.014"	6	155
3.0	100	0.014"	6	137
3.0	100	0.014"	6	155
3.5	100	0.014"	6	137
3.5	100	0.014"	6	155
4.0	100	0.014"	6	90
4.0	100	0.014"	6	137
4.0	200	0.014"	6	137
5.0	100	0.014"	6	90
5.0	100	0.014"	6	137
5.0	200	0.014"	6	137
6.0	100	0.014"	6	90
6.0	100	0.014"	6	137
6.0	200	0.014"	6	137

III. INDICATIONS

The AngioSculpt PTA Scoring Balloon Catheter is intended for dilatation of lesions in the iliac, femoral, ilio-femoral, popliteal, infra popliteal, and renal arteries, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. Not for use in the coronary neuro-vasculature.

IV. CONTRAINDICATIONS

None known for Percutaneous Transluminal Angioplasty (PTA) procedures.

V. WARNINGS

This device is intended for single (one) patient use only. Do not resterilize and/or reuse, as this can potentially result in compromised device performance and increased risk of inappropriate reesterilization and cross contamination.

The inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis, in order to reduce potential vessel damage.

When the catheter is exposed to the vascular system, it should be manipulated while under high quality fluoroscopic observation.

Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.

Balloon pressure should not exceed the rated burst pressure (RBP). Refer to product label for device specific information. The RBP is based on results of in-vitro testing. At least 99.9% of the balloons (with a 95% confidence level) will not burst at or below their RBP. Use of a pressure monitoring device is recommended to prevent over-pressurization.

Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate the balloon.

Proceed cautiously when using the AngioSculpt catheter in a freshly deployed bare metal or drug eluting stent. The AngioSculpt catheter has not been tested for post-dilation of stents or in lesions distal to freshly deployed stents in clinical studies. Bench testing has shown no additional risk when inserting or withdrawing the AngioSculpt catheter through stents (no interference with stent struts, no retention of or damage to the AngioSculpt catheter).

Use the catheter prior to the "Use Before" (expiration) date specified on the package.

VI. PRECAUTIONS

A thorough understanding of the principles, clinical applications and risks associated with PTA is necessary before using this product.

Any use for procedures other than those indicated in these instructions is not recommended.

The device is not recommended for use in lesions which may require inflation pressures higher than those recommended for this catheter.

Do not use if package is opened or damaged.

Prior to angioplasty, the catheter should be examined to verify functionality, device integrity and to ensure that its size and length are suitable for the specific procedure for which it is to be used.

During and after the procedure, appropriate anti-coagulants, anti-platelet agents and vasodilators should be administered to the patient according to institutional practice for peripheral angioplasty of similar arteries.

Pass the AngioSculpt catheter through the recommended introducer sheath size or minimum size guiding catheter indicated on the product label.

VII. ADVERSE EFFECTS

Possible adverse effects include, but are not limited to, those listed in Table 2 below:

Table 2: Possible Adverse Effects

- Total occlusion of the treated artery
- Arterial dissection or perforation
- Arterial spasm
- Pseudo-aneurysm
- Re-stenosis of the dilated artery
- Embolism
- Thrombus
- Retained device components
- Hemorrhage or hematoma
- Arteriovenous fistula

VIII. MATERIALS REQUIRED FOR USE WITH THE ANGIOSULPT® CATHETER:

WARNING: Use single use items only. Do not resterilize or reuse.

- Femoral introducer sheath and/or guiding catheter (GC):
 - ≥ 6F introducer sheath and/or ≥ 7F GC
- Hemostatic valve

- Radiographic contrast medium diluted ~1:1 with normal saline
- Sterile heparinized normal saline
- 10-cc and 20-cc syringes for flushing and balloon prep
- Inflation device (indeflator)
- Guide wire: 0.014"
- Guide wire introducer
- Guide wire torque device
- Manifold (for pressure monitoring and contrast injection), extension pressure tubing
- Sterile gauze

IX. INSTRUCTIONS FOR USE

Prior to use of the AngioSculpt, examine carefully for damage and device integrity. Do not use if the catheter has bends, kinks, missing components or other damage. Do not use if inner package is open or damaged.

1. Premedicate patients with anti-coagulants, anti-platelet agents and vasodilators according to institutional protocol for PTA procedures.
2. Perform peripheral angiogram in the view best demonstrating the target lesion prior to device deployment.
3. Utilizing standard fluoroscopic technique, position an appropriately sized guide wire of choice beyond the target lesion. Use an exchange length 300 cm guide wire.
4. Using sterile technique, remove an appropriately sized AngioSculpt catheter from the sterile package and place on the sterile field.
 - Size the AngioSculpt a minimum of 0.5 mm LESS than the reference vessel diameter (RVD).
5. Remove the protective tubing from the balloon. Inspect the balloon catheter to ensure that all components are intact.
6. Flush the guide wire lumen by connecting the proximal hub to a 10-cc syringe and injecting heparinized saline into the proximal guide wire lumen until droplets emerge from the distal end.
7. Attach a stopcock to the catheter's balloon inflation port.
8. Attach a 20-cc syringe filled with 2-3 cc of 1:1 mixture of radiographic contrast and normal saline to the stopcock.
9. Open the stopcock to the syringe, aspirate/remove air from the catheter balloon lumen using the 20-cc syringe filled with 2-3 cc of radiographic contrast and leave on vacuum for 30 seconds.
10. Close the stopcock to the catheter balloon inflation port and remove the syringe.
11. Attach inflation device (indeflator), filled with 1:1 mixture of radiographic contrast and normal saline, to the stopcock by creating a meniscus. Avoid introducing air bubbles into the catheter balloon lumen.
12. Open the stopcock to the inflation device and aspirate using the inflation device, locking in vacuum.

NOTE: All air must be removed from the balloon and displaced with contrast medium prior to inserting into the body (repeat steps 8-12, if necessary).
13. To activate the hydrophilic coating, wet the catheter's coated area or wipe a wet sterile gauze over the coated

area immediately prior to insertion into the sheath (see Figure 1).

14. Advance the AngioSculpt catheter over the guide wire (through either a previously placed and appropriately sized hemostatic introducer sheath or guiding catheter) and position at the target lesion utilizing standard fluoroscopic technique. If there is difficulty advancing the catheter through the hemostatic valve of the introducer sheath due to lubricity, gently grasp the catheter with a wet sterile gauze.

NOTE: When backloading the catheter onto the guide wire, the catheter should be supported, ensuring that the guide wire does not come in contact with the balloon. Do not advance or retract the AngioSculpt catheter over the floppy portion of the guide wire. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
15. Inflate the AngioSculpt balloon per the following recommended protocol:
 - increase the inflation pressure by 2 atmospheres every 10-15 seconds until full balloon inflation is achieved
 - do not exceed the rated burst pressure (RBP) printed on the package label
16. Apply negative pressure to the inflation device and confirm that the balloon is fully deflated prior to removing the AngioSculpt catheter. It may be necessary to hold vacuum for up to 1 minute to ensure the balloon is fully deflated prior to retrieval into the sheath.
17. Remove the AngioSculpt catheter. If there is difficulty removing the catheter through the hemostatic valve of the introducer sheath due to lubricity, gently grasp the catheter with a wet sterile gauze.

NOTE: Do not rotate the catheter shaft in excess of 180 degrees when the tip is constrained. Do not rotate the catheter luer hub in excess of five (5) turns during use. Catheter manipulation, including advancement and retraction, should be performed by grasping the shaft.
18. Perform peripheral angiogram (in the same view(s) as step 2) of the target lesion following completion of AngioSculpt treatment.
19. Inspect all components to ensure that the catheter is intact. Follow institutional procedures for disposal of biohazards. If device malfunction occurs or any defects are noted on inspection, flush the guide wire lumen and clean the outer surface of the catheter with saline, store the catheter in a sealed plastic bag, and contact Spectranetics® for further instructions.
20. Complete any additional interventions as clinically indicated (e.g. stent placement).
21. Remove the guide wire and perform peripheral angiography (in the same view(s) as step 2) of the target lesion following completion of all interventions.
22. Remove all catheters and manage the arterial access site according to institutional protocol.
23. Continue treatment with anti-coagulants, anti-platelet agents and vasodilators according to institutional

protocol for PTA procedures.

X. REFERENCES

The physician should consult recent literature on current medical practice regarding balloon dilatation and PTA procedures.

XI. DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

There is no express or implied warranty, including any implied warranty of merchantability or fitness for a particular purpose, on the Spectranetics® product(s) described in this publication. Under no circumstances shall Spectranetics® be liable for any direct, incidental or consequential damages other than as expressly provided by specific law. No person has the authority to bind Spectranetics® to any representation or warranty except as specifically set forth herein.

Descriptions or specifications in Spectranetics® printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Spectranetics® assumes no liability with respect to instruments reused, reprocessed or resterilized.

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Spectranetics Corporation
Spectranetics International B.V.