Quick-Cross
EXTREME
Support Catheter

Instructions for Use

Instructions for Use

English / English

Instructions for Use - Sections by Language

Table of Contents

Pg	Language	Instructions For Use
3	English	Instructions for Use
6	Russian	Инструкции по применению

Instructions for Use

English / English

Table of Contents

1. DESCRIPTION.....	3
2. INDICATIONS FOR USE.....	3
3. CONTRAINDICATIONS.....	3
4. WARNINGS.....	3
5. PRECAUTIONS.....	3
6. ADVERSE EFFECTS.....	3
7. HOW SUPPLIED.....	3
8. COMPATIBILITY.....	4
9. PROCEDURE SET-UP AND DIRECTIONS FOR USE.....	4
10. MANUFACTURER'S LIMITED WARRANTY.....	4
11. NON-STANDARD SYMBOLS.....	4
12. APPLICABLE STANDARDS.....	5
13. DISPOSAL.....	5

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician

1. DESCRIPTION

The Spectranetics Quick-Cross Extreme Support Catheters are intravascular catheters. These catheters are available in a variety of lengths and tip configurations. All models have 3 radiopaque markers spaced equally along the distal shaft to aid in estimating geometry within the vascular system. The distal radiopaque marker is positioned within 3 mm of the distal catheter tip. A standard female luer is placed on the proximal end of each model. The catheter is coated with a lubricious, hydrophilic coating.

2. INDICATIONS FOR USE

Quick-Cross Extreme Support Catheters are intended to guide and support a guidewire during access of the coronary or peripheral vasculature, allow for wire exchanges and provide a conduit for the delivery of saline solutions or diagnostic contrast agents.

3. CONTRAINDICATIONS

None known.

4. WARNINGS

- Maximum Infusion Pressure: **300 psi**, for 0.014" & 0.018" catheters and **500 psi** for 0.035" catheters.
- The catheter is designed and intended for one time use only. **Do not re-sterilize and/or reuse.**
- The catheter should only be used by physicians qualified to perform percutaneous vascular interventions.

5. PRECAUTIONS

- DO NOT resterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing. Reuse of this single use device could lead to serious patient injury or death and voids manufacturer warranties.
- Do not use if device or packaging is damaged.
- Use the catheter prior to the "Use By" date (Expiration Date) specified on the package.

6. ADVERSE EFFECTS

Vascular catheterization and/or vascular intervention may result in complications including but not limited to:

- Vessel dissection, perforation, rupture or total occlusion
- Infection
- Hematoma
- Unstable angina
- Embolism
- Hypo/hypertension
- Acute myocardial infarction
- Arrhythmia, including ventricular fibrillation
- Death

7. HOW SUPPLIED

This device has been sterilized using Ethylene Oxide and is supplied **STERILE**. The devices are designated and designed for **SINGLE USE ONLY** and must not be resterilized and/or reused.

7.1 Sterilization

- The sterility of the product is guaranteed only if the package is unopened and undamaged. Before use, visually inspect the sterile package to ensure that the seals have not been broken. Do not use the catheter if the integrity of the package has been compromised. Do not use catheter if its "Use Before Date," found on package labeling, has been passed.

7.2 Transportation and Storage

Keep dry. Store in a cool dry place. Protect from direct sunlight and high temperature (greater than 55°C or 131°F).

7.3 Inspection Prior to Use

- Before use, examine carefully for defects, all of the equipment to be used. Do not use any equipment if it is damaged or unintentionally opened.

Instructions for Use

English / English

8. COMPATIBILITY

Product Features & Model Numbers Table

Model (Ref)	Working Length	Tip Shape / Type	Hydrophilic Coating Length	GW Compatibility	Minimum Guide I.D.	Minimum Sheath I.D.	Maximum Catheter O.D.	Marker Band Spacing
518-084	135 cm	Straight	100 cm	0.014" 0.36 mm	5F	4F	3.2F	15 mm
518-086	150 cm	Straight					0.042" 1.07 mm	
518-088	90 cm	Straight	60 cm	0.018" 0.46 mm	5F	4F	3.4F	15 mm
518-090	135 cm	Straight	100 cm				0.044" 1.12 mm	
518-092	150 cm	Straight	0.035" 0.89 mm	N/A	5F	4.5F	50 mm	
518-076	65 cm	Straight				45 cm		0.059" 1.50 mm
518-078	90 cm	Straight	70 cm					
518-080	135 cm	Straight	115 cm					
518-082	150 cm	Straight	130 cm					

9. PROCEDURE SET-UP AND DIRECTIONS FOR USE

Note: Follow instructions for use for all equipment to be used with the Quick-Cross Extreme Support Catheter.

1. Preparation: Using sterile technique, open the sterile package. Gently remove the protective hoop with the catheter from the pouch. Fill a sterile standard luer-lock syringe with sterile saline. Before removing the catheter from the hoop, connect the syringe to the catheter proximal luer fitting, flush the catheter and allow the saline to fill the hoop. Set catheter in hoop aside until ready for use.
2. Insertion: Through a previously inserted, appropriately sized guiding catheter or introducer sheath, introduce the catheter over an appropriately sized guidewire (see specifications) using standard technique.
3. Advancement: Use fluoroscopic guidance when advancing the catheter to the desired location within the vasculature.
4. Removal: Gently withdraw the catheter using standard technique, being careful to maintain guidewire position if the guidewire is to remain in place.

Infusion: To perform infusion, withdraw the guidewire and reference the specifications for maximum infusion pressure.

Note: Do not exceed the maximum infusion pressures.

After use, dispose of all equipment in accordance with applicable requirements relating to hospital waste, and potentially bio-hazardous materials.

10. MANUFACTURER'S LIMITED WARRANTY

Manufacturer warrants that the Quick-Cross Extreme is free from defects in material and workmanship when used by the stated "Use By" date and when package is unopened and undamaged immediately before use. Manufacturer's liability under this warranty is limited to replacement or refund of the purchase price of any defective Quick-Cross Extreme. Manufacturer will not be liable for any incidental, special, or consequential damages resulting from use of the Quick-Cross Extreme. Damage to the Quick-Cross Extreme caused by misuse, alteration, improper storage or handling, or any other failure to follow these Instructions for Use will void this limited warranty. **THIS LIMITED WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** No person or entity, including any authorized representative or reseller of Manufacturer, has the authority to extend or expand this limited warranty and any purported attempt to do so will not be enforceable against Manufacturer.

11. NON-STANDARD SYMBOLS

Distal Marker Spacing	
Shape	
GW Compatibility	
Maximum OD	
Maximum PSI/kPa	
Sheath Compatibility	
Working Length	
Quantity	QTY
CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.	Rx ONLY
Importer	

Instructions for Use

English / English

12. APPLICABLE STANDARDS

The standards identified in Table 12.1 were applied to the development of the device.

Table 12.1: Standards applied to the development of the Quick-Cross Extreme.

Standard/Norm and Date	Title
Biocompatibility:	
ISO 10993-1	Biological evaluation of medical devices – Part 1: Evaluation of testing
Clinical Study:	
EN ISO 14155	Clinical Investigation of medical devices for human subjects
Design:	
EN ISO 10555-1	Sterile, single-use intravascular catheters – Part 1: General requirements
US 21 CFR 820	Quality System Regulation
Environmental:	
ISO 14644-1	Cleanrooms and associated controlled environments – Part 1: Classification of air cleanliness.
ISO 14644-2	Cleanrooms and associated controlled environments – Part 2: Specification for testing and monitoring to prove continued compliance with ISO 14644-1.
Labeling:	
EN 556-1	Sterilization of Medical Devices. Requirements for medical devices to be designated "STERILE". Requirements for terminally sterilized medical devices.
EN 980	Graphical symbols for use in labeling of Medical Devices.
EN 1041	Information supplied by the manufacturer with medical devices
ISO 15223	Medical devices Symbols to be used with medical device labels, labeling and information to be supplied
Packaging:	
EN ISO 11607-1	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO 11607-2	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes
ISO 780	Packaging - Pictorial marking for handling of goods
Quality:	
EN ISO 13485	Medical devices — Quality management systems
Risk Management:	
EN ISO 14971	Medical Devices – Application of risk management to medical devices
EN 62366	Medical Devices - Application of usability engineering to medical devices
Sterilization:	
EN ISO 11135	Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices

13. DISPOSAL

Use and disposal must be in accordance with generally accepted medical practice and applicable local, state, and federal laws and regulations.

Products may pose a potential biohazard after use.

If the packaging is damaged and / or the expiration date expires, the product is to be disposed of in accordance with the rules for the treatment of medical waste class A according to SanPiN 2.1.7.2790-10.

The used product is disposed of in accordance with the rules for the management of medical waste class B according to SanPiN 2.1.7.2790-10.

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Spectranetics Corporation
PHILIPS Limited Liability Company