The logo for Thrombocath XT features a thick, black, curved line above the text. The text "Thrombocath" is in a bold, italicized sans-serif font, followed by a registered trademark symbol (®). The "XT" is in a larger, bold, italicized sans-serif font, with a diagonal slash through the "T". Below this, the words "Thrombus Removal System" are written in a smaller, plain sans-serif font.

Thrombocath[®] XT
Thrombus Removal System

INSTRUCTIONS FOR USE

Spectranetics[®]

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to observe warnings and precautions may result in complications. Any recommendations within these instructions are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

Instructions for Use

Table of Contents

Pg	Language		Instructions For Use
3	English	English	Instructions For Use
9	Bulgarian	български език	Указания за употреба
16	Croatian	hrvatski	Upute za uporabu
22	Czech	Česky	Návod k použití
28	Danish	Dansk	Brugsanvisning
34	Dutch	Nederlands	Gebruiksaanwijzing
40	Estonian	Eesti	Kasutusjuhis
46	Finnish	Suomi	Käyttöohjeet
52	French	Française	Mode d'emploi
58	German	Deutsch	Gebrauchsanweisung
64	Greek	Ελληνικά	Οδηγίες χρήσης
70	Hungarian	Magyar	Használati útmutatóaz
76	Icelandic	Íslenska	Leiðbeiningar um notkun
82	Italian	Italiano	Istruzioni per l'uso
88	Latvian	latviešu valoda	Lietošanas pamācība
94	Lithuanian	Lietuvių kalba	Naudojimo instrukcija
100	Norwegian	Norsk bokmål	Bruksanvisning
106	Polish	Polski	Instrukcja obsługi
112	Portuguese	Português	Instruções para Utilização
118	Romanian	Română	Instrucțiuni de utilizare
124	Russian	Русский язык	Инструкции по применению
131	Serbian	српски језик	Uputstvo za upotrebu
137	Slovak	Slovenčina	Návod na použitie
143	Slovenian	Slovenščina	Navodila za uporabo
149	Spanish	Español	Instrucciones de uso
155	Swedish	Svenska	Bruksanvisning
161	Turkish	Türkçe	Kullanım Talimatları

INSTRUCTIONS FOR USE

DEVICE DESCRIPTION	3
INDICATIONS / INTENDED USE	3
CONTRAINDICATIONS	3
WARNINGS	3
PRECAUTIONS	3
POTENTIAL ADVERSE EVENTS	3
VISUAL ALARMS / OPERATING MODES	4
INSTRUCTIONS FOR USE	4
HOW SUPPLIED	5
MANUFACTURER'S LIMITED WARRANTY	7

DEVICE DESCRIPTION

The ThromCat® XT is a single-use, disposable device designed for removing thrombus from vessels of the arterial system. The device consists of a flexible catheter containing a helix designed to track over a 0.014" guidewire. The catheter provides extraction flow at approximately 38 ml/min to remove thrombus. The catheter is advanced and retracted within the vessel to disrupt and remove thrombus. The entire system is supplied sterile and non-pyrogenic. The device is intended for single-use by physicians trained and experienced in diagnostic and interventional techniques.

DEVICE COMPONENTS

The ThromCat XT consists of the following components:

- **Catheter:** Constructed from a nylon blend and stainless steel braid with working length of 150 cm. The catheter operates on a 0.014" guidewire via 6F sheath or 7F guide catheter (≥0.078" minimum ID). The catheter encases a stainless steel helix designed to break up and remove thrombus. The catheter design prevents direct contact between the internal helix and vessel wall. The distal tip of catheter contains an atraumatic non-latex flexible tip.
- **Control Unit:** Contains a motor, which powers the extraction helical pump; an Extraction Line with a 500 ml extraction bag and a power cord that connects the Control Unit to the power supply.
- **Power Supply:** A non-sterile power adaptor that connects the Control Unit power cord to an electrical outlet.

INDICATIONS / INTENDED USE

The ThromCat XT is indicated for percutaneous transluminal thrombus removal in native coronary arteries and infrainguinal arteries.

CONTRAINDICATIONS

- Vessel or lesion that cannot be crossed with a 0.014" diameter guidewire
- Vessel with reference diameter by visual estimation less than (<) 2.5 mm or greater than (>) 7mm
- Severe vessel tortuosity
- Patients with uncontrolled coagulation disorders.
- Infection

WARNINGS

- Do not operate device without a guidewire, as vessel injury may result.
- Do not attempt to advance catheter against resistance until cause of resistance has been determined by fluoroscopy or other means.
- Do not use device if catheter has been kinked, as catheter damage and/or vessel injury may result.
- If excessive slack or a loop in the guidewire is observed between the guide catheter and the monorail segment of the catheter during the procedure, the guidewire may become kinked within the vessel during catheter advancement or retraction. Remove the slack or loop in the guidewire before advancing or retracting the catheter to avoid catheter and/or vessel damage.
- Do not use device in access catheters or vessels with excessive curvatures (<2.5 cm diameter).
- Do not allow catheter to remain in a stationary position while operating, as catheter damage and/or vessel injury may result.
- Confirm access catheter placement to insure firm seating and alignment with target vessel, as catheter damage and/or vessel injury may result.
- Do not use ThromCat XT in vessels with an existing dissection, as vessel damage may result.

PRECAUTIONS

- Manipulation of catheter, while in vivo, should be performed under fluoroscopic guidance.
- Do not operate catheter over guidewire floppy tip as kinking of floppy tip may result.
- Do not force, rapidly advance, or excessively torque catheter as catheter damage may result.
- Do not use the device in calcified lesions as catheter damage may result.
- Do not over tighten Tuohy Borst valve as catheter may become damaged or inoperable.
- Do not use the ThromCat XT catheter in access catheters smaller than 0.078" inner diameter as catheter damage may result.
- Portable and mobile RF communications equipment can affect Medical Electrical Equipment.
- Do not re-sterilize, re-process, or re-use device.
- DO NOT re-sterilize or reuse this device, as these actions can compromise device performance or increase the risk of cross-contamination due to inappropriate reprocessing.
Reuse of this single use device could lead to serious patient injury or death and voids manufacturer warranties.

POTENTIAL ADVERSE EVENTS

Potential adverse events associated with use of the device may include, but are not limited to, the following:

- Access site bleeding/hematoma
- Abrupt closure or total occlusion of treated vessel
- Distal embolization of debris resulting in myocardial infarction or limb ischemia
- Infection

- Vasospasm
- Vessel dissection or perforation
- Adverse reaction to contrast medium
- Emergent surgery
- Death

VISUAL ALARMS / OPERATING MODES

ThromCat XT incorporates controls to maintain safe and effective performance. A green light illuminates when power is supplied to device. A red LED light illuminates when device is inactive. During the device priming procedure, the red LED light will illuminate until the catheter is fully primed. Additionally, during device operation, the red LED light will illuminate if catheter becomes kinked or clogged.

INSTRUCTIONS FOR USE

DEVICE PREPARATION:

1. Using sterile technique, carefully remove the ThromCat XT system from packaging.
CAUTION: Do not use device if catheter has been kinked.
NOTE: A slight curvature in the catheter is normal due to packaging and will not impact device performance or safety.
2. Orient device on procedure table so that catheter is toward patient.
3. Remove and discard paper bands from catheter and extraction bag.
4. Hand off attached extraction bag and power cord to an assistant outside sterile field.
5. Hang extraction bag from IV pole.
6. Connect device power cord to provided power supply.
7. Connect power supply cord to nearby electrical outlet. The green light will illuminate indicating electrical power is being delivered to the device, at the same time the red LED light will illuminate indicating the catheter is inactive.

PRIME PROCEDURE:

1. Submerge distal tip of catheter approximately 5 cm into a container of sterile saline.
NOTE: To successfully prime device, distal tip must be completely submerged during prime sequence.
2. Depress and hold control switch in the ON position, designated by "-", until red LED light turns off (approximately 20 seconds) indicating device is primed.
NOTE: An audible change in speed will occur as the red LED light turns off.
3. Release control switch. The switch will automatically return to its center position and device will continue to run.
4. Depress control switch to OFF position, designated by "O". The red LED light will reilluminate indicating the catheter has reset. **Device is now fully primed and ready for use.**

CATHETER DELIVERY:

1. Place an appropriately sized guide catheter or sheath into vessel using standard percutaneous techniques.
2. Cross the lesion with a 0.014" diameter guidewire.
3. Secure guidewire position; insert guidewire through flexible distal tip and rapid exchange lumen of ThromCat XT catheter (Figure 1).
4. Advance catheter just proximal to thrombus.
5. Standard techniques for placement of vascular access sheaths, angiographic catheters and guidewires may be employed.

Figure 1

WARNING: Incorrect catheter loading may cause catheter damage.

WARNING: Do not attempt to advance catheter against resistance until cause of resistance has been determined by fluoroscopy or other means.

CAUTION: Do not force or excessively torque catheter as this may result in deformation of distal tip, or kinking of catheter.

CAUTION: Do not over tighten Tuohy Borst valve as catheter may become damaged or inoperable.

NOTE: If any problems occur with ThromCat XT, please refer to "TROUBLESHOOTING" section for guidance.

DEVICE OPERATION:

1. With the catheter positioned proximal to lesion site, depress and hold control switch to the ON position, ("."), to activate device.
2. Release control switch when red LED light turns off.
CAUTION: Do not operate device while catheter distal tip is inside guide catheter or sheath as this may result in catheter damage or helix failure.
CAUTION: If red LED light does not turn off, remove catheter from patient before attempting to restart. Attempt to restart device following "PRIME PROCEDURE". If device fails to restart, discard device.
3. Advance catheter approximately 2 mm/sec through lesion beginning proximally and continuing distally. For optimal results, use short back and forth strokes. Repeat as necessary.
CAUTION: Under normal running conditions, the helix will not be visible under fluoroscopy due to its high rotational speed. If helix appears stationary (not rotating) while device is running, turn off device and remove catheter from patient. Discard device.

CAUTION: If device shuts down during use, remove catheter from patient. Attempt to restart device following "PRIME PROCEDURE". If device fails to restart or helix does not rotate, discard device.

- Retract catheter back until just outside of guide catheter or sheath and depress control switch to OFF position ("O").

CAUTION: Do not retract catheter distal tip inside guide catheter or sheath while operating device as this may result in catheter damage or helix failure.

- Repeat steps 1 through 4 of "DEVICE OPERATION" as necessary to achieve desired results.

CAUTION: If necessary to perform multiple procedures, and catheter is removed from patient, catheter should be flushed with saline according to "PRIME PROCEDURE" immediately following removal from patient. Continue to run device during priming sequence until effluent entering extraction line becomes clear.

- When procedure is complete, remove catheter from patient.
- Discard entire device, including power supply and cords, according to hospital procedures.

POST-PROCEDURE CARE

After the procedure, hospital standard of care should be followed for removing sheath and providing hemostasis to prevent bleeding at vascular access site.

HOW SUPPLIED

Supplied sterilized by gamma radiation in peel-open package. Intended for one-time use; do not resterilize, reprocess or reuse. Sterile if package is unopened and undamaged. Do not use the device if there is doubt as to whether product is sterile. Upon removal from package, inspect device to ensure no damage has occurred.

ELECTRICAL STANDARDS COMPLIANCE (EN 60601-1-2:2001, sub-clause 6.8)

- The ThromCat XT device requires installation into service in accordance to the supplied instructions to ensure safety (compliance with electrical standards is provided in Tables 2 through Table 5).
- The following supplied cables were utilized during electrical standard testing and are considered compliant when used with the ThromCat XT device:
 - Power Cable (from device connecting to power supply) with maximum length of 5.5 feet
 - Power Supply with maximum length of 5.0 feet
 - Power Cord with maximum length of 6.5 feet
- Usage of accessories, transducers, and cables other than those provided, may result in increased emissions or decreased immunity of the ThromCat XT device.
- The ThromCat XT device should not be used adjacent to or stacked with other equipment. If adjacent or stacked usage is necessary, the ThromCat XT device should be observed to verify normal operation in the configuration that it will be used.

STORAGE CONDITIONS

Store in a cool, dry location.

TROUBLESHOOTING

Table 1

Problem	Possible Cause(s)	Recommended Action(s)
Device Preparation		
Noticeable damage while still in packaging	Packaging/shipping	DO NOT USE; replace with new device
Device damaged anytime before use	Packaging/shipping	DO NOT USE; replace with new device
	Inappropriate handling	DO NOT USE; replace with new device
Priming Procedure		
Device will not prime; Red LED light will not go out; device will not reach RUN speed	Switch not held long enough	Hold switch until Red LED light turns off or hold for up to 30-40 seconds
	Catheter tip not submerged in saline	Submerge catheter tip in bowl of saline and repeat priming sequence
	Device damaged	DO NOT USE; replace with new device
Catheter Delivery or Removal		
Catheter meets resistance during advancement through guide catheter or sheath	Catheter damaged or kinked	DO NOT USE; replace with new device
	Guide catheter or sheath diameter not compatible with ThromCat XT	Replace guide catheter or sheath with ≥ 0.078 " diameter guide catheter or sheath.
	Guide catheter or sheath damaged or kinked	Replace guide catheter or sheath; repeat PRIME PROCEDURE to verify helix rotation
	Device improperly loaded on guidewire	Verify device is properly loaded on guidewire
Catheter meets resistance during advancement in vessel	Vessel size too small	Verify vessel size is ≥ 2.5 mm
	Device improperly loaded on guidewire	Verify device is properly loaded on guidewire. Remove any slack in guidewire
	Distal tip advanced beyond guidewire	Remove device and guide wire together; Re-deliver guidewire and catheter
Catheter meets resistance during retraction	Device not properly loaded on guidewire	Verify device is properly loaded on guidewire. Remove any slack in guidewire.
	Distal tip advanced beyond guidewire	Remove device and guidewire together
Device Operation		
Device will not run	No power to device	Check power connection; verify LED lights are on
	Switch not fully depressed into ON position	Depress switch fully to ON position
	Switch not held in ON position long enough	Hold switch in ON position until Red LED light turns off
	Attempting to restart device too quickly after turning off	Wait for Red LED light to illuminate before attempting to restart
Catheter helix not rotating when device is operational	Extraction helix damaged	DO NOT USE; replace with new device
Effluent not flowing into collection bag	Extraction pathway clogged	Remove device from patient; repeat PRIME PROCEDURE. If device does not successfully prime, replace with new device
	Extraction helix fractured	DO NOT USE; replace with new device
Device shuts off during use	Extraction pathway clogged	Remove device from patient; repeat PRIME PROCEDURE. If device does not successfully prime, DO NOT USE; replace with new device
	Extraction helix fractured	DO NOT USE; replace with new device
	Device inadvertently deactivated	Verify switch position and operation. If switch does not operate properly, DO NOT USE; replace with new device
Device pulsates or ratchets	The control switch was not kept long enough in the off position	Place the control switch in the off position and wait until the red LED light turns on

Table 2

Guidance and manufacturer's declaration – electromagnetic emissions		
The ThromCat XT is intended for use in the electromagnetic environment specified below. The the user of ThromCat XT should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The ThromCat XT uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause and interference in nearby electronic equipment. The ThromCat XT is suitable for all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage Fluctuations & Flicker Emissions IEC 61000-3-3	Complies	

Table 3

Guidance and manufacturer's declaration – electromagnetic immunity			
The ThromCat XT is intended for use in the electromagnetic environment specified below. The User of the ThromCat XT should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8kV air	+/- 6 kV contact +/- 8kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 1 kV line(s) to line(s) +/- 2 kV line(s) to earth	+/- 1 kV line(s) to line(s) +/- 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95% dip in U_T) for 0,5 cycles 40 % U_T (60% dip in U_T) for 5 cycles 70 % U_T (30% dip in U_T) for 25 cycles <5 % U_T (>95% dip in U_T) for 5 s	<5 % U_T (>95% dip in U_T) for 0,5 cycles 40 % U_T (60% dip in U_T) for 5 cycles 70 % U_T (30% dip in U_T) for 25 cycles <5 % U_T (>95% dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of ThromCat XT requires continued operation during power mains interruptions, it is recommended that the ThromCat XT be powered from an uninterruptible power supply or a battery.
Power frequency (50/60/ Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Table 4

Guidance and manufacturer's declaration – electromagnetic immunity			
The ThromCat XT is intended for use in the electromagnetic environment specified below. The Customer or the user of the ThromCat XT should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the ThromCat XT device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. $d = 1.17 \sqrt{P}$ $d = 1.17 \sqrt{P}$ $d = 2.33 \sqrt{P}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ThromCat XT is used exceeds the applicable RF compliance level above, the ThromCat XT should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the ThromCat XT.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strength should be less than 3 V/m.</p>			

Table 5


Recommended separation distances between portable and mobile RF communications equipment and ThromCat XT			
The ThromCat XT is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of ThromCat XT can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and ThromCat XT as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of the transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = [3.5/\nu_{\text{r}}] \sqrt{P}$	80 MHz to 800 MHz $d = [3.5/E_{\text{r}}] \sqrt{P}$	800 MHz to 2.5 GHz $d = [7/E_{\text{r}}] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.70	11.70	23.33
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			

MANUFACTURER'S LIMITED WARRANTY

Manufacturer warrants that the ThromCat® XT is free from defects in material and workmanship when used by the stated "Use By" date and when package is unopened and undamaged immediately before use. Manufacturer's liability under this warranty is limited to replacement or refund of the purchase price of any defective ThromCat® XT. Manufacturer will not be liable for any incidental, special, or consequential damages resulting from use of the ThromCat® XT. Damage to the ThromCat® XT caused by misuse, alteration, improper storage or handling, or any other failure to follow these Instructions for Use will void this limited warranty. **THIS LIMITED WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** No person or entity, including any authorized representative or reseller of Manufacturer, has the authority to extend or expand this limited warranty and any purported attempt to do so will not be enforceable against Manufacturer.

NON-STANDARD GRAPHICAL SYMBOLS

Table 6

<p style="text-align: center;">-</p> <p style="text-align: center;">ON</p>	<p style="text-align: center;"></p> <p style="text-align: center;">OFF</p>	<p style="text-align: center;"></p> <p style="text-align: center;">Wire Compatibility</p>
<p style="text-align: center;"></p> <p style="text-align: center;">Working Length</p>	<p style="text-align: center;"></p> <p style="text-align: center;">Vessel Diameter</p>	<p style="text-align: center;">Rapid Exchange (RX)</p>
<p style="text-align: center;"></p> <p style="text-align: center;">Sheath Compatibility</p>	<p style="text-align: center;"></p> <p style="text-align: center;">Guide Catheter Compatibility</p>	<p style="text-align: center;">NON-PYROGENIC</p> <p style="text-align: center;">Non-Pyrogenic</p>

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©2011 Spectranetics Corporation
Spectranetics Corporation
Spectranetics International B.V.