



Philips Sonicare AirFloss
Rechargeable AirFloss

AirFloss



HX8111

An easier way to floss

Gently and effectively improves oral health

Fall in love with flossing. Switch to Sonicare.

Fast: cleaning in 60 seconds

- Simple one-button use; as easy as a mouse-click

Simple: point & press

- Microburst technology cleans at the touch of a button
- Charging after two weeks
- Slim, angled nozzle for easy to even the hard-to-reach areas
- Guidance tip for easy placement
- Easy-to-fill reservoir
- Slim, easy to hold handle
- Fits standard Sonicare chargers

Proven improvement of oral health

- Removes up to 99% more plaque between teeth*

PHILIPS
sonicare
the sonic toothbrush

PHILIPS
sense and simplicity


Highlights

Microburst technology




Microburst technology applies a quick burst of pressurized air and micro-water droplets to clean deep between teeth where a toothbrush can't reach

Slim, angled nozzle



Slim, angled nozzle for easy to even the hard-to-reach areas

Slim, easy to hold handle




Easy to manoeuvre around the mouth.

Charging after two weeks



The AirFloss lasts two weeks between charging.

Guidance tip




Simply glide the guidance tip along the gum-line until you feel it settle between the teeth. Provides easy positioning even on those hard-to-reach back teeth.

Fits Sonicare chargers




You only need one charger to charge your Sonicare products.

Simple one-button use



Simple one-button use for easy cleaning in 60 seconds.

Easy-to-fill reservoir



Easy-to-fill reservoir uses only a tiny amount of water for less mess. Fill with mouthwash for enhanced oral health benefits.

Removes up to 99% more plaque

Removes up to 99% more plaque between teeth than brushing with a manual toothbrush alone; Proven to improve oral health effectively.

Specifications

Ease of use

- Handle: Slim ergonomic design for easy handling.
- Battery indicator: Light in power button blinks to indicate when to recharge

Technical specifications

- Battery: Rechargeable
- Charging time: 24 hour(s)
- Operating time (full to empty): 2 weeks or 14 full-

- mouth uses
- Operating time: Replace nozzle every 6 months

Power

- Voltage: Multi-voltage charger

Service

- Warranty: 2-year limited warranty



Issue date 2011-06-30

Version: 9.2.1

© 2011 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com

* Compared to manual toothbrushing alone