



Giving providers what they need

Philips Respironics System One REMstar SE

Designed around your needs

Philips Respironics took input from home care providers and developed a device to meet the demands of an evolving Sleep industry. As an extension of the proven System One family, the REMstar SE rounds out the System One line by providing a robust set of features in one affordable unit. Now, most of your patients can enjoy many of the same comfort benefits, including:

- Flex pressure relief
- Enhanced humidification with rain out protection and Heated Tube option
- Modem and SD card capable

Simplicity at its best

System One REMstar SE offers users simplicity and ease of use through its intuitive, icon-based interface. Philips provides cost savings to home care providers by offering one family of devices with common components such as filters, tubing, humidifiers and more.

PHILIPS

RESPIRONICS

sense and simplicity

Flex and System One Humidity Control for enhanced comfort

Philips Respironics Flex technology has become our clear standard for pressure relief therapy.

A more natural way of breathing

Flex technology helps make PAP therapy more natural by sensing when the patient is exhaling and reducing pressure at those times to make breathing more natural and comfortable.

More comfort with our rainout protection

System One Humidity Control delivers high humidity and ultimate comfort to the patient, while also reducing rainout, even in changing environmental conditions.

A heated hose option is also available to provide more humidity even in cold environments.



Ordering Information	
Devices	Part Number
REMstar SE	IN221S
REMstar SE with Heated Humidifier	IN221HS
REMstar SE with Heated Tube Humidifier	IN221TS
System One Heated Humidifier	IN6H
System One Heated Tube Humidifier	IN6T
Accessories	
SD cards, 10 per pack	1063859
SD card mailers, 10 per pack	1065146

Specifications	
Pressure range	4 to 20 cm H ₂ O
Flex pressure relief	0 to 3
Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	Patient adjustable
Dimensions	7" L x 5.5" W x 4" H 18 cm L x 14 cm W x 10 cm H
Weight	3lbs; 1.36 kg (without humidifier)
Filters	Foam and optional ultra-fine
Device set-up	LED/keypad
Data Storage capacity	Compliance 7-30-day averages, total therapy hours, blower hours SD Card > 1 year On-board > 1 year
Compliance meter	Breathing detection
Altitude compensation	Low, medium or high
Electrical requirements	100-240VAC, 50/60 Hz
Humidification	Integrates with REMstar System One heated humidifier and REMstar System One heated tubing humidifier
Warranty	Two years
Miscellaneous	Fine Pressure Adjustment
Leak compensation	Blower control algorithm

Respironics, REMstar, System One, and Flex are trademarks of Koninklijke Philips Electronics N.V. and its affiliates. All rights reserved.

Please visit www.philips.com/systemone



© 2012 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Broudy NA 7/24/12 MCI 4104925 PN 1101230

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics International Headquarters
+33 1 47 28 30 82
Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9947 0440
Philips Respironics United Kingdom
+44 800 1300 845

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respironics