

PHILIPS
RESPIRONICS

**Premium Alice 6 LDxS
PSG Sleep System**

Upgrade your PSG. Downsize your budget.

Philips Alice 6 upgrade bundle

Upgrade any existing PSG system by purchasing our **premium Alice 6 LDxS PSG Sleep System** with premium accessory package for only **\$16,499!**



Premium Alice 6 LDxS PSG Sleep System includes:

- Base station and headbox (amplifier)
- Sleepware G3 software with Somnolyzer**
- Unlimited core sleepware/software updates
- Unlimited Sleepware/software site licensing
- Installation and on-site customer training
- Four hours of custom report writing per site
- Live 24/7/365 technical/product support
- Loaner units available while under warranty
- Two-year product warranty
- Excellent leasing packages available

\$17,999 – \$1,500* = \$16,499

Alice 6 convenience bundle price

Upgrade discount



Premium accessory package includes:

- Dell OptiPlex 7020 or better desktop computer with three-year pro support warranty
- 24" wide-screen color monitor
- Premium digital pan-tilt zoom video camera package
- Full Pro-Tech sensor kit, including two zRIP Durabelts
- Nonin Xpod oximeter and sensor

For more information about Alice 6 or to receive a demo, please contact your representative or submit your information to sleep.diagnostic@philips.com:

Name:

Telephone:

Email:

To obtain a quote email us at sleep.diagnostic@philips.com

www.philips.com/sleepdx

*Customer must relinquish their existing PSG system to receive the \$1500 upgrade discount.

**Somnolyzer is an optional, per use, auto-scoring feature

Alice 6 system, software and service advantages

Alice 6 LDxS system advantages	Benefits	Sleepware software advantages	Benefits
Fully integrated sleep laboratory system	<ul style="list-style-type: none"> All functions (data, video, audio, & titration control of all RI lab devices) via no more than 2 CAT-5 cables 	Software updates included	<ul style="list-style-type: none"> Substantial value without additional service agreements No yearly add-on charges No hidden costs Eliminate unanticipated fees or licenses
Full international 10-20 EEG inputs (19 + 2 references)	<ul style="list-style-type: none"> Continuous real-time impedance display Standard 1.5 inputs 		
Increased resolution and inputs	<ul style="list-style-type: none"> 2000 Hz sampling rate Up to 500 Hz recording 16-bit sample depth 4 physical ECG leads, including "V" lead EOG independent referencing Pulse Transit Time (PTT) 	Unlimited site licenses included	<ul style="list-style-type: none"> Add sites as your business grows *No additional fees
Real time impedance display	<ul style="list-style-type: none"> Don't lose sight of your patient data with real-time signal quality Excellent training tool for new sleep techs 		
Automatic chin referencing	<ul style="list-style-type: none"> 3 chin EMG inputs EMG technology that continuously monitors all three chin EMG inputs and selects the best pair in real-time 		
Intuitive headbox	<ul style="list-style-type: none"> Built-in amplifier Integrated RIP module Reduces excess wiring and batteries Single plug "key hole" sensor inputs 		
Dedicated differential pressure transducer	<ul style="list-style-type: none"> Powered by headbox (no batteries needed) Snore channel 		
Full Pro-Tech sensor kit, including two zRIP Durabelts	<ul style="list-style-type: none"> Select between adult or pediatric kits An additional \$2,300 value 		
2-year system warranty	<ul style="list-style-type: none"> An additional \$1,500 value over more traditional 1-year warranties Loaner unit while under warranty 		

