

Lumify v1.8 Device Compatibility List



Supported Devices

Device Name	Screen Size	Connectivity Options	Internal Storage	Approximate Continuous Scan Time	USB Type	Supports Reacts
Samsung Galaxy Tab S3 9.7" (Model SM-T820, SM-T825)	9.7 in	Wi-Fi / LTE	32 GB	4 h 30 min	C	Yes
Samsung Galaxy S8 (SM-G950U1)	5.8 in	Wi-Fi / LTE	32 GB	3 h 30 min	C	Yes
Samsung Galaxy S8+ (SM-G955U1)	6.2 in	Wi-Fi / LTE	64 GB	4 h	C	Yes
Panasonic Toughpad FZ-A2	10.1 in	Wi-Fi/LTE	32 GB	3 h 30 min	C	Yes
Panasonic Toughpad FZ-B2 (USB C models only)	7 in	Wi-Fi/LTE	32 GB	3 h 30 min	C	Yes
Samsung Galaxy Note 8 (SM-N950)	6.3 in	Wi-Fi/LTE	64 GB	3 h 30 min	C	Yes
Samsung Galaxy S7 Edge (SM-G935)	5.5 in	Wi-Fi/Cellular	32 GB	2 h 15 min	B	Yes
Samsung Galaxy Tab S2 8" (Only Model SM-T713)	8 in	Wi-Fi	32 GB	2 h 30 min	B	Yes
Samsung Galaxy Tab S2 9.7" (Only Model SM-T813)	9.7 in	Wi-Fi	32 GB	3 h	B	Yes
Samsung J7 (Model SM-J710)	5.5 in	Wi-Fi/ Cellular	16 GB	4 h 15 min	B	Yes
Nokia 6 (Model TA-1003)	5.5 in	Wi-Fi/Cellular	32 GB	2 h 30 min	B	Yes
Huawei Mediapad M3 (Model BTV-W09, BTV-DL09)	8.4 in	Wi-Fi	32 GB	3 h 30 min	B	Yes



Philips Ultrasound, Inc.

22100 Bothell Everett Hwy

Bothell, WA 98021-8431, USA

www.healthcare.philips.com/ultrasound

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.

Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.
239034 Rev O February 2018

Lumify v1.8 Device Compatibility List



No Longer Supported

Device Name	Screen Size	Connectivity Options	Internal Storage	Approximate Continuous Scan Time	USB Type
Samsung Galaxy View	18.4 in	Wi-Fi/LTE	32 GB	5 h	B
Google Nexus 6 (not Nexus 6P)	6 in	Wi-Fi/Cellular	32/64 GB	1 h 45 min	B
NVIDIA SHIELD Tablet (V1)	8 in	Wi-Fi	16 GB	1 h 30 min	B
Samsung Galaxy Tab Active	8 in	Wi-Fi	16 GB	3 h 30 min	B
Samsung Galaxy Tab S	8.5 in	Wi-Fi	16 GB	2 h 15 min	B
Nvidia Shield Tablet K1	8 in	Wi-Fi	16 GB	2 h 30 min	B
Samsung Galaxy Tab S2 8" (Model SM-T710) – Not Compatible with S4-1	8 in	Wi-Fi	32 GB	2 h 30 min	B
Samsung Galaxy Tab S2 9.7" (Model SM-T810) – Not Compatible with S4-1	9.7 in	Wi-Fi	32 GB	3 h	B
Samsung Galaxy Note5	5.7 in	Wi-Fi/Cellular	32/64/128 GB	3 h 15 min	B
Samsung Galaxy Tab S	8.5 in	Wi-Fi	16 GB	2 h 15 min	B
Google Nexus 9	9 in	Wi-Fi/LTE	16/32 GB	2 h	B
Lenovo Yoga Tab 3 (Model YT3-850F, P/N#ZA090094US)	8 in	Wi-Fi	16 GB	4 h	B
Lenovo Yoga Tab 3 (Model YT3-X50F P/N#ZA0H0022US, ZA0H0064US)	10 in	Wi-Fi	16 GB	4 h	B

Tablet Compatibility is subject to change, and Philips reserves the right to make changes without notice. Devices listed under “No Longer Supported” are not tested for use with the current version of Lumify by Philips.



Philips Ultrasound, Inc.
 22100 Bothell Everett Hwy
 Bothell, WA 98021-8431, USA
www.healthcare.philips.com/ultrasound

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.
 Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.
 239034 Rev O February 2018