

# Lumify v1.8 Device Compatibility List



## Supported Devices\*

Device Name	Screen Size	Connectivity Options	Internal Storage	Approximate Continuous Scan Time	USB Type	Supports Reacts
Samsung Galaxy Tab S4 (Model SM-T830)	10.5 in	Wi-Fi	64 GB	5 H 30 min	C	Yes
Samsung Galaxy S9 (SM-G960U1)	5.8 in	Wi-Fi / LTE	32 GB	3 h 45 min	C	Yes
Samsung Galaxy S9+ (SM-G965U1)	6.2 in	Wi-Fi / LTE	64 GB	4 h	C	Yes
Samsung Galaxy Note8 (SM-N950U1)	6.3 in	Wi-Fi / LTE	64 GB	3 h 30 min	C	Yes
Samsung Galaxy S8 (SM-G950U1)	5.8 in	Wi-Fi / LTE	32 GB	3 h 30 min	C	Yes
Samsung Galaxy S8+ (SM-G955U1)	6.2 in	Wi-Fi / LTE	64 GB	4 h	C	Yes
Google Pixel 2 XL	6.0 in	Wi-Fi/LTE	64 GB	3 h	C	Yes
Panasonic Toughpad FZ-A2	10.1 in	Wi-Fi/LTE	32 GB	3 h 30 min	C	Yes
Huawei Mediapad M5 (Model SHT-W09, SHT-AL09)	8.4 in	Wi-Fi or Wi-Fi/LTE	64 GB	3 h 30 min	C	Yes
Panasonic Toughpad FZ-B2 (USB C models only)	7 in	Wi-Fi/LTE	32 GB	3 h 30 min	C	Yes
Samsung Galaxy S7 Edge (SM-G935)	5.5 in	Wi-Fi/LTE	32 GB	2 h 15 min	B	Yes
Samsung Galaxy Tab S2 8" (Only Model SM-T713)	8 in	Wi-Fi	32 GB	2 h 30 min	B	Yes
Samsung Galaxy Tab S2 9.7" (Only Model SM-T813)	9.7 in	Wi-Fi	32 GB	3 h	B	Yes
Samsung J7 (Model SM-J710)	5.5 in	Wi-Fi/LTE	16 GB	4 h 15 min	B	Yes
Nokia 6 (Model TA-1003)	5.5 in	Wi-Fi/LTE	32 GB	2 h 30 min	B	Yes
Huawei Mediapad M3 (Model BTV-W09, BTV-DL09)	8.4 in	Wi-Fi or Wi-Fi/LTE	32 GB	3 h 30 min	B	Yes

\*Devices on the supported list satisfy all of the following criteria: 1. Passed software testing on the latest Lumify software, 2. Passed Philips hardware testing, 3. Is readily available from the device manufacturer. Philips cannot test every device available on the market, so this list represents a range of devices to satisfy a variety of user needs with respect to form factor, price, battery life and connectivity options



**Philips Ultrasound, Inc.**  
 22100 Bothell Everett Hwy  
 Bothell, WA 98021-8431, USA  
[www.healthcare.philips.com/ultrasound](http://www.healthcare.philips.com/ultrasound)

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.  
 Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.  
 239034 Rev R September 2018

# Lumify v1.8 Device Compatibility List



## No Longer Supported\*\*

Device Name	Screen Size	Connectivity Options	Internal Storage	Approximate Continuous Scan Time	USB Type
Samsung Galaxy Tab S3 (Model SM-T820,SM-T825)	9.7 in	Wi-Fi or Wi-Fi/LTE	32 GB	4 h 30 min	C
Samsung Galaxy View	18.4 in	Wi-Fi/LTE	32 GB	5 h	B
Google Nexus 6 (not Nexus 6P)	6 in	Wi-Fi/Cellular	32/64 GB	1 h 45 min	B
NVIDIA SHIELD Tablet (V1)	8 in	Wi-Fi	16 GB	1 h 30 min	B
Samsung Galaxy Tab Active	8 in	Wi-Fi	16 GB	3 h 30 min	B
Samsung Galaxy Tab S	8.5 in	Wi-Fi	16 GB	2 h 15 min	B
Nvidia Shield Tablet K1	8 in	Wi-Fi	16 GB	2 h 30 min	B
Samsung Galaxy Tab S2 8" (Model SM-T710) – Not Compatible with S4-1	8 in	Wi-Fi	32 GB	2 h 30 min	B
Samsung Galaxy Tab S2 9.7" (Model SM-T810) – Not Compatible with S4-1	9.7 in	Wi-Fi	32 GB	3 h	B
Samsung Galaxy Note5	5.7 in	Wi-Fi/Cellular	32/64/128 GB	3 h 15 min	B
Samsung Galaxy Tab S	8.5 in	Wi-Fi	16 GB	2 h 15 min	B
Google Nexus 9	9 in	Wi-Fi/LTE	16/32 GB	2 h	B
Lenovo Yoga Tab 3 (Model YT3-850F, P/N#ZA090094US)	8 in	Wi-Fi	16 GB	4 h	B
Lenovo Yoga Tab 3 ( Model YT3-X50F P/N#ZA0H0022US, ZA0H0064US)	10 in	Wi-Fi	16 GB	4 h	B

\*\*The devices on this list represent devices that at one point supported the Lumify system but no longer do for one or more of the following reasons: 1. Did not pass software testing on the latest Lumify software, 2. Did not pass Philips hardware testing, or 3. Is no longer available from the manufacturer. Lumify may still run on these devices, but due to the above criteria, Philips can no longer verify Lumify compatibility



# Lumify v1.8 Device Compatibility List



Tablet Compatibility is subject to change, and Philips reserves the right to make changes without notice. Devices listed under “No Longer Supported” are not tested for use with the current version of Lumify by Philips.



**Philips Ultrasound, Inc.**  
22100 Bothell Everett Hwy  
Bothell, WA 98021-8431, USA  
[www.healthcare.philips.com/ultrasound](http://www.healthcare.philips.com/ultrasound)

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.  
Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.  
239034 Rev R September 2018