

*“If there is a problem, it can be addressed sooner, rather than waiting a few weeks for an appointment”*  
David C., patient on EncoreAnywhere



*“Remotely accessing data via modems enables us to be more efficient, getting patients onto effective therapy more quickly”*  
Sara Parsons, Respiratory Physiologist  
St. George's Hospital, United Kingdom



Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

© 2015 Koninklijke Philips N.V. All rights reserved.  
Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Broudy CDL 1/19/15 MCI 4106435 PN 1119728

How to reach us:  
[www.healthcare@philips.com](http://www.healthcare@philips.com)

Product information:  
[www.philips.com/respironics](http://www.philips.com/respironics)



Connecting you  
**with your patients**

A web-based patient management system that supports the monitoring of home ventilation and sleep therapy patients...

# Introducing EncoreAnywhere

In one centralized and comprehensive database, easily manage your patients' records, upload and access patient therapy data, and generate clinical and business reports.

Wireless modems collect therapy data direct from patients' homes each day, giving you access to the most up-to-date therapy information, which aids proactive and confident clinical decision making.

EncoreAnywhere is accessible anytime, from any internet-enabled location, securely connecting a care team with its patients and promoting collaborative care.

EncoreAnywhere falls under the medical device directive.

## Manage your day with 'My Day'

The 'My Day' page automatically prioritizes the patients that require your clinical attention.

By defining the therapy and time parameters on which patients are closely monitored, only those not meeting these parameters are shown on your 'My Day' page – allowing you to review and manage any interventions from one centralized page.

And with wireless modem technology, up-to-date therapy data can be monitored quickly and remotely, allowing for more timely interventions.



Further enhancing our sleep and home ventilation offering, EncoreAnywhere is compatible with:

### SLEEP

System One REMstar



### SLEEP

System One BiPAP autoSV Advanced



### NIV

System One BiPAP S/T



### NIV

BiPAP A30



### NIV

BiPAP A40



## A low maintenance database complemented with high security

### Low maintenance

EncoreAnywhere's web-based technology removes the need for hospital server space and local installations.

Updates are automatic and discreet, removing any disruption to clinics and clinician access.

### High security

The EncoreAnywhere database is hosted by a secure data center owned by an independent European service provider. The data center is fully compliant with EU data protection laws with strict privacy and security measures. It maintains ISO 27001 certification and is subject to audit.

The infrastructure is operated and maintained in accordance with an Information Security Management System (ISMS).

## EncoreAnywhere packages

### Essentials

- Manual therapy data upload
- Unlimited access to the EncoreAnywhere database
- Unlimited number of clinical reports
- Unlimited number of business reports
- Unlimited server space
- Automatic software updates

### Enhanced

- All the benefits of Essentials but with wireless modem technology. Additional benefits include:
  - Automatic daily therapy data upload
  - Remote prescription changes
  - Instant messaging to devices\*

\*System One REMstar 60 Series

For more information, please speak to your Philips Respironics representative.

### Make confident clinical decisions

EncoreAnywhere is simple to use yet comprehensively covers your clinical requirements.

Therapy data can only be accessed by authorized members of the patient's care team, from any internet-enabled location.

## Sleep

- Cheyne-Stokes respiration
- Clear airway apnea
- Obstructive airway apnea
- Hypopnea

## NIV

- Tidal volume
- Patient-triggered breaths

### For Sleep and NIV

- Daily details
- Patterns of use
- Short-term trends
- Long-term trends
- Compliance and adherence
- Therapy statistics
- Pressure support
- Hours of usage
- Pressure
- AHI
- Breaths per minute
- Leak

### Generate clinical and business reports with ease

EncoreAnywhere offers customizable reports for cross-functional team use and auditing:

#### Clinical

- Download device usage and therapy compliance statistics
- View patient therapy set-up and device information

#### Auditing

- Review sleep/NIV service growth over time
- Monitor therapy equipment allocation and modem activity
- Download patient demographic and therapy status details


Additionally, data can be easily exported to Microsoft Excel and adapted for more specific auditing.



### How therapy data is uploaded to EncoreAnywhere


#### Automatic therapy data upload

The device automatically uploads therapy data each day. Prescription changes and instant messages\* can be sent remotely to the device.



#### Manual therapy data upload

Therapy data stored on the device's data card is manually uploaded via a computer.



\*System One REMstar 60 Series