

Lumify Device Compatibility



Device Name	Screen Size	Type	Approximate Continuous Scan Time	Connectivity Options	Internal Storage	USB Type
Samsung Galaxy Tab S3 9.7" (Model SM-T820)	9.7 in	Tablet	4 h 30 min	Wi-Fi	32 GB	C
Samsung Galaxy S8 (SM-G950U1)	5.8 in	Phone	3 h 30 min	Wi-Fi / LTE	32 GB	C
Samsung Galaxy S8+ (SM-G955U1)	6.2 in	Phone	4 h	Wi-Fi / LTE	64 GB	C
Panasonic Toughpad FZ-A2	10.1 in	Tablet	3 h 30 min	Wi-Fi/LTE	32 GB	C
Panasonic Toughpad FZ-B2 (USB C models only)	7 in	Tablet	3 h 30 min	Wi-Fi/LTE	32 GB	C
Samsung Galaxy S7 Edge	5.5 in	Phone	2 h 15 min	Wi-Fi/Cellular	32 GB	B
Samsung Galaxy Tab S2 8" (Only Model SM-T713)	8 in	Tablet	2 h 30 min	Wi-Fi	32 GB	B
Samsung Galaxy Tab S2 9.7" (Only Model SM-T813)	9.7 in	Tablet	3 h	Wi-Fi	32 GB	B
Samsung Galaxy View	18.4 in	Tablet	5 h	Wi-Fi/LTE	32 GB	B
Samsung J7 (Model SM-J710)	5.5 in	Phone	4 h 15 min	Wi-Fi/ Cellular	16 GB	B
Nokia 6 (Model TA-1003)	5.5 in	Phone	2 h 30 min	Wi-Fi/Cellular	32 GB	B
Huawei Mediapad M3 (Model BTV-W09, BTV-DL09)	8.4 in	Tablet	3 h 30 min	Wi-Fi	32 GB	B



Philips Ultrasound, Inc.
 22100 Bothell Everett Hwy
 Bothell, WA 98021-8431, USA
www.healthcare.philips.com/ultrasound

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.
 Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.
 239034 Rev N October 2017

No Longer Supported

Device Name	Screen Size	Type	Approximate Continuous Scan Time	Connectivity Options	Internal Storage	USB Type
Google Nexus 6 (not Nexus 6P)	6 in	Phone	1 h 45 min	Wi-Fi/Cellular	32/64 GB	B
NVIDIA SHIELD Tablet (V1)	8 in	Tablet	1 h 30 min	Wi-Fi	16 GB	B
Acer Iconia Tab 10 A3-A30	10 in	Tablet	2 h 30 min	Wi-Fi	16 GB	B
Samsung Galaxy Tab Active	8 in	Tablet	3 h 30 min	Wi-Fi	16 GB	B
Samsung Galaxy Tab S	8.5 in	Tablet	2 h 15 min	Wi-Fi	16 GB	B
Google Nexus 9	9 in	Tablet	2 h	Wi-Fi/LTE	16/32 GB	B
Nvidia Shield Tablet K1	8 in	Tablet	2 h 30 min	Wi-Fi	16 GB	B
Samsung Galaxy Note5	5.7 in	Phone	3 h 15 min	Wi-Fi/Cellular	32/64/128 GB	B
Lenovo Yoga Tab 3 (Model YT3-850F, P/N#ZA090094US)	8 in	Tablet	4 h	Wi-Fi	16 GB	B
Lenovo Yoga Tab 3 (Model YT3-X50F P/N#ZA0H0022US, ZA0H0064US)	10 in	Tablet	4 h	Wi-Fi	16 GB	B
Samsung Galaxy Tab S2 8" (Model SM-T710)	8 in	Tablet	2 h 30 min	Wi-Fi	32 GB	B
Samsung Galaxy Tab S2 9.7" (Model SM-T810)	9.7 in	Tablet	3 h	Wi-Fi	32 GB	B

Tablet compatibility is subject to change, and Philips reserves the right to make changes without notice.

