

PHILIPS

Respiratory drug delivery

Choosing an **aerosol delivery system** for your drug development program

For more information please call us at **+44(0) 800-130-0857** or for the **US dial +44 800 130 0857***

**This number is not toll free if dialed from the US. Regular charges will apply.*

Does your drug have any issues with solubility, osmolality, or biological activity?

NO

Consider a **jet nebulizer**

YES

Consider a **mesh device**

Do you need high efficiency, portability or device differentiation for your drug?

YES

Does your drug have challenges with high cost, narrow therapeutic window or patient adherence?

YES

Do you need features such as adherence monitoring, compliance monitoring, feedback or training ?

NO

YES



SideStream Plus jet nebulizer

Breath-enhanced, continuous and designed to improve drug delivery and reduce waste.



InnoSpire Go portable mesh nebulizer

Highly efficient, continuous delivery of medication with virtually silent operation.



Micro breath-activated mesh nebulizer**

For delivery of small volumes of medication over a few breaths



I-neb breath-activated mesh nebulizer

Precision dosing with patient compliance tracking and adherence monitoring.

**Non-CE marked/Investigational Device. Limited to investigational use.