

GD4220



Article Codes: **GD4220** For 2 Monitors up to 19" for Portegra2 Suspension Arms
GD4220-L For 2 Monitors up to 21/24" for Portegra2 Suspension Arms
GD4220-XL For 2 Monitors up to 24/26" for Portegra2 Suspension Arms


CE
 Made in Germany

Portegra2 Dual Monitor Holder – GD4220 Models

Monitor Support for Two Flat Screen Displays

The GD4220 can easily hold up to two 26" monitors. The large curved handlebar protects the monitor and provides a means for quick & easy positioning. Combined with the mobility and safety of the Portegra2 suspension arm, the GD4220 offers the user a worry-free and flexible application.

Type:	Monitor holder with handle and VESA 100 adapter for two flat panels, external cabling	
Weight Load:	GD4220	15.0 kg / 33.0 lbs
(Max.)	GD4220-L	14.5 kg / 32.0 lbs
	GD4220-XL	14.2 kg / 31.0 lbs
Weight:	GD4220	3.0 kg / 6.6 lbs
	GD4220-L	3.5 kg / 7.7 lbs
	GD4220-XL	3.8 kg / 8.4 lbs
Colour:	RAL 9002	
Test mark:	CE, UL	



* (Dimensions for L and XL versions)

MAVIG GmbH
 X-Ray Protection
 & System Solutions

Stahlgruberring 5
 81829 Munich
 Germany

PO Box 82 03 62
 81803 Munich
 Germany

Phone: +49 (0) 89 / 420 96 - 0
 Fax: +49 (0) 89 / 420 96 - 200
 e-Mail: info@mavig.com

www.mavig.com