

# X-Ray Protective Shield OT50001 & OT90001



Accessories, also available:



**STEA-OT5** Disposable, gamma sterile, complete shield cover.



**OTSTE1** Disposable, gamma sterile, partial cover for lower third.

**OT50001**  
**OT90001**  
**E-OT25B05**  
**STEA-OT5** or **OTSTE1**

Shield complete with Portegra2 Suspension Arm (Extension arm length 75 cm)  
Shield complete with Portegra2 Suspension Arm (Extension arm length 95 cm)  
Replacement Shield and Connection Bar  
Complete (STEA-OT5) or partial shield (OTSTE1) disposable, sterile cover


**CE**, IEC 60601-1  
Made in Germany

## MAVIG X-Ray Protective Shield – OT50001 & OT90001

MAVIG's products are designed to offer solid protection and optimized functionality. The high quality acrylic shield allows for an excellent transmittance and unimpeded view. All edges and corners have been rounded and shatter-proof materials are used to minimize risk of injury.

MAVIG's patented Portegra2 suspension arm provides its industry known safety features and movement flexibility. The centrally guided shield with a special patient cut-out can be easily positioned over the patient and quickly removed in case of an emergency.

- **Lead Acrylic Shield:** 60 x 76 cm (23.6 x 29.9 in)
- **Lead Equivalent:** Pb 0.50 mm
- **OT50001:** Extension / Spring Arm  
75 cm / 91 cm (29.5 in / 35.8 in)  
Colour RAL 9002
- **OT90001:** Extension / Spring Arm  
95 cm / 91 cm (37.4 in / 35.8 in)  
Colour RAL 9002



**MAVIG GmbH**  
X-Ray Protection  
& System Solutions

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone: +49 (0) 89 420 96 - 0  
Fax: +49 (0) 89 420 96 - 200  
e-Mail: info@mavig.com

[www.mavig.com](http://www.mavig.com)