


RA614 Thyroid Gland and Sternum Protection



Hygienic Cover in Colour Ocean

- RA614CC*** Thyroid/Sternum Shield RA614 with NovaLite X-Ray Protective Material
RA614CC*AK Thyroid/Sternum Shield Similar to RA614 with an Exchangeable Velcro Band
RA614F Washable hygienic cover made of ComforTex®-MF material in colour Ocean

RA614 – Thyroid Gland and Sternum Protection

Today we know that the thyroid gland is particularly sensitive to radiation, and requires special protection. It has therefore become common place to wear appropriate additional protection for high dosage X-ray examinations.

Our thyroid gland and sternum protection shows itself through outstanding features such as first-class design and ideal fit, while keeping in mind all aspects of radiation protection. The well conceived idea behind ensures that one size fits all neck sizes from 32 to 48 cm. An uncomplicated velcro fastener, additionally available as exchangeable, gives plenty of range for individual adjustments.

In order to keep up with the hygienic standards, if the shield is used by several persons, it is recommended to order the collar with a special washable hygienic cover made of ComforTex®-MF material in colour Ocean.

Thyroid gland and sternum protection with snap buttons are for the apron model series RA631 and RA632.

- Lead Equivalent Pb 0.50 mm
- Edging All edges covered with a seam band
- Outer Cover ComforTex® HPMF

(*CC - ColourCode: 1 = Curacao, 2 = Orchidee, 3 = Regatta, 4 = Indian Summer, 5 = Limette, 7 = Ocean)

MAVIG GmbH
X-Ray Protection
& System Solutions

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone: +49 (0) 89 420 96 – 0
Fax: +49 (0) 89 420 96 – 200
e-Mail: info@mavig.com

www.mavig.com