

<i>Differences between IPL devices</i>	<b>Lumea Comfort</b>	<b>Lumea Essential</b>			<b>Lumea Precision</b>			
	SC1981/ SC1982	SC1991	SC1992	SC1996	SC2003	SC2004	SC2006	SC2008
Number of flashes when fully charged	corded, not applicable	corded, not applicable	corded, not applicable	corded, not applicable	370 flashes	370 flashes	370 flashes	370 flashes
Total amount of flashes	100.000 flashes	100.000 flashes	100.000 flashes	100.000 flashes	100.000 flashes	100.000 flashes	140.000 flashes	140.000 flashes
Lifetime lamp (based on total amount of flashes)	4 years	5 years	5 years	5 years	5 years	5 years	5 years	5 years
Speed (based on treating 2 lower legs)	22 minutes	15 minutes	15 minutes	15 minutes	15 minutes	15 minutes	15 minutes	15 minutes
Pouch included	no	yes	yes	yes	yes	no	yes	yes
Window size of attachment	2.7 cm <sup>2</sup>	4 cm <sup>2</sup>	4 cm <sup>2</sup>	4 cm <sup>2</sup>	4 cm <sup>2</sup>	4 cm <sup>2</sup>	4 cm <sup>2</sup>	4 cm <sup>2</sup>
body	n.a.	n.a.	n.a.	2 cm <sup>2</sup>	2 cm <sup>2</sup>	n.a.	2 cm <sup>2</sup>	2 cm <sup>2</sup>
face	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	2 cm <sup>2</sup>
bikini								
Slide&Flash	Yes	No	Yes	Yes	Yes	Yes	Yes	Yes
Corded/Cordless	Corded	Corded	Corded	Corded	Cordless	Cordless	Cordless	Cordless
Skin Tone Sensor	Yes	No	No	Yes	No	No	No	No
Skin types	I - IV	I - IV	I - IV	I - IV	I - V	I - V	I - V	I - V