


support

for you and your child



The Wisp pediatric mask has a fun, child-friendly giraffe print and soft, fabric materials. To help ease children into therapy and feel

more comfortable with the mask, we've created a fantasy world featuring our cartoon giraffe character, Jacky, and friends Sami the Seal and Tucker the Turtle.

We have also created a storybook to help the youngest children with breathing problems understand the benefit of wearing a mask.



Jacky's adventures and instructional videos to help you and your child get comfortable with therapy can be found on Philips Respironics YouTube channel: philips.to/wisppediatric.



PHILIPS
RESPIRONICS

Wisp

Pediatric nasal mask



Rediscover dreams

© 2016 Koninklijke Philips N.V. All rights reserved.
Specifications are subject to change without notice.

www.philips.co.uk/sleep

Clay Comm JW 6/9/16 MCI 4107612

Minimal contact, **maximum comfort**

Wisp pediatric nasal mask— designed specifically for children

Wisp pediatric applies a child-friendly fabric pattern, a modified cushion curvature, and family support tools to provide a positive experience for your child.

- The cushion shape is designed to fit the curvature of an infant's face and nose.
- An open line of sight allows children to comfortably wear the mask during their normal daily activities.
- A Leak Correction Dial lets you adjust the mask without removing it or waking the child.
- Its three cushion sizes are designed specifically for young children.



Supporting your family with our family of products

The Wisp family of masks

Serves infants to adults, providing continuity and consistency in their therapy experience.



Wisp pediatric



Wisp youth



Wisp adult

Trilogy series of ventilators

Sensitive delivery of non-invasive ventilation therapy in the home or wherever it's needed.



For more information about Wisp masks or Trilogy ventilators, visit www.philips.co.uk/respironics.

Avent products

Makes feeding, soothing and sleeping more natural and more comfortable for parent and child.

- Baby bottles
- Pacifiers and Soothies
- Breastfeeding supplies
- Bottle warmers and sterilisers
- Toddler drinking cups
- Baby monitors



For more information about Avent products, visit www.philips.co.uk/avent.

Sonicare for Kids

A rechargeable sonic toothbrush with dynamic action helps kids brush better, longer, and more regularly.



For more information about Sonicare for Kids, visit www.philips.co.uk/sonicare.