

Pro-Tech ezRIP PSG compatibility chart

PSG w/1.5mm connectors	Pro-Tech ezRIP adult kit, Alice 5/LE, 1.5mm, 84"	Pro-Tech ezRIP adult kit, 1.5mm, 84"	Pro-Tech ezRIP pediatric kit, 1.5mm, 84"	SLP/ SleepSense RIP inductive interface	AMBU, SleepMate RIPmate
Alice 5/LE	✓			✓	✓
Cadwell		✓	✓	✓	✓
Compumedics E & S		✓	✓	✓	✓
Compumedics Graef		✓	✓	✓	✓
Embla Artisan Amp		✓	✓	✓	✓
Embla Lamont Amp		✓	✓	✓	✓
Embla Monet Amp		✓	✓	✓	✓
Grass Technologies		✓	✓	✓	✓
Nicolet		✓	✓	✓	✓
Nihon Kohden		✓	✓	✓	✓
Sandman SD20		✓	✓	✓	✓
Sandman SD32		✓	✓	✓	✓
Sensormedics Cephalo/Omega		✓	✓	✓	✓
SomnoStar z4		✓	✓	✓	✓
Stellate		✓	✓	✓	✓
Xltek		✓	✓	✓	✓

PSG w/249 connectors	Pro-Tech ezRIP adult kit, 249, 84"	Pro-Tech ezRIP pediatric kit, 249, 84"	Pro-Tech ezRIP adult kit, 249, abd 16"/ thx 20"	Pro-Tech ezRIP pediatric kit, 249, abd 16"/ thx 20"	SLP/ SleepSense RIP inductive interface cable	AMBU, SleepMate RIPmate
Nihon Kohden/Nomad			✓	✓		
Nihon Kohden/Sphinx			✓	✓		
Embletta			✓	✓	✓	
Embla S4000/S4500	✓	✓	✓	✓	✓	
Embla N7000	✓	✓	✓	✓	✓	

HST w/harwin connector	Pro-Tech ezRIP Stardust II kit, harwin
Stardust II	✓

Respironics and ezRIP are trademarks of Koninklijke Philips Electronics N.V. and its affiliates. All rights reserved. All other trademarks are the property of their respective owners.

©2012 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

SB 03/21/12 MCI 4104765

PHILIPS

RESPIRONICS
sense and simplicity