



PHILIPS

Trumpf Medical

Making the difference **with Live Image Guidance**

Seamlessly perform open and minimally invasive procedures in a single room

As procedures evolve, surgeons are performing an increasingly broad range of procedures. A Hybrid operating room (OR) is designed to provide the optimal environment for open surgery, hybrid, as well as endovascular procedures. To create this exceptional care environment, a well-designed Hybrid OR enhances the patient and staff experience. Integration is crucial and Philips understands that interoperability with your partner of choice is the key to creating an exceptional care environment.

A truly multifunctional room

In Philips Hybrid OR, the combination of the Philips Allura with the Trumpf TruSystem 7500 Hybrid plus seamlessly interfaces leading edge interventional X-ray with a leading edge OR table. A reliable investment that will provide you with a truly multifunctional room where you can perform interventional, hybrid, and traditional surgical procedures. It supports high utilization which supports a high return on your investment.



Fast and confident positioning

The Allura and Trumpf OR table are interfaced to work smoothly and flexibly with each other. The interface between the Allura stand and the Trumpf OR table facilitates the stop of gantry and table movements after a Bodyguard signal so you can move the Allura at high speeds with confidence. C-arc functions, such as automatic position control and store recall, can be used with the Trumpf OR table.

Smooth and easy control

Intuitively pan and free-float the Trumpf tabletop with the Trumpf SensorControl handle. The Allura stand and Trumpf OR table are controlled via their individual control modules.



Transfers made simple

Use the motorized shuttle to easily exchange tabletops and transfer patients.

Excellent flexibility and efficiency

When a leading edge interventional X-ray system works side-by-side with a leading edge OR table, your medical teams benefit from excellent flexibility and efficient workflow. The Trumpf OR table provides extreme adjustment angles and positions to support the full range of surgical procedures.

The Trumpf TruSystem 7500 Hybrid plus OR table system has a fixed column in the floor and can be outfitted with different tabletops:

- The Trumpf TruSystem 7500 Carbon FloatLine Imaging Tabletop. The special Carbon Fiber tabletop supports a wide range of procedures in endovascular and cardiovascular surgery, as well as in interventional cardiology, radiology, orthopedics, and traumatology.
- The Trumpf TruSystem 7500 SQ14-XTRA Imaging Tabletop supports neuro and cranial procedures. Carbon segments are permanently mounted at the head and foot end of the tabletop.
- The Trumpf TruSystem 7500 Universal Tabletop for open surgical procedures can be positioned in many different ways, following the requirements of the specific procedure. In these cases, the Allura C-arc should be parked away.



*The product is not available in all countries. Consult your Philips representative for more information.

© 2015 Koninklijke Philips N.V. All rights reserved.
Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
(Royal Philips) or their respective owners.



How to reach us
Please visit www.philips.com/healthcare
healthcare@philips.com