

# Lumify v1.8 Device Compatibility List



## Supported Devices\*

Device Name	Screen Size	Connectivity Options	Internal Storage	Approximate Continuous Scan Time	USB Type	Supports Reacts
Samsung Galaxy Tab S6 (Model SM-T860)	10.5 in	Wi-Fi	128 GB	TBD	C	Yes
Samsung Galaxy Tab S4 (Model SM-T830, SM-T835, SM-T837)	10.5 in	Wi-Fi or Wi-Fi/LTE	64 GB	5 h 30 min	C	Yes
<b>Samsung Galaxy Tab A w/ S-Pen (SM-P200, SM-P205)</b>	<b>8.0 in</b>	<b>Wi-Fi / LTE</b>	<b>32 GB</b>	<b>3 h 16 min</b>	<b>C</b>	<b>Yes</b>
Samsung Galaxy Tab A- 2018 (Model SM-T387)	8.0 in	Wi-Fi / LTE	32 GB	3 h 29 min	B	Yes
<b>Lenovo Tab4 8</b>	<b>8.0 in</b>	<b>Wi-Fi / LTE</b>	<b>16 GB</b>	<b>4 h 26 min</b>	<b>B</b>	<b>Yes</b>
Samsung Galaxy S10e (SM-G970U)	5.8 in	Wi-Fi / LTE	128 GB	4 h 7 min	C	Yes
Samsung Galaxy S10 (SM-G973U)	6.1 in	Wi-Fi / LTE	128 GB	4 h 3 min	C	Yes
Samsung Galaxy S10+ (SM-G975U)	6.4 in	Wi-Fi / LTE	128 GB	3 h 34 min	C	Yes
Samsung Galaxy S9 (SM-G960U)	5.8 in	Wi-Fi / LTE	32 GB	3 h 45 min	C	Yes
Samsung Galaxy S9+ (SM-G965U)	6.2 in	Wi-Fi / LTE	64 GB	4 h	C	Yes
Samsung Galaxy Note10+ (SM-N975U, SM-N976V)	6.8 in	Wi-Fi / LTE / 5G	256 GB	4 h 52 min	C	Yes
Samsung Galaxy Note10 (SM-970U)	6.3 in	Wi-Fi / LTE	256 GB	TBD	C	Yes
Samsung Galaxy Note9 (SM-N960U)	6.25 in	Wi-Fi / LTE	128 GB	3 h 26 min	C	Yes
Samsung Galaxy Note8 (SM-N950U)	6.3 in	Wi-Fi / LTE	64 GB	3 h 30 min	C	Yes
Samsung Galaxy S8 (SM-G950U)	5.8 in	Wi-Fi / LTE	32 GB	3 h 30 min	C	Yes
Samsung Galaxy S8+ (SM-G955U)	6.2 in	Wi-Fi / LTE	64 GB	4 h	C	Yes
<b>Google Pixel 4</b>	<b>5.7 in</b>	<b>Wi-Fi / LTE</b>	<b>64 GB</b>	<b>3 h 5 min</b>	<b>C</b>	<b>Yes</b>
Google Pixel 3a / Pixel 3a XL	5.6 in/ 6.0 in	Wi-Fi / LTE	64 GB	3 h 26 min	C	Yes
Google Pixel 3 / Pixel 3 XL	5.5 in/ 6.3 in	Wi-Fi	64 GB	3 h 43 min	C	Yes
Samsung Galaxy Tab S2 8" (Only Model SM-T713)	8 in	Wi-Fi	32 GB	2 h 30 min	B	Yes
Panasonic Toughpad FZ-A2	10.1 in	Wi-Fi/LTE	32 GB	3 h 30 min	C	Yes



**Philips Ultrasound, Inc.**

22100 Bothell Everett Hwy

Bothell, WA 98021-8431, USA

[www.healthcare.philips.com/ultrasound](http://www.healthcare.philips.com/ultrasound)

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.

Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.  
239034 Rev X January 2020

# Lumify v1.8 Device Compatibility List



Panasonic Toughpad FZ-B2 (USB C models only)	7 in	Wi-Fi/LTE	32 GB	3 h 30 min	C	Yes
--	------	-----------	-------	------------	---	-----

\*Devices on the supported list satisfy all of the following criteria: 1. Passed software testing on the latest Lumify software, 2. Passed Philips hardware testing, 3. Device is readily available from the device manufacturer. Philips cannot test every device available on the market, so this list represents a range of devices to satisfy a variety of user needs with respect to form factor, price, battery life and connectivity options.



**Philips Ultrasound, Inc.**  
22100 Bothell Everett Hwy  
Bothell, WA 98021-8431, USA  
[www.healthcare.philips.com/ultrasound](http://www.healthcare.philips.com/ultrasound)

© Koninklijke Philips N.V. 2017. All rights reserved. Printed in USA.  
Reproduction or transmission in whole or in part, in any form or by any means, electronic, mechanical, or otherwise, is prohibited without the prior written consent of the copyright owner.  
239034 Rev X January 2020