



Philips
Ear thermometer

Accurate temperature reading
Fast measurement in 2 seconds
Syncs with Philips health app



DL8740/37

Connected temperature tracking

Syncs with free Philips HealthSuite health app

The connected ear thermometer allows you to quickly and accurately take readings and store them direct on your smart device using the free Philips HealthSuite health app so you can conveniently monitor your temperature over time.

Delivers fast and accurate results

- Accurate temperature reading in 2 seconds
- Provides accurate temperature measurements

Personalized tracking in the Philips health app

- Automatically syncs to the Philips health app via Bluetooth®
- Monitor your temperature over time

Comfortable and easy to use

- Suitable for people of all ages
- Gives results in Fahrenheit and Celsius
- Red LED display ring for rapid high temperature indication



asimpleswitch.com

PHILIPS

Ear thermometer

Accurate temperature reading Fast measurement in 2 seconds, Syncs with Philips health app

DL8740/37

Highlights

Phone compatibility



Automatically transfer your measurements to the Philips HealthSuite health app on your smart phone, so you can monitor changes over time. Compatible with iOS8 or higher/Android 4.4 or higher. For an overview of (in)compatible devices please have a look at www.philips.com/personalhealthcompatibility

Fast



Specially designed for speed and ease. Just one press of the button and you will have a reading in 2 seconds.

Easy and comfortable to use



The thermometer has a simple one-touch button and a probe tip that is comfortable to insert and suitable for use by all ages (newborn and up) No replacement caps required. Just clean the thermometer tip before usage.

Track your results over time



Determine your normal body core temperature so you can track and spot changes quickly and easily in the Philips HealthSuite health app.

Accurate



Body core temperature can be measured in five locations - the mouth, the ear, the armpit, the forehead and the rectum. In ear measurement (also known as tympanic measurement), along with rectum, is the most accurate form. The Philips ear thermometer uses this form to offer precise readings you can be confident are accurate. The thermometer offers an accurate measurement of $\pm 0.4^{\circ}\text{F}$ ($\pm 0.2^{\circ}\text{C}$) within the range of 90.3°F - 109.3°F (32.4°C - 42.9°C).



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO₂ emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Syncing

- Phone compatibility: iOS8 or higher / Android 4.4 or higher
- Data transfer technology: Bluetooth® Smart

Design

- Color: White
- Dimensions: 5.35(L)x1.36(W)x2.24(D)in
- Net. Weight: About 0.12lb

Ease of use

- Fast: 2 second measurement
- Clear LCD display
- Probe tip for all ages
- Can be used as stand-alone
- Switch to Fahrenheit & Celsius

Technical

- Accuracy: $\pm 0.4^{\circ}\text{F}$ ($\pm 0.2^{\circ}\text{C}$)
- Measurement range: 90.3°F - 109.3°F (32.4°C - 42.9°C)

Power

- Power source: 6V (2xAAA alkaline batteries)

Accessories

- 2x AAA batteries

Service

- 2-year guarantee



Issue date 2017-01-20

Version: 4.2.1

© 2017 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com