

PHILIPS

Home Care

Oxygen simplified

- reliable and cost-effective oxygen therapy

Philips EverFlo is a forward-thinking, efficient, compact and easy to handle oxygen product that save you time and money over the course of ownership.



AVAILABLE BY SUBSCRIPTION

innovation  you

Ideal oxygen concentrator for your patients



Philips Respiratory Home Care

Respiratory Therapist: 12 visits / year | 24x7 Technical support available

- Observation & Assessment • Pulmonary Rehabilitation • Patient and Family Education
- Regular Reporting of clinical to prescribing physician

Product Specifications

Input voltage	120 VAC +/- 10%
Avg power consumption	350 watts
Oxygen concentration* (at 5 lpm)	93% (+/- 3%)
Liter flow	0.5 to 5 liters per minute
Weight	31 lbs (14 kg)
Sound level	45 dBA typical
OPI alarm levels Low oxygen:	82%, very low oxygen: 70%
Operating humidity	Up to 95% relative humidity
Outlet pressure	5.5 PSI
Operating altitude	0 to 7,500' (0 to 2,286 m)

Product Benefits

- Compact and easy to move • Light weight: 14 Kg
- 93% purity at all flow rates from 0.5 to 5 LPM • Cost effective & Durable

Why Philips

Easy to store, transport, and deliver

Humidifier bottle holder compatible with all bottle styles

Internal tubing storage to keep tubing conveniently available for use

Metal cannula connector for greater durability

Interface panel with on/off switch and alarm indicators

Recessed flow meter reduces the risk of breakage



Philips Home Care Services India Pvt. Ltd.
9th Floor, DLF 9-B, DLF Cyber City,
DLF Phase - 3, Sector 25, Gurgaon - 122 002.

© 2009 Koninklijke Philips Electronics N.V.
All rights are reserved.

For more information contact
our call center : 18004197979