



The first choice for comfort and fit

Respironics PerforMax Full Face Mask for single-patient use (EE, SE)

The Respironics PerforMax Full Face Mask improves comfort, minimizes skin breakdown, and eliminates nasal bridge challenges by sealing around the perimeter of the face where, patients have less pressure sensitivity and smoother facial contours. The PerforMax Full Face Mask also simplifies mask selection and speeds setup by fitting a wide range of patients with only two sizes.

Key advantages

- Standard adult size fits most patients.
- Perimeter facial contact maximizes comfort and simplifies fitting by sealing around the flat surfaces of the face and avoiding the sensitive nasal bridge.
- The clear, unobstructed view provided by the PerforMax Full Face Mask helps reduce patient claustrophobia.

PHILIPS

Respironics PerforMax

The Respironics PerforMax Full Face Mask for single-patient use has additional advantages to support patient comfort and a quick, reliable fit.

- Interchangeable elbows accommodate quick transitions between single- and dual-limb ventilator circuits.
- With a choice of the standard elbow (SE) for dual-limb systems and the entrainment elbow (EE) for single-limb systems, clinicians can use the PerforMax Mask with any noninvasive ventilator.
- Snap clips require minimal effort to attach and remove headgear.
- EE radial diffusion leak ports reduce noise and deflect exhaled gas away from clinicians.
- Velcro® grab tabs allow quick, easy removal and adjustment of headgear.
- Integrated sizing gauge provides a convenient way to reduce cross-contamination.
- Low cost replacement headgear is available in two sizes.
- Digital Auto-Trak compatibility enhances NIV performance.
- Optional sealing pad minimizes air leaks and localized skin breakdown from nasogastric (NG) tubes.



Entrainment elbow (EE)



Standard elbow (SE)

Description	Quantity	Part Number
PerforMax Full Face Mask EE, adult	10	1052536
PerforMax Full Face Mask EE, small adult	10	1052535
PerforMax Full Face Mask SE, adult	10	1052548
PerforMax Full Face Mask SE, small adult	10	1052547
PerforMax entrainment elbow (EE) replacement kit	1	1045313
PerforMax standard elbow (SE) replacement kit	1	1043210
PerforMax small headgear	5	1051857
PerforMax small headgear	9	1045325
PerforMax large headgear	5	1051816
PerforMax large headgear	9	1045314
Accessory swivel replacement kit (22mm male-male)	1	7041
Nasogastric tube sealing pad	10	452040



© 2009 Koninklijke Philips Electronics N.V. All rights are reserved. PerforMax and Respironics are trademarks of Respironics, Inc., dba Philips Respironics, and its affiliates.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Healthcare is part of Royal Philips Electronics

www.philips.com/healthcare
healthcare@philips.com
fax: +31 40 27 64 887

MCI 4102353 PN 1062760
4522 962 51611 *JUL 2009

Philips Healthcare
Global Information Center
P.O. Box 1286
5602 BG Eindhoven
The Netherlands