

For every pressure you face

Philips non-invasive blood pressure solutions



More than blood pressure

The pressure for clinical excellence and efficiency

Philips latex-free blood pressure cuffs help you deliver customized care, including reusable and disposable cuffs with universal connectivity to Philips and non-Philips monitors.

Our high-quality blood pressure supplies are thoughtfully designed, validated, and manufactured to precise specifications to simply optimize performance throughout your operation to help you realize full accuracy, precise measurements, and smoother workflows. With updated industry standards for NIBP devices, you can stock fewer cuffs and adapters to meet all of your needs.

The pressure for improved patient care

We help you overcome the challenges you're facing, building in convenience, ease of use, and comfort to meet the needs of each patient. From color-coding to special artwork to extra-long indexing, we've made it easy for you to pick the right cuff for the right patient the first time, saving time and enhancing care.

With Philips, you have a choice of materials to balance patient comfort, cost, and durability, as well as a range of cuff sizes, pediatric through large adult, and DEHP-free and EVA materials on selected models.

Philips cuffs are designed with rounded edges and in a range of materials, including extra soft components to enhance patient comfort. We offer solutions to help prevent hospital-acquired infections (HAIs). Our reusable cuffs are easy to clean and validated with many different types of disinfectants. Our range of disposable cuffs also reduces the chance of spreading infection.

The pressure for cost-effective performance

High performance is cost-effective. Philips supplies and accessories are 100% guaranteed for quality, reliability, and durability, and tested well beyond the average life of a cuff, with up to a three-year warranty on certain cuffs.

Philips cuffs are easy to order and are backed with extensive clinical and technical support. The ease and assurance of a single, trusted source for all of your supply needs means that you can buy, stock, and standardize on one brand.





Backward compatible (with connector upgrade kits) to allow you to meet new safety requirements using older equipment.

Making safe connections

More hospitals trust Philips monitoring equipment with the care of their neonatal patients then all other manufacturers combined. As the leader in neonatal monitoring, we've introduced a line of neonatal blood pressure solutions designed specifically to prevent tragic misconnection errors that might harm your most vulnerable patients.

These neonatal blood pressure cuffs, hoses, and upgrade kits do not use Luer connectors, which the FDA and other leading organizations have said can pose a significant

safety hazard. The solution was created to prevent misconnections between blood pressure devices and other medical accessories, such as feeding tubes, catheters, IV tubing, and needleless IV ports.

Ideal for today's NICU

The neonatal blood pressure solution is intuitive and easy to use. Clinicians can deliver a high standard of care without the bright lights and loud noises that interfere with premature infants' neurological development. The cuffs are made of comfortable materials to be sensitive to the needs of these most fragile patients.

Quality cuffs for every clinical application and patient

Philips blood pressure solutions are simple to use, perform as expected, and save time.

Our reusable cuffs are easy to care for with a broad line of disinfectants, helping to prevent the spread of HAIs.

Easy Care for an easier day

With rounded edges for greater comfort, and designed for use in most areas of the hospital, reusable Easy Care cuffs are durable, come in a wide range of sizes, and are treated with Micropel $^{\text{TM}}$, an antimicrobial agent which prevents the growth of fungal, bacterial and algal organisms on the cuff.¹

Traditional two-piece solution

Ideal for general care, spot checks, and step-down units, Philips reusable Traditional cuffs feature a two-piece design with removable bladder for easy cleaning.

Comfort that doesn't quit

Quality reusable Comfort cuffs provide a practical alternative to traditional multi-patient cuffs and are a good choice for busy clinical units that use cuffs frequently. The cuffs feature soft material for outstanding patient comfort, while offering the economic benefits of multi-patient use with an easy-to-clean, waterproof design that resists stains and odors.



Philips reusable cuffs include Easy Care, Traditional, and Comfort.

¹The mechanism of action of this agent does not contribute to the development of drug-resistant microbial strains. Reduction in colonization or microbial growth on the device has not been shown to correlate with a reduction in infections in patients. Clinical studies to evaluate reduction in infection have not been performed.



Philips Soft and Traditional disposable cuffs offer greater comfort, durability and ease of use.

Philips disposable cuffs help prevent HAIs, and are designed to last well beyond the average patient stay so that clinicians don't have to apply more than one for each patient, which can allow more time to focus on patient care.

Greater comfort and ease

A soft, felt-like cover brings exceptional comfort to economical Soft disposables. These cuffs come with an integrated bladder, making them ideal for use in the OR, ER, cath lab, isolation, or other applications that call for single-patient use. The cuffs feature color-coding and exterior size markings for ease of use, and come in a full range of AHA-compatible adult, pediatric, and extralong sizes. Philips Soft disposable cuffs are also easier on the environment because they do not contain PVC.

Durable and long-lasting

Philips Traditional single-patient cuffs are an ideal choice for areas where risk of cross-contamination is high, or where cuffs are subject to heavy soiling, such as the OR and ER. Available in a full range of AHA-compatible sizes, you can count on the durability of these cuffs for the entire length of a patient's stay.

Universal performance

The Philips system of cuffs and easily secured adapters allow for universal connectivity to other manufacturers' NIBP monitors. Philips blood pressure cuffs include both single- and double-hose models to help you buy, stock, and standardize on one brand with confidence.



Philips gives the you ease and assurance of a single, trusted source for your supply needs.

Philips Healthcare is part of Royal Philips Electronics

How to reach us

www.philips.com/healthcare healthcare@philips.com

Asia +49 7031 463 2254

Europe, Middle East, Africa +49 7031 463 2254

Latin America +55 11 2125 0744

North America +1 425 487 7000 800 285 5585 (toll free, US only)

Please visit www.philips.com/healthcarestore



© 2011 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Printed in The Netherlands. 4522 962 69561 * JAN 2011