






# Rigorous testing for reliability you can count on.

Philips is a trusted name in oxygen therapy, providing the highest quality solutions for oxygen patients for more than 35 years.

So you can expect the Philips **SimplyGo Mini** to be an incredible, reliable machine.

In fact, every Philips **SimplyGo Mini** is subjected to the most rigorous testing, including impact, vibrations and temperature extremes. That means.....drop it, shake it, accidentally leave it in the car on a hot summer day....Philips made sure it can stand up to the abuses of everyday life.





## Portable oxygen therapy the way you want it.

The Philips **SimplyGo Mini** is sleek and lightweight... a tiny powerhouse that's easy to use, and simplifies your life by making its own oxygen from the air around you! And because it's a tank-free delivery system with no big, awkward tanks to hassle with, no oxygen deliveries to schedule..... the Philips **SimplyGo Mini** is virtually worry-free.

This machine features a sleek, non-clinical design, and its stylish carrying case is so good-looking that it won't ever draw attention to your oxygen use. That's freedom!

And, whether you decide to carry it by hand or over your shoulder, you can do it with complete comfort.

# \*Battery duration times are based on new, fully charged batteries used with a new Philips SimplyGo Mini system. Battery duration times will degrade with battery age, environmental operating conditions, use over time, and operational condition of the concentrator. Call us Toll-Free at 1-800-258-7678

## Battery-powered to support your life – wherever it takes you.

Now, a long plane flight is within reach.....and all-day outings with the grandchildren can be worry-free.....because Philips **SimplyGo Mini** offers easy-to-change lithium ion battery options with standard AC and DC power cords that add to its convenience and flexibility.

The manufacturer of this POC has determined this device conforms to all applicable FAA requirements for POC carriage and use on board aircraft.



So, you can choose the power option that works best for your life and how you want to live it!

#### **Available battery options:**

- Standard battery for battery duration up to 4.5 hours\* (pulse setting of 2 at 20 BPM). Resulting in total Philips SimplyGo Mini weight of 5 lbs.
- Extended battery for battery duration up to 9 hours\* (pulse setting of 2 at 20 BPM). Resulting in total Philips SimplyGo Mini weight of 6 lbs.





## SimplyGo Mini at a glance

Want to learn more about Philips SimplyGo Mini and see it in action? Simply visit www.philips.us/mySimplyGoMini there, you'll find helpful videos and more information.

Call us Toll-Free at 1-800-258-7678





### Meet the SimplyGo

For those who prefer continuous flow oxygen delivery in a portable device..... The SimplyGo is Philips' first ever oxygen concentrator that offers a continuous flow mode.

In fact, this model offers three different oxygendelivery modes to suit your needs:



**Pulse mode** – detects breath intake and delivers the necessary volume of oxygen.



**Continuous flow** – uninterrupted oxygen delivered from 0.5 to 2 LPM.



**Sleep mode** – Pulse mode with increased trigger sensitivity for use during sleep. For peace of mind, SimplyGo switches automatically to 2 LPM Continuous flow if no breath intake is detected.



Call us Toll-Free at 1-800-258-7678

## Which Device is Right for Me?

Buying an oxygen concentrator doesn't have to be complicated, especially when we're here to help. Take a look to compare the key benefits and features of the Philips SimplyGo Mini and the Philips SimplyGo. Still not sure? Give us a call at 1-800-258-7678.



#### SimplyGo Mini

The Philips SimplyGo Mini is right for you if you have a prescription for pulse flow oxygen



#### **SimplyGo**

The Philips SimplyGo is right for you if you require continuous flow oxygen delivery or need a higher oxygen flow

#### Weight

5 lbs. (2.3 kg) with standard battery installed 6 lbs. (2.7 kg) with extended battery installed 10 lbs. (4.5 kg) with battery installed

#### **Typical Battery Duration**

Standard battery lasts up to 4.5 hours<sup>1</sup> (pulse setting of 2, 20 BPM)
Extended battery lasts up to 9 hours<sup>1</sup> (pulse setting of 2, 20 BPM)

3 hours<sup>1</sup> (pulse dose setting of 2) 0.9 hours<sup>1</sup> (continuous flow 2 LPM)

#### **Settings**

Five pulse dose settings

Three different oxygen delivery modes

#### **Travel**

Meets FAA requirements for travel

Meets FAA requirements for travel

#### **Accessory Options<sup>2</sup>**

- ✓ Shoulder Bag
- ✓ Accessory Bag
- √ Backpack
- ✓ Standard Battery
- ✓ Extended Battery
- ✓ External Battery Charger
- ✓ AC Power Supply
- ✓ DC Power Supply

- ✓ Mobile Cart
  - ✓ Accessory Bag
  - ✓ Carrying Case
  - ✓ Standard Battery
  - ✓ External Battery Charger
- ✓ AC Power Supply
- ✓ DC Power Supply

Call us Toll-Free at 1-800-258-7678

<sup>1</sup>Battery duration times are based on new, fully charged batteries used with a new system. Battery duration times will degrade with battery age, environmental operating conditions, use over time, and operational condition of the concentrator.

## **Accessories**Available In All Bundles

#### Designed for travel

Be simply confident on the go. Designed to make traveling with your new oxygen concentrator even easier.



#### **Shoulder Bag**

- · Lightweight and durable
- · Carry handle and adjustable padded shoulder strap

SimplyGo Mini shoulder bag (Brown)

SimplyGo Mini shoulder bag (Black)



#### **Accessory Bag**

· Weather resistant and durable

SimplyGo Mini accessory bag (Brown)

SimplyGo Mini accessory bag (Black)

### You have the power

You lead an active and busy life. Sometimes that means you need even more battery life to keep you on the go.



#### Standard battery

• For battery duration up to 4.5 hours\*\* (pulse setting of 2 at 20 BPM). Resulting in total SimplyGo Mini weight of 5 lbs.

Standard lithium ion battery

## Charge wherever you are

No more moving the cords from place to place. Keep an extra in your bedroom, living room, car or RV.

AC power supply







SimplyGo Mini AC power supply

SimplyGo Mini DC power supply



#### Extended battery\*

• For battery duration up to 9 hours\*\* (pulse setting of 2 at 20 BPM). Resulting in total SimplyGo Mini weight of 6 lbs.

#### Extended lithium ion battery

\* Included in Gold and Platinum packages

\*\* Battery duration times are based on new, fully charged batteries used with a new SimplyGo Mini system. Battery duration times will degrade with battery age, environmental operating conditions, use over time, and operational condition of the concentrator.

#### Cannula Starter Kit

The Cannula Starter Kit comes with three 4-foot cannulas, three 7-foot cannulas, and a Quick Start Guide for reference.



Cannula Starter Kit

Call us Toll-Free at 1-800-258-7678

## **Accessories**Available In All Bundles

## Ready for a stroll

This lightweight fold-up cart features oversized wheels for a smoother roll and excellent maneuverability over challenging terrain.

Sliding handle allows personalized adjustment and can hold the SimplyGo accessory bag\*

\*Accessory bag sold separately

Oversized 6" wheels for a smoother roll and excellent maneuverability



SimplyGo mobile cart

### Designed for travel

Wear it over the shoulder or attach easily to the mobile cart.



- · Lightweight and durable
- · Padded shoulder strap
- Zippered pouch to store an extra battery or other accessories

SimplyGo carrying case

#### Cannula Starter Kit

The Cannula Starter Kit comes with three 4-foot cannulas, three 7-foot cannulas, and a Quick Start Guide for reference.



Cannula Starter Kit

### You have the power

Our powerful lithium ion battery is fully charged from empty in only 2 to 3 hours, and weighs just 1-1/2 pounds.

 For battery duration up to 3 hours\* (pulse dose setting at 2)

#### SimplyGo lithium ion battery

\*Battery duration times are based on new, fully charged batteries used with a new SimplyGo system. Battery duration times will degrade with battery age, environmental operating conditions, use over time, and operational condition of the concentrator.

## Charge wherever you are

No more moving the cords from place to place. Keep an extra in your bedroom, living room, car or RV.

AC power supply DC power supply





SimplyGo AC power supply

SimplyGo DC power supply



© 2019 Koninklijke Philips N.V All rights are reserved. Trademarks are the property of Koninklijke Philips N.V. or their respective owners. Respironics Inc. is a Philips company.

Specifications are subject to change without notice.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Caution: US federal law restricts these devices to sale by or on the order of a physician.

1010 Murry Ridge Lane, Murrysville, PA 15668 800-345-6443 • 724-387-4000

www.philips.com/simplygomini

MCI 4108455 POC Portfolio Patient Brochure Dig, US Eng