





Coronary Laser Atherectomy Catheter

## GS1 BARCODE SHEET

### GS1-128 Barcode Format




(01)00813132020743(17)160910(10)XXXXXXXXXX




- 1 Application Identifier, AI (01), indicates that the GS1 AI data field contains a GTIN
- 2 GTIN-14
- 3 AI (17) indicates that the GS1 AI data field contains expiration date
- 4 Expiration Date (YYMMDD, example: 160910)
- 5 AI(10) indicates that the GS1 AI data field contains Lot number
- 6 Lot number (XXXXXXXXXX, 9 characters)

### GTIN-14 Barcode Format



00813132020743













- 1 Indicator Digit for different levels of product hierarchy (For product with package of 1 unit per box: the 0 indicator is for both pouch and box; For product with package of multiple units: 0 is for pouch and 2 is for box)
- 2 GS1 Company Prefix (SPNC Prefix: 081313202)
- 3 Item Reference Number (example: 074)
- 4 Calculated Check Digit Reference (calculation based on previous 13 digits to ensure global uniqueness).





Coronary Laser Atherectomy Catheter

## GS1 BARCODE SHEET

Model #	Configuration	GTIN Portion of the GS1 Barcode
		GTIN-14 BOX BARCODE
110-001	0.9mm OTW	 00813132020002
110-001-A	0.9mm OTW	 00813132025489
110-002	0.9mm X-80 OTW	 00813132020019
110-002-A	0.9mm X-80 OTW	 00813132025472
110-003	0.9mm RX	 00813132020026
110-003-A	0.9mm RX	 00813132025465
110-004	0.9mm X-80	 00813132020033

Model #	Configuration	GTIN Portion of the GS1 Barcode
		<b>GTIN-14 BOX BARCODE</b>
110-004-A	0.9mm X-80	 00813132025458
114-009	1.4mm RX	 00813132020040
114-009-A	1.4mm RX	 00813132025441
117-016	1.7mm RX	 00813132020057
117-016-A	1.7mm RX	 00813132025434
117-205	1.7mm E* RX	 00813132020064
117-205-A	1.7mm E* RX	 00813132025427
120-008	2.0mm RX	 00813132020071
120-008-A	2.0mm RX	 00813132025410

Model #	Configuration	GTIN Portion of the GS1 Barcode
		<b>GTIN-14 BOX BARCODE</b>
120-009	2.0mm RX	 00813132020088
120-009-A	2.0mm RX	 00813132025403

**Corporate Headquarters**








The Spectranetics Corporation  
 9965 Federal Drive  
 Colorado Springs, CO 80921  
 Tel: 719-447-2000 Fax: 719-447-2022  
 Customer Service: 800-231-0978










Coronary Laser Atherectomy Catheter

## GS1 BARCODE SHEET

Model #	Configuration	GTIN Portion of the GS1 Barcode
		<b>GTIN-14 POUCH BARCODE</b>
110-001	0.9mm OTW	 00813132020002
110-001-A	0.9mm OTW	 00813132025489
110-002	0.9mm X-80 OTW	 00813132020019
110-002-A	0.9mm X-80 OTW	 00813132025472
110-003	0.9mm RX	 00813132020026
110-003-A	0.9mm RX	 00813132025465
110-004	0.9mm X-80	 00813132020033

Model #	Configuration	GTIN Portion of the GS1 Barcode
<b>GTIN-14 POUCH BARCODE</b>		
110-004-A	0.9mm X-80	 00813132025458
114-009	1.4mm RX	 00813132020040
114-009-A	1.4mm RX	 00813132025441
117-016	1.7mm RX	 00813132020057
117-016-A	1.7mm RX	 00813132025434
117-205	1.7mm E* RX	 00813132020064
117-205-A	1.7mm E* RX	 00813132025427
120-008	2.0mm RX	 00813132020071
120-008-A	2.0mm RX	 00813132025410

Model #	Configuration	GTIN Portion of the GS1 Barcode
		<b>GTIN-14 POUCH BARCODE</b>
120-009	2.0mm RX	 00813132020088
120-009-A	2.0mm RX	 00813132025403

**Corporate Headquarters**

The Spectranetics Corporation  
 9965 Federal Drive  
 Colorado Springs, CO 80921  
 Tel: 719-447-2000 Fax: 719-447-2022  
 Customer Service: 800-231-0978

