


# Turbo-Power™


Laser Atherectomy Catheter

## GS1 AND HIBC BARCODE SHEET

### GS1-128 Barcode Format




(01)00813132020743(17)160910(10)XXXXXXXXXX




- 1 Application Identifier, AI (01), indicates that the GS1 AI data field contains a GTIN
- 2 GTIN-14
- 3 AI (17) indicates that the GS1 AI data field contains expiration date
- 4 Expiration Date (YYMMDD, example: 160910)
- 5 AI(10) indicates that the GS1 AI data field contains Lot number
- 6 Lot number (XXXXXXXXXX, 9 characters)

### GTIN-14 Barcode Format



00813132020743



- 1 Indicator Digit for different levels of product hierarchy (For product with package of 1 unit per box: the 0 indicator is for both pouch and box; For product with package of multiple units: 0 is for pouch and 2 is for box)
- 2 GS1 Company Prefix (SPNC Prefix: 081313202)
- 3 Item Reference Number (example: 074)
- 4 Calculated Check Digit Reference (calculation based on previous 13 digits to ensure global uniqueness).



The Spectranetics Corporation

9965 Federal Drive, Colorado Springs, CO 80921  
Tel: 719-447-2000 • Fax: 719-447-2022

# **Turbo-Power™**

Laser Atherectomy Catheter

## **GS1 AND HIBC BARCODE SHEET**

<b>Model #</b>	<b>Configuration</b>	<b>GTIN Portion of the GS1 Barcode</b>
		<b>GTIN-14 BOX BARCODE (Updated)</b>
423-050	2.3 mm catheter OTW 0.018 in. compatible	 <b>00813132022983</b>
420-050	2.0mm Turbo-Power OTW, Manual	 <b>00813132024833</b>