

Philips Vacuum Accessories

Model	Accessory Kit Retail stores	Authorized Service Agents - Optional Individual Accessories			
FC6141/01		9965 100 57121 Inner filter			
FC6162/02		4322 004 94361 Freudenberg filter			
FC6168/01		4322 004 94361 Freudenberg filter			
FC6402/01	FC8063/01 KIT	4322 004 94361 Freudenberg filter	4322 004 94311 Pad	4322 004 27471 Pads	
FC6409/01		4322 004 94361 Freudenberg filter	4322 004 94312 Pad	4322 004 27471 Pads	
FC7020/01	FC8056/01 KIT	4322 004 93841 D.Calc filter	4322 004 93831 Cleaning pad		
FC7025	FC8057/01 KIT	4322 004 93841 D.Calc filter	4322 004 25641 Cleaning pad		
FC7028/01	FC8057/01 KIT	4322 004 93841 D.Calc filter	4322 004 25641 Cleaning pad		
FC8383/01	FC8074/01 KIT	4322 004 94331 Inlet filter	4322 004 94341 Super clean air filter	8838 021 03010 FC8021/2	FC8019 Classic S-bag (paper bag)
FC8471/01		4322 004 93801 HEPA10	4322 004 93811 Sponge	4322 004 93911 Exhaust Filter Sponge	FC8058 KIT
FC8588/01	FC8059 KIT	4322 004 94021 HEPA 13 filter	4322 004 93931 Inlet filter	8838 022 04010 FC8022 Clinic S-bags 4 pcs	FC8019 Classic S-bag (paper bag)
FC8671/01		4322 004 93801 HEPA10	4322 004 93811 Sponge	4322 004 93821 Micro filter	FC8058 KIT
FC9713/01		4322 005 33153 Inlet filter assy	4322 004 93791 HEPA10	4322 004 93350 FC8038 HEPA 13	

Philips Accessory Kit

4 microfiber pads

Compatible with PowerPro Aqua
Replace every 6 months

FC8063/01



Philips Performer Active Replacement Kit

6 x dust bags (s-bag® CLP)

1x Allergy filter
1 x triple-layer motor filter

FC8059/01



Philips Replacement Kit

3x dust bags (s-bag® CLP)

1x Super Clean Air filter
1x motor inlet filter
1x Lemon scent (s-fresh®)

FC8074/01



Philips Replacement Kit

2 microfiber pads

1 active calc filter

FC8057/01



Philips Steam Plus Sweep and Steam Cleaner

2 washable micro-fiber brushes

1 active calc filter
Replace every 6 months
Compatible with FC7020...21

FC8056/01

