



2015, April 16th

Beato Angelico at Museo Diocesano in Turin



During the Holy Shroud exposition, from 16 April to 30 June 2015, the Beato Angelico famous painting, Lamentation of the Dead Christ (tempera on panel, 152 x 200 cm, 1436), in extraordinary loan from the Museum of San Marco in Florence. is exhibited in a special display at Diocesan Museum of Turin. For the illumination of this important masterpiece, Ilti Luce used the new shaper, Kalibro M HP track GEN.2, an evolution in terms of design and lighting performance from the previous version.



Kalibro M HP track gen.2