


MASTER TL-D Food Pro (Butcher)

The colour 79 butcher tube. Next Generation triphoshor tube designed specifically for the presentation of fresh meat and poultry. Lamps emphasis the deep red colours such as those of meat, without the lamp having a red appearance.



Product ID	A max.	B max.	C max.	D max.
18W	589.8	594.5	604.0	28
30W	894.6	899.3	908.8	28
36W	1199.4	1204.1	1213.6	28
36W	970.0	974.7	984.2	28
58W	1500.0	1504.7	1514.2	28

Dimensions in mm

MASTER TL-D Food Pro (Butcher)

Shortcode	Wattage (W)	Cap Base	Luminous Flux Lamp (lm)	Average Lifetime EL (hrs)	Colour Rendering Index (Ra)	Colour Temperature (K)	Min Order Quantity	Carton Quantity
TLD18W79	18	G13	740	20000	74	3800	25	25
TLD30W79	30	G13	1300	20000	74	3800	25	25
TLD36W79	36	G13	1750	20000	74	3800	25	25
TLD36W791M	36	G13	1600	20000	74	3800	25	25
TLD58W79	58	G13	2800	20000	74	3800	25	25