

## Nguyen Lam Dang, Germany - October 16th 2017

I started my new tank with the Philips Coralcare on June 6th. The plan was to set up a nano reef tank stocked with colorful Zoanthus and LPS corals only. I started filling the tank with RODI water, mixed in the salt and after 4 hours the tank was fully stocked with many corals and two fishes.

Please find the main specs of the tank below:

- 95x45x45cm (including integrated technical compartment on the right side)
- Philips CoralCare
- Aqua Medic DC 1.1
- Coral Box s150

I have tried several other LED fixtures on my other tanks before, so LED is not a new thing to me. The Philips CoralCare is heavy and massive. Installation and programming were easy and self explanatory. After almost half year I am still very pleased and impressed by the Philips CoralCare. The light distribution and strength of the LEDs are great. Even with 10% total power coral growth is great.

